



Specification

Consolidated VOXP

© 2007 by VIASYS® Respiratory Care Inc.

All rights reserved. No part of this document may be reproduced, stored, or transmitted in any form or by any means without the written permission of VIASYS Respiratory Care Inc., except as specifically permitted under the copyright laws of the United States of America.

AVEA is a registered trademark of VIASYS Respiratory Care Inc. VIASYS and the VIASYS logo are registered trademarks of VIASYS Healthcare Inc. in the United States and other countries. Microsoft and Windows are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. Other brands and their products are trademarks of their respective holders and should be noted as such.

U.S.A.

VIASYS Respiratory Care Inc.
1100 Bird Center Drive
Palm Springs, California 92262-8099

Telephone: (800) 231-2466
(1) (714) 283-2228
Fax: (1) (714) 283-8493

Authorized European Representative

VIASYS Healthcare GmbH
Leibnizstrasse 7
97204 Hoechberg
Germany

Telephone: (49) (931) 4972-0
Fax: (49) (931) 4972-423

www.viasyshealthcare.com

Contents

Introduction.....	1
1.1 Purpose	1
1.2 Scope	1
1.3 Definitions/Abbreviations	1
1.4 References	1
1.5 Overview / Design Considerations.....	2
2 Use Cases	2
2.1 Immediate Patient Care	2
2.2 Link Control	2
2.3 Service / Maintenance	2
2.4 Asset Management.....	3
2.5 Manufacturing.....	3
3 Physical Requirements	3
4 General Characteristics	3
5 Message Format.....	3
5.1 Profile	4
5.2 Profile Info	8
5.3 Configuration	9
5.4 Link Control	10
5.5 Data Exchange	12
5.6 Data Snapshot.....	14
6 State Diagrams and Message Sequence.....	16
6.1 Protocol States	16
6.2 Discovery.....	22
6.3 Keep-Alive Exchange	23
6.4 Passive Mode Exchange	23
6.5 Active Mode Exchange.....	25
6.6 Log Query, Passive Mode	26
6.7 Query, Active Mode	27
6.8 Update Profile.....	28

6.9	Update Configuration	28
6.10	Link Restart.....	28
7	Nomenclature	30
7.1	Operational Settings.....	30
7.2	Alarm Limit Settings	31
7.3	Session Settings	32
7.4	Patient Information	32
7.5	Digital Monitors	32
7.6	Trended Parameters	34
7.7	Alarms	36
7.8	Scalars	42
7.9	System Information	43
8	Constant Definitions	43
9	Appendix A – CRC Definition and Algorithm	44
10	VIASYS VOXP Compliant Ventilators and Supported Data Set	45
10.1	Ventilator Configuration / Setup.....	45
10.2	VOXP Message, Syntax, and Sequence Description	46
10.3	VOXP Data Support.....	50
10.4	Info Class.....	51
10.5	Setting Class.....	54
10.6	Monitor Class.....	83
10.7	Alarm Class	101
10.8	Scalar Class.....	141
11	Foreign Language Support.....	145
12	VOXP Profile, Config, Data, and Link Messages.....	145
12.1	Profile Message – AVEA Adult	145
12.2	Profile Message – AVEA Pediatric.....	156
12.3	Profile Message – AVEA Neonate	167
12.4	Profile Message – VELA	179
12.5	Config Message.....	191
12.6	Data Message – AVEA	195
12.7	Data Message - VELA	197
12.8	Link Commands.....	198

13 Revision History..... 199

13.1 Rev A 199

13.2 Rev B 199

13.3 Rev C 207

Consolidated VOXP

Introduction

1.1 Purpose

This document provides the specifications of a flexible, extensible data communication protocol for interchange of digital data with VIASYS Healthcare ventilators. It is intended for software engineers developing software for VIASYS ventilators and any third party who may wish to develop host applications which exchange digital data with VIASYS ventilators.

1.2 Scope

The information contained in this document applies to any development, ventilator or host, which uses the VIASYS Comm (VOXP) Protocol for digital communication.

1.3 Definitions/Abbreviations

This section identifies terms, abbreviations and acronyms which are used in this document that are not commonly known or otherwise require clarification.

XML: Extensible Markup Language.

VOXP: VIASYS Open XML Protocol used to communicate information from a VIASYS Respiratory Care medical device to an external host system.

Big-endian: Big-endian and little-endian are terms that describe the order in which a sequence of bytes are stored in computer memory. Big-endian is an order in which the "big end" (most significant value in the sequence) is stored first (at the lowest storage address).

10BaseT: 10Mbps baseband data transmission over twisted-pair copper wire.

CRC: Cyclic Redundancy Check is a number derived from, and stored or transmitted with, a block of data in order to detect corruption.

ER: Engineering Report.

MIB: Medical Information Buss – see IEEE 1073 Medical Information Buss.

RJ: Registered Jack (e.g. RJ-45, RJ-11, etc.).

RS: Recommended Standard (e.g. RS-232, RS-422, etc.)

RxD: Receive data.

STP: Shielded Twisted Pair.

TxD: Transmit data.

VOXP: VIASYS Open XML Protocol

1.4 References

Identified below are documents which are referenced throughout or which must be used to interpret this document:

Extensible Markup Language (XML) V1.0

ISO 8601:2004 - Representation of Dates and Times

1.5 Overview / Design Considerations

This document defines digital data communication between a VIASYS ventilator and a host computer system. This specification is intended to provide the flexibility required to support any ventilator model and configuration as well as any host application. To achieve this flexibility, the protocol is based on a point to point, stream oriented exchange of data using only the printable ASCII character set. Messages are formatted using a subset of XML to provide a standardized method of structuring data. Instrument capability is negotiated with the host upon connection to allow generalization of the host compatibility across ventilator product models.

The document presents the specification of the protocol in a number of sections. First, a section describing the probable host applications will help ensure that the protocol will meet the needs of the uses identified. Secondly, the characteristics of the physical connection are described followed by a section with definitions that apply to multiple other sections. The syntax of properly formatted messages is defined along with the sequence that must be exchanged between the ventilator and host. Finally, a lexicon of nomenclature is presented that provides the identification and description of data entities that may be exchanged using the protocol.

NOTE: CLIENT SYSTEMS THAT SUPPORT THE VOXP PROTOCOL SHOULD BE DESIGNED IN SUCH A WAY THAT THESE AND FUTURE ENHANCEMENTS DO NOT RESULT IN FAILURE OF THE SYSTEM TO FUNCTION PROPERLY.

For example, some ventilators may support the “query-all” command and others may not, or some may support the EtCO₂ semantics and others not. The presence or absence of these protocol enhancements should not necessarily cause failure of other client system capabilities.

Also, the commands and semantics specified in this document were not supported in previous versions of the VOX Protocol.

2 Use Cases

2.1 Immediate Patient Care

- Integration of patient data with a vital signs monitor at the bedside.
- Central station monitoring of patient data.
- Remote Alarm monitoring – outside room, central station, mobile with care provider.

2.2 Link Control

- Automatic recording of Flow Sheet.
- Record of Alarm Log in central patient record.
- Record of Event Log in central patient record.
- Locally recorded patient chart information.

2.3 Service / Maintenance

- Fault Log.

2.4 Asset Management

- Asset location tracking.
- Asset use tracking.
- Service tracking.

2.5 Manufacturing

- Automated test.

3 Physical Requirements

This protocol is intended not to be constrained to any particular physical link. Any point-to-point connection from a ventilator to a single host, logical or physical, that can support the General Characteristics, should be adequate. This may include, for example, RS-232, TCP/IP, USB, etc.

4 General Characteristics

This section defines characteristics of this communication protocol which apply generally.

Only printing ASCII characters shall be transmitted in either direction between the ventilator and host. However, whitespace (SP, TAB, CR, and LF) shall be allowed, but should be ignored by the receiver. Also, XON/XOFF shall be allowed for flow control only. Note, however, that the use of XON/XOFF while the keep-alive mechanism is activated can cause the link to disconnect and return the session to the discovery state. This will occur if XOFF lasts longer than the expected 5 second link.ack() reply message.

The transmission of characters shall be stream oriented with XML element tags being the delimiters of a message and values shall be transmitted in big-endian order.

Consumption of communication channel bandwidth is the responsibility of the host application. It should not configure more data access than the channel will support.

Each message shall assume but not include the following XML Declaration statement:
`<?xml version="1.0" encoding="UTF-8" ?>`

Protocol extensions (e.g., additional data classes or types) may be added in the future. Though “backwards compatibility” is the goal of all protocol additions, client applications must be designed in a manner that allows them to ignore content that is syntactically correct but not recognized; otherwise, a new data class, message type or alarm level, for example, may cause the application to fail because it isn’t recognized.

5 Message Format

Messages consist of a well-formed XML Element, which may contain child elements. Whitespace is optional and should be ignored in a received message. The following examples have ample whitespace to more clearly depict the message structure.

Four primary message types are defined: Profile, Config, Link and Data. Profile and Config types are used to negotiate the capabilities of the ventilator and capabilities the

host intends to use. Data messages transport the requested data between ventilator and host. Link messages stimulate, acknowledge or control in some way the exchange of all the messages of other types.

Each message type (i.e. Profile, Config, Link, and Data) may contain a message identifier (msgID) element attribute. The msgID shall be a 16-bit value encoded as four ASCII HEX characters that the message exchange originator can use to uniquely identify a response message containing the echoed msgID. It may be implemented as a fixed value for each exchange sequence or, may be implemented as a simple counter across all messages that are sent (i.e. one counter for all message types). In the case where only one exchange is active at a time, the simpler fixed msgID approach may be used. However, when multiple exchanges may be active in parallel, the msgID mechanism provides the means for sorting out which response (esp. ACK and NAK messages) are associated with which command. A message with a negative acknowledgement (NAK), shall echo the msgID in its message to the originator. For example, if a link (send-profile) message originated by the Host and has an msgID="5A5A", and the ventilator identifies that the message has an error, the ventilator will respond with a link(nak) with an msgID="5A5A". The table below summarizes the VOXP message exchange between various command types:

VOXP Message Exchange					
Version Released	Command Type	Response Type	Originator*	Responder	Status
1.0	link(ping)	link(ack)	M/V	M/V	Mandatory
1.0	link(send-profile)	profile	M	V	Mandatory
3.1	link(send-profile-info)	profileInfo	M	V	Optional
3.1	link(set)	link(ack)	M	V	Optional
1.0	config	link(ack)	M	V	Mandatory
1.0	link(query)	data	M	V	Mandatory
3.1	link(query-all)	dataSnapshot	M	V	Optional
1.0	data**	link(ack)	V	M	Mandatory
1.0	link(restart)		M/V		Mandatory

* M=Manager (or external host); V=Ventilator

** Active mode only

5.1 Profile

The profile message identifies the capabilities of the ventilator to the host system and is signified by the XML element tag "<profile>". The profile element contains a number of child elements which describe the attributes of each data unit available for exchange and are signified by the XML element tag "<unit>". Unit elements are presented in Attribute Normal Form. The order in which the unit elements appear establishes the order in which the ventilator will send data. For sets of data items that have a priority relationship (e.g., alarm indicators), higher priority units precede lower priority units.

Units that represent enumerated values have a data type of ENUM and may contain one or more child elements that map a specific enumerated value to a textual label, signified by an <enum> XML element tag. Enum elements are presented in Attribute Normal Form.

Profile Element Attributes

- | | | | | | |
|----------------|--|----------|---|----------|-------------------|
| model | Identifies the basic ventilator name and defines the context for the subsequent profileVersion attribute value. For example: model="Avea". The value of this attribute is defined in the ventilator's VOXP "SysInfoModel" parameter. | | | | |
| profileVersion | The version number (major.minor) associated with the configuration of the subsequently defined information profile. Updates to the major number indicate either significant additions and modifications of content or backwards compatibility issues. Minor version updates indicate relatively small additions or modifications to the information content. (See also voxpVersion attribute.) The value of this attribute is defined in the ventilator's Major/Minor Version Number Assignments tables (see Section 10.2). | | | | |
| textEncoding | <p>Specifies the character encoding for the TEXT attribute data type (see below). This attribute shall comply to the character set encodings supported by XML, as defined by IANA (cf. www.iana.org/assignments/character-sets). Supported character sets include:</p> <table border="0" style="margin-left: 40px;"> <tr> <td style="padding-right: 20px;">"UTF-16"</td> <td>16-bit characters (default, if not otherwise specified)</td> </tr> <tr> <td>"UTF-32"</td> <td>32-bit characters</td> </tr> </table> <p>Unicode character sets are per Unicode Version 4.0.0 or more recent.
For example: textEncoding="UTF-32"</p> | "UTF-16" | 16-bit characters (default, if not otherwise specified) | "UTF-32" | 32-bit characters |
| "UTF-16" | 16-bit characters (default, if not otherwise specified) | | | | |
| "UTF-32" | 32-bit characters | | | | |
| voxpVersion | Version number (major.minor) of the baseline VOX Protocol specification is used by this device. Major number changes typically indicate syntactic changes to the protocol that may require parser modifications on the client application receiving the message. Minor number changes typically indicate content changes, such as additional nomenclature terms that are available. The value of this attribute is defined in the ventilator's Major/Minor Version Number Assignments tables (see Section 10.2). For example: voxpVersion="1.2" | | | | |

Unit & Enum Element Attributes

- class:** “setting” | “monitor” | “alarm” | “scalar” | “info” | “trend” | “log” | “pragma”
 The “class” attribute provides a grouping of unit data by similar characteristics. All unit data available for a class, except “trend” and “log” classes are grouped together and transmitted in a single message. Data in “trend” and “log” is record oriented and transmitted one record at a time. Any unit in the “pragma” class is generally for special functions or is instrument specific and should be ignored by the host.
- ID:** May be any of the Parameter ID’s from the Nomenclature (see §7).
 Parameter identification establishes a common reference of meaning and intent for unit data between the ventilator and host.
- type:** BOOL | BYTE | UBYTE | WORD | UWORD | INT | UINT | FLOAT | TEXT | ENUM

BOOL: A Boolean value represented as “0” | “1”.

BYTE: Eight bit signed integer, represented in two's complement hexadecimal “00”..”FF”.

UBYTE: Eight bit unsigned integer, represented in hexadecimal “00”..”FF”.

WORD: Sixteen bit signed integer, represented in two's complement hexadecimal “0000”..”FFFF”.

UWORD: Sixteen bit unsigned integer, represented in hexadecimal “0000”..”FFFF”.

INT: Thirty-two bit signed integer, represented in two's complement hex. “00000000”..”FFFFFFFF”.

UINT: Thirty-two bit unsigned integer, represented in hexadecimal “00000000”..”FFFFFFFF”.

FLOAT: Single precision IEEE floating point, represented in hexadecimal “xxxxxxx”.

TEXT: A NULL terminated character string represented as a sequence of ASCII Hex digits. The number of ASCII Hex digits per character depends on the character encoding specified for the profile (See “textEncoding” attribute above). For example, if UTF-16 is specified, then each character requires 16 bits: “ABC” = “0041004200430000”.

ENUM: A UWORD type with an associated label (see <enum> element).

scale: "E-99" .. "E+99"

"E" signifies the exponent of 10 by which the transmitted data is multiplied. If this attribute is omitted, the scaling is assumed to be 10^0 . E.g. "E+2" is $10^2 = 100$. In this case, transmitted data is represented x100, therefore should be divided by 100 to return to actual scale.

range: A pair of values, separated by ":", of the type and scale specified by the "type" and "scale" attributes respectively. The left value indicates the minimum valid value of the unit data and the right value indicates the maximum valid value of the unit data transmitted. E.g. "0064:03E8" specifies a valid range of 100..1000 for a WORD parameter.

resolution: A value of the type and scale specified by the "type" and "scale" attributes respectively. The value indicates the coarseness of the unit data values. E.g. "0005" specifies a WORD parameter should only have values in multiples of 5 (0, 5, 10, 15, etc.). This value is only relevant to the "settings class". The "monitor class" resolution is one (1).

epoch: Specifies the time duration of an epoch of scalar data in milliseconds. Scalar data is continuous, but must be transmitted at discrete intervals, so samples are collected over a time period (epoch) and transmitted as a unit. The value is expressed as a hexadecimal value. E.g. "01F4" specifies each epoch of scalar data contains 500ms of data samples.

size: Specifies the number of data samples per epoch of scalar data. Expressed as a hexadecimal value.

level: HIGH | MED | LOW | ALERT

Indicates the severity of Alarm unit data¹. Priority is determined by the order in which Alarms appear in the profile, highest priority first.

units: A string that may be displayed to identify Units of Measure of the associated data. Represented as type TEXT. E.g. "mL" would be represented as: units="006D004C".

label: A string that may be displayed to capture the associated data. Represented as type TEXT. E.g. "Volume" would be represented as: label="0056006F006C0075006D0065".

value: Identifies a specific enumerated value to be associated with a text label.

¹ Consult the ventilator's operator documentation for the proper interpretation of the severity levels.

Example

```

<profile voxpVersion="3.0" profileVersion="1.2" textEncoding="UTF-16" msgID="5A5A">

  <unit class="setting" ID="SetFlowBias" type="WORD" scale="E+1"
    resolution="0001" range="0004:0032" units="004C002F006D0069006E"
    label="004200690061007300200046006C006F0077"/>

  <unit class="setting" ID="SetPatSize" type="ENUM"
    label="00500041005400490045004E0054002000530049005A0045002000530
    045004C004500430054">
    <enum value="0000" label="004E0065006F"/>
    <enum value="0001" label="005000650064"/>
    <enum value="0002" label="004100640075006C0074"/>
  </unit>

  <unit class="monitor" ID="MntrVte" type="INT" scale="E+8"
    range="00000000:3B8B87C0" units="004C" label="005600740065"/>

  <unit class="scalar" ID="WavePaw" type="WORD" range="FDA8:04B0"
    epoch="01F4" size="0032" units="0063006D00480032004F"
    label="005000610077"/>

  <unit class="info" ID="SysInfoTimeSession" type="UINT"/>

  <unit class="trend" ID="TrendVte" type="WORD" scale="E+1" range="0000:2344"
    units="006D004C" label="005600740065"/>

  <unit class="alarm" ID="AlarmPpeakHigh" type="BOOL" level="HIGH"
    label="004800490047004800200050007000650061006B"/>

  <unit class="log" ID="LogAlarm" type="TEXT"/>

  <unit class="pragma" ID="Password" type="INT"/>
</profile>

```

5.2 Profile Info

The “profileInfo” message is identical to the Profile message with the exception that it only provides basic profile information and not detailed data <unit> specifications. This message is sent in response to a link(send-profile-info) command. The main purpose for this message is to allow manager systems to quickly determine the VOXP profile supported by the ventilator without having to retrieve the entire Profile message (a difference of a few dozen characters vs. a 25k byte or greater message length).

The syntax for this message is the same as that for the <profile> message described above. For example,

```
<profile voxpVersion="3.1" profileVersion="1.3" textEncoding="UTF-16"
msgID="5A5A"/>
```

or

```
<profile voxpVersion="3.1" profileVersion="1.3" textEncoding="UTF-32"
msgID="5A5A"/>
```

5.3 Configuration

The “config” message is the host response to the profile. It is the means by which the host informs the ventilator of the data units it intends to use. Only those data units configured will be available. To configure a data unit, its “unit” element from the profile is returned with only the class and ID attributes.

In addition, the config message establishes the pattern of interaction between host and ventilator. In the ACTIVE mode, the ventilator sends new data as soon as it becomes available. In the PASSIVE mode, the host must request all data transactions pertinent to the profile. The mode applies to classes “setting”, “monitor”, “alarm” and “scalar”. All other classes are by host request only.

As in any communication between two or more devices, there may be the potential for insufficient bandwidth or delays of various types of messages due to the specifics in the transmission scheme for each communication device (i.e. ventilator).

Config Attributes

mode: “ACTIVE” | “PASSIVE”

ACTIVE: Ventilator sends new data when available.

PASSIVE: Ventilator remains quiet until host requests data.

Unit Attributes

class: Identifies the class of unit data that will be configured for availability to the host system. Must be of a class appearing in the profile.

ID: Identifies the specific data units that will be configured for availability to the host system. Must correspond to a unit appearing in the profile.

Example

This example configures all data specified in the profile example above.

```
<config mode="PASSIVE" msgID="5A5A">
  <unit class="setting" ID="SetVol"/>
  <unit class="monitor" ID="MntrVte"/>
```

```

<unit class="scalar" ID="WavePaw"/>
<unit class="info" ID="SysInfoTimeSession"/>
<unit class="trend" ID="TrendVte"/>
<unit class="alarm" ID="AlarmPpeakHigh"/>
<unit class="log" ID="LogAlarm"/>
</config>

```

5.4 Link Control

Link messages stimulate, acknowledge or control in some way the exchange of all the messages of other types. Link messages are always transmitted in Attribute Normal Form. In ACTIVE Mode (where the ventilator originates data messages) when a data message is sent by the ventilator, a timeout period is activated and a reply message with the msgID is expected within that timeout period. If the reply message is not received within the timeout period, the ventilator will send out a link(nak-timeout) message with the same msgID. The timeout period is defined as the following:

timeout period = (message size (bytes) * byte TxTime) + 2seconds

This timeout period will allow longer timeout period for slower BAUD rates. The host system will effectively have approximately 2 seconds to respond to the message.

Link Attributes

cmd: "ack" | "nak" | "ping" | "send-profile" | "send-profile-info" | "query" | "query-all" | "restart" | "set" |

ack:

Acknowledge successful receipt of a data message.

nak:

Negative acknowledges receipt of data message. Errors, etc.

ping:

Discovery or connection keep-alive.

send-profile:

Host request for ventilator to send profile message.

send-profile-info:

Allows a Host to request for only the profile version information; not the complete profile specification as provided by the "send-profile" message.

query:

Host request for data.

query-all:

Allows the host to request all data elements in every class configured to return in a single "dataSnapshot" message response.

restart:

Indication that link will change characteristics, i.e. disconnect, change in profile, change in config. Either the ventilator or host can initiate a link restart message. During this phase the ventilator will need some time to

reinitialize its system parameters, and the host system must wait until the Discovery phase has started (i.e. wait for the link ping messages subsequent to initiation of the link restart message).

set:

Allows the host to request the target device (ventilator) to configure operational link parameters such as the keep-alive period.

error: Associated with “nak”. Indicates type of error. Valid error values are as follows:

- “nos” – Non Otherwise Specified.
- “config” – Query before valid configuration.
- “crc” – CRC of data doesn’t match attribute value.
- “na” – Data not available.
- “seq” – Message not expected or out of sequence.
- “timeout” – Complete message not received before timeout end of period.
- “syntax” – There was a general protocol syntax error in a received message.
- “too many elements” – A received message contained too many XML elements.
- “too many element attributes” – A received message element contained too

many XML

Attributes:

- “message too long” – The overall length of a received message was too long
- “invalid message type” – The message type (e.g., “config”) is unsupported
- “invalid link command” – A “link” message command type was unrecognized
- “not configured” – a query for data that is not configured.
- “invalid config unit element” – one or more data items are not recognized as a valid config unit element.
- “value out of range” – a specified value (e.g., a target keep-alive period) is not in the allowable range.

class: Associated with “query”. Identifies class of data host is requesting. May be any class value specified in a profile.

ID: Required only to differentiate between “log” data units configured from the profile.

time: Indicates a time after which log records are reported. If omitted, all records are reported. Format follows the ISO 8601 standard for representing dates and times and should include full field separators (e.g. “2001-11-01T16:35:45”).

keepAlivePeriod:

Specifies the time period (in seconds) for the keep-alive interval once a connection is established. The default period is 10 seconds. The maximum period settable by this command is 999 seconds or 16 minutes, 39 seconds. The minimum is 5 seconds.

Examples

Following are examples of a variety of different Link messages.

```

<link cmd="ack" msgID="5A5A"/>
<link cmd="nak" msgID="5A5A"/>
<link cmd="nak" error="crc" msgID="5A5A"/>
<link cmd="ping" msgID="5A5A"/>
<link cmd="send-profile" msgID="5A5A"/>
<link cmd="send-profile-info" msgID="5A5A"/>
<link cmd="query" class="log" msgID="5A5A"/>
<link cmd="query" class="log" ID="LogAlarm" msgID="5A5A"/>
<link cmd="query" class="log" ID="LogAlarm" time="20041127T235807" msgID="5A5A"/>
<link cmd="query-all" msgID="5A5A"/>
<link cmd="restart"/>
<link cmd="set" keepAlivePeriod="999" msgID="5A5A />
<link cmd="ack" keepAlivePeriod="999" msgID="5A5A"/>

```

5.5 Data Exchange

The essence of digital communication is that data can be exchanged between the ventilator and a host system. Data messages are the transport mechanism that achieves this exchange. Other message types influence the content and format of data messages. Except for Log and Trend class data, all unit data for a configured class is appended one after the next in the body of the data message in the order they appear in the profile. For Scalar class data, all samples in an epoch for a data unit are appended, then the samples for the next data unit, etc. For Log and Trend classes, each data message includes one record from the identified unit.

Data Attributes

- class:** Identifies the class data contained in the body of the message. Must be of a class configured in the config message.
- ID:** Required only to differentiate between "log" data units in response to a query.
- time:** Indicates the time that a log record was recorded. Format follows the ISO 8601 standard for representing dates and times.
- rec:** Identifies individual record from a log. The format of the attribute value is "*n/m*" where *n* and *m* are both hexadecimal values and indicate the *nth* record of *m* records being reported.
- code:** If the data being transmitted is discrete, there may be an associated numerical code that is an alias to the data. E.g., an "Intubation" event may alternatively be referred to as "Event 128" as follows:
- ```

<data class="log", ID="LogEvent", time="20041127T235807", code="0080",
 rec="01/03">0049006E007400750062006100740069006F006E</data>

```

- crc:** Optional attribute for error detection in the data stream. CRC is of type UWORD indicating the 16-bit CRC of the ASCII characters representing the data in the body of the message. Its calculation is performed on the message encoded character content. For example, given the message `<data...>1234567890ABCDEF</data>`, the CRC will be calculated for the ASCII characters between the ‘>’ and ‘<’ characters. See Appendix A for a detailed description of the CRC algorithm used for this protocol.
- seq:** The seq attribute is a 16-bit ASCII HEX encoded value which identifies the scalar/wave epoch update sequence that is being returned. The seq increments monotonically once a connection is established and allows polled data retrieval to determine whether the same epoch has been retrieved or whether epoch updates have been missed (possibly due to the slow polling frequency).
- units:** A string that may be displayed to identify Units of Measure the associated data. Represented as type TEXT. E.g. “mL” would be represented as: `units="006D004C"`.
- label:** A string that may be displayed to capture the associated data. Represented as type TEXT. E.g. “Volume” would be represented as: `label="0056006F006C0075006D0065"`.

### Examples

Following are examples of a variety of different Data messages.

*Four settings of type WORD have been configured:*

```
<data class="setting">012E08A961F03010</data>
```

*Four settings of type WORD have been configured, with CRC:*

```
<data class="setting" crc="9F37">012E08A961F03010</data>
```

*Two monitors of type WORD and two of type INT have been configured:*

```
<data class="monitor">012E08A9000361F000203010</data>
```

*Two scalar waveforms and the scalar ‘Wave Metric’ of type WORD with 500ms epochs of 50 samples each have been configured:*

```
<data class="scalar" crc="CB01" seq="0070"
 msgID="0EBE">0021003E01220402060E06EB0749076C07710760074407240
 71106F406DB06BE06A00682066C064C062C061505F805DC05C705AF05970
 57E0568055305380523050D04F304D604B3049704840472045A04440428033
 D002EFE5DFE11FD45FC95FC00FB8902730275028202A70324033C037C03
 8403A703B003C503CD03DD03E803F60402040C04180422042B0438044404
 4E0454045D0468046E04770483048C0494049D04A804AF04B704BE04C404
 CB04D204D804DF04E204E704AD048F0495047A045D044304260032000200
 08002700400066007F00A600BF00E500FD01210139015C0174019601AC01C
 D01E3020302170236024A0268027B029802AB02C702D902F4030603200331
 034A035A03710381039803A703BD03CB03E003E803DF03D303C803BF03A
 B039B0383</data>
```

*An Alarm log has been configured:*

```
<data class="log" ID="LogAlarm" time="20041127T235807" code="0135"
 rec="03/1F">0053006F006D00650020004C006F006700200045006E0074007
 20079</data>
```

The above alarm log example contains a code="0135", which represents "Alarm 309" in the list of potential alarms associated with this device. The "03/1F" after "rec=". This represents the 3<sup>rd</sup> of 31 individual log entries that are being transmitted as part of the Alarm log. The information, there after represents the description for this log message (i.e. alarm log message) in ASCII text.

## 5.6 Data Snapshot

When the manager system seeks to capture a "snapshot" of the currently configured data – that is a single data update set – it can issue a link(query-all) and receive back the "dataSnapshot" message. This simplifies the communication requirements in that a single command is required vs. multiple query requests for each desired data class.

This message simply wraps one or more "data" messages as defined above, within a single "dataSnapshot" message element. The syntax for the contained "data" elements is identical to what would be returned from the ventilator if they were retrieved using single link(query) requests. The only difference is that if the triggering link(query-all) contained a message ID (msgID), this will be echoed back in the "dataSnapshot" element, and not in the contained "data" elements. All "crc" attributes shall be contained in the individual "data" class elements (see the example below).

Note that though this does result in a larger single message size (than individual link(query) commands), the data is guaranteed to be internally consistent and cohesive, addressing potential issues that may result from, for example, ventilator setting changes taking effect during multiple data class poll requests (i.e., if it takes 30 seconds to poll all the various classes of data, the monitored values may not match those indicated by the retrieved settings).

Example

```
<dataSnapshot msgID="0007" >
 <data class="setting" crc="67EA">
 0015001400000019000100C80000000F000F0006000300060000000C0005000
 50005002800000000004B00140000004B000A001E03E800000BB80000000000
 000014000300280008004B01F400320BB800000000000370096100000000D00
 000000007402F8004F0003</data>
 <data class="monitor" crc="3A99">
 8000800080000067800005C580008000800080000015FF9D000A80008000800
 00029001E000604218000023F00080000002480008000800000080008000800
 0800000608000800002A9004502FCC25500000000013F02A43E5D0260170B8
 000800000000000000002F52A7C8000800080000098</data>
 <data class="alarm" crc="059C">
 00
 10000000000011000000000010000</data>
 <data class="info" crc="CA65">
 00FF000100300030003000370041003100300030003500390042004200340037
 00300030000000410042005600300031003000310037000000320000000260BF
 0000</data>
 <data class="scalar" crc="B501" seq="0391">
 00
 00
 00
 000200
 00
 0000000000000010002
 0000000000000000FEFCFEF9FEF9FEFAFEF9FEFAFEFAFEFAFEF9FEFAFEF
 BFEF9FEF8FEFBFEFAFEFBFEFAFEFFFEFBFEF7FEF8FEF8FEF8FEFBFEFB
 FEFBFEFBFEFAFEF9FEF8FEF9FEF8FEF8FEFDFEFCFEFBFEFBFEF8FEF9F
 EF7FEF9FEFAFEF9FEF9FEF9FEF8FEF8FEFAFEFAFEF70110011701160115
 01170110010B0112011601150111010D0114010E0114011401100112010E011
 401190117011501140110011501110112011201110114011001150111010C011
 5011401150115010D011201170117010E011001110117011601150114FFFDFF
 FDFFEFFFEFFDFDFFEFFDFDFBFFF8FFFEFFDFDF8FFF8FFFCFFFCFF
 CFFDFDFFEFFBFFF6FFFDFFDFBFFDFDFCFFFEFFDFDFEFFFEB
 FFFDFDFBFFFCFFF9FFFDFFFEFFFEFFDFFEFFFCFFF2FFFAFFFAFFFEF
 FFEFFFEFFFEFFDFFEFFFE0022002200220022002200220022002200
 2200220022002200220022002200220022002200220022002200220022002200
 2200220022002200220022002200220022002200220022002200220022002200
 2200220022002200220022002200220022002200220022002200220022002200
 0129012D0132012E012A012C012D012C012C012D012E01290129012E012801
 28012B012D012B012C012E012B012C012E012D012A012B012E012B0129012
 B012D012801270129012E012D012C012B012A012B012A012B000000000000
 00
```

[illegible]

## 6 State Diagrams and Message Sequence

This section describes the sequence of messages necessary to manage the connection and exchange of data between ventilator and host. Once a logical connection is established, communication is effected with a succession of transactions, each consisting of a stimulus and response. Each transaction must complete before the next commences. Some common examples of transactions are:

```
Vent: link.ping() replied by Host: link.ack()
Vent: data(monitor) replied by Host: link.ack()
Host: link.query(settings) replied by Vent: data(settings)
```

The sequence diagrams, starting in section 6.2, assume some data is configured in all classes.

## 6.1 Protocol States

The next two sub-sections displays the various states of the VOX Protocol (Figure 1 VOXP Protocol State Model) as well as the states while monitoring a connection (Figure 2 Connection Monitoring State Model).




Figure 1 VOXP Protocol State Model

## Notes:

- Processing Error – Could be error indication from the communication port or a timeout in receiving an ACK/NAK to a sent message.

**VOX Protocol State Definitions**

<b>Inactive</b>	Start-up state for the protocol, in which all connections with the transport and port hardware are broken and another protocol may be running.
<b>Activating VOXP</b>	Performs all initialization necessary to prepare the protocol stack to receive and process commands from an external system. This includes data member initialization and configuration of the Comm I/O transport.
<b>Establishing Connection</b>	The VIConnection object is started and begins to try and establish a connection with an external host. See the object's description for more details
<b>Waiting for Profile Request</b>	Once a connection has been established, the protocol waits for the host to request the device's communications profile, a list of all available (and configurable) information and services.
<b>Waiting for Config Command</b>	The system is waiting for a configuration message, ignoring all other messages save keep-alive pings and ACK's.
<b>Configuring Link</b>	When a valid "configuration" message is received, the protocol stack configures itself accordingly (e.g., active or passive mode), and transitions to the appropriate operational state.
<b>Passive Mode</b>	In this mode, the system responds to requests but does not automatically send updates.
<b>Active Mode</b>	In this mode, the system both responds to requests, as well as automatically sends data updates as they become available.
<b>Restarting Link</b>	The connection is restarting and being reconfigured as a result of changes that may, for example, change the information that may be communicated to the external system.
<b>Shutting Down VOXP</b>	All VOXP communications activities are terminated and any acquired resources (e.g., communication port) released.

**VOX Protocol Event Definitions**

<b>Activate</b>	The protocol has been activated, either during startup initialization or dynamically by the operator.
<b>Commlo Ready</b>	The communication port has been configured for the VOXP protocol and is ready to send and receive application level messages.



<b>Connection Established</b>	The VIConnection object has successfully established a connection with a remote system.
<b>Profile Request Received</b>	A “profile request” message has been received.
<b>Configuration Command Received</b>	A “configuration command” message has been received.
<b>Configure Passive</b>	The external system has selected passive mode operation.
<b>Configure Active</b>	The external system has selected active mode operation.
<b>Configure</b>	A “configuration command” has been received.
<b>Connection Broken</b>	Indicates that an active connection has been broken.
<b>Processing Error</b>	A problem was encountered either processing VOXP messages or data update messages, in either case requiring the connection to be broken and reestablished.
<b>Mode Changed</b>	A new language, patient size, or communication parameters (e.g., BAUD rate) needs to be supported on the VOXP connection, requiring cycling of the connection.
<b>Link Restart</b>	A “link restart” command message has been received.
<b>Restart Complete</b>	The transport and protocol have been initialized and reconfigured.
<b>Protocol Shutdown</b>	A different protocol has been selected, requiring VOXP processing to be shutdown.




Figure 2 Connection Monitoring State Model

### Connection Monitoring State Definitions

<b>Inactive</b>	The state when the VOXP Protocol is also inactive.
<b>Monitoring Connection</b>	State during which the object is trying to establish and maintain a connection.
<b>Discovery</b>	For 2 minutes the object sends a “ping” message every few seconds until an external system responds with an “ack” to the “ping”.
<b>Waiting for Ping</b>	After the 2 minute discover period, the object stops sending “pings” and waits until it receives a “ping” message from an external system.
<b>Waiting for Link Config</b>	Once a connection has been established, the ventilator waits for the link to be configured (i.e., exchange of a “profile” and

“config” message pair, at which point it transitions to the keep-alive state. If 60 seconds elapses without receipt of any messages, the connection is reset. Depending on the selected BAUD rate and the length of the “config” message, completion of link configuration can require a significant period of time.

**Keeping Connection Alive** Once the link has been configured, the object sends periodic keep-alive pings, and waits to receive the “ack”. If an “ack” isn’t received in time, the object initiates a resetting of the connection.


**Resetting Connection** The object stops the timer, signals that the connection is closed, and transitions automatically to the Inactive state.

### **Connection Monitoring Event Definitions**

<b>Start Discovery</b>	Event from the VOX Protocol indicating that connection discovery and monitoring should be initiated.
<b>Send 5 second ping</b>	Timer event to send another ping.
<b>Ping received</b>	A link.ping() message was received.
<b>Ack received</b>	A link.ack() message was received.
<b>Reset connection</b>	Either an ack wasn’t received for the ping in the required time or the VOX Protocol has determined that the connection needs to be reset (e.g., a change in mode).
<b>After: X seconds</b>	If a message has not been received from the manager system within “X” seconds, a link.ping() message is issued. “X” defaults to 10 seconds but is configurable using the link.set_ka_period() function.


## 6.2 Discovery

When the ventilator is first powered on or when the VIASYS Comm (VOXP) protocol is selected as the communication protocol, the ventilator will attempt to discover a host. Every 5 seconds for the first 2 minutes, or until a reply is received, a link message with the ping attribute will be transmitted. Any reply will establish the link, but queries will be responded with link.nak() until configuration is complete. If connected to a host after the initial two minute period, the host may initiate the connection by sending a link.ping().




### 6.3 Keep-Alive Exchange


When the ventilator is first powered on or when the VIASYS Comm (VOXP) protocol is selected as the communication protocol, the ventilator will attempt to discover a host. Every 5 seconds for the first 2 minutes, or until a reply is received, a link message with the ping attribute will be transmitted. Any reply will establish the link, but queries will be responded with link.nak() until configuration is complete. If connected to a host after the initial two minute period, the host may initiate the connection by sending a link.ping().




### 6.4 Passive Mode Exchange



## 6.5 Active Mode Exchange




## 6.6 Log Query, Passive Mode







## 6.7 Query, Active Mode



## 6.8 Update Profile




## 6.9 Update Configuration



## 6.10 Link Restart

The following diagrams illustrate that a `link.restart()` command is exchanged (in either direction), the vent will go through a reconfiguration/initialization phase and eventually restart the connection discovery phase. The manager system should wait until the vent starts sending `link.ping()` messages before reconfiguring the link.



## 7 Nomenclature

This section defines several groups of identifiers in the nomenclature which are used in the VOXP protocol. However, the implementation specifics relating to any identifier in this nomenclature can be found in the target ventilator specification complying with the VOXP protocol. These details in the ventilator specifications include, but are not limited to, data type and domain values.

### 7.1 Operational Settings

	Parameter ID	Description
1	SetCO2MonitoringEnable	Enable or Disable CO2 Monitoring
2	SetEndTidalCO2Average	Number of breaths EtCO2 is averaged over.
3	SetFiO2	FiO2
4	SetFiO2AutoControlEnable	Activates and de-activates the Automatic FiO2 Control system.
5	SetFiO2IncreaseActive	Increase FiO2 Activated (SetFiO2 + SetIncrFiO2)
6	SetFlowBias	Bias Flow
7	SetFlowCycle	Flow Cycle % of PIFR
8	SetFlowCyclePsv	PSV Flow Cycle
9	SetFlowDemand	Intrabreath Demand Flow OFF/ON
10	SetFlowInsp	Inspiratory Flow
11	SetNebulizerActive	Nebulizer Active / Inactive
12	SetPanelLockActive	Panel Lock Active / Inactive
13	SetPauseInsp	Inspiratory Pause
14	SetPresHigh	APRV/BiPhasic High Phase Pressure
15	SetPresInsp	Inspiratory Pressure
16	SetPresInspNPPV	Inspiratory Pressure, Non-Invasive
17	SetPresLow	APRV/BiPhasic Low Phase Pressure
18	SetPresNasalCPAP	Nasal CPAP pressure level
19	SetPresPeep	Positive End-Expiratory Pressure
20	SetPresPsv	Pressure Level, Pressure Support Ventilation
21	SetPresPsvNPPV	Pressure level, Pressure Support Ventilation, Non-Invasive
22	SetRate	Breath Rate
23	SetRiseInsp	Pressure Control Rise Time
24	SetRisePsv	PSV Rise Time
25	SetRiseVsync	Vsync Rise Time
26	SetSpO2AlarmDelay	Determines the duration that a High or Low SpO2alarm or "SpO2 Invalid" alarm must persist before the alarm is activated.
27	SetSpO2AveragingTime	Determines the time period for averaging the output of the SpO2 readings.
28	SetSpO2Enable	Enable or Disables Pulse Oximeter related functions.
29	SetSpO2Mode	Configures the pulse oximeter module to a specified operating mode.

	Parameter ID	Description
30	SetTimeHigh	APRV/BiPhasic High Phase Time
31	SetTimeHighPsv	APRV/BiPhasic High Phase PSV OFF/ON
32	SetTimeHighSync	APRV/BiPhasic Low Phase Sync Window % of Time High
33	SetTimeInsp	Inspiratory Time
34	SetTimeLow	APRV/BiPhasic Low Phase Time
35	SetTimeLowSync	APRV/BiPhasic Low Phase Sync Window % of Time Low
36	SetTmaxPsv	PSV Maximum Inspiratory Time
37	SetTrigFlow	Flow Trigger Sensitivity
38	SetTrigPres	Pressure Trigger Sensitivity
39	SetVCO2Average	Patient's exhaled minute volume of CO2
40	SetVol	Tidal Volume
41	SetVolAssured	Machine Volume
42	SetVolLimit	Volume Limit
43	SetVolSigh	Sigh Volume Breath OFF/ON
44	SetVolWave	Volume Control Waveform.
45	SetVsync	Vsync mode OFF/ON

## 7.2 Alarm Limit Settings

	Parameter ID	Description
1	LimitApnea	Apnea Interval
2	LimitEndTidalCO2High	High End Tidal CO2 limit
3	LimitEndTidalCO2Low	Low End Tidal CO2 limit
4	LimitFiO2AutoHigh	Upper bound for the High Auto FiO2 Limit to trigger alarm.
5	LimitFiO2AutoLow	Lower bound for the Low Auto FiO2 Limit to trigger alarm.
6	LimitFiO2BaselineHigh	Upper bound for the Baseline FiO2 to trigger an alarm.
7	<i>LimitLowPpeak</i>	<i>Low Peak Airway Pressure – <b>DEPRECATED<sup>2</sup></b> (See LimitPpeakLow)</i>
8	LimitPeepLow	Low PEEP
9	LimitPpeakHigh	High Peak Airway Pressure
10	LimitPpeakLow	Low Peak Airway Pressure
11	LimitPulseRateHigh	Upper bound of pulse rate to trigger an alarm
12	LimitPulseRateLow	Lower bound of pulse rate to trigger an alarm
13	LimitRateHigh	High Breath Rate
14	LimitSpO2High	Upper bound SpO2 to trigger an alarm
15	LimitSpO2Low	Low bound SpO2 to trigger an alarm
16	LimitVeHigh	High Minute Volume
17	LimitVeLow	Low Minute Volume
18	LimitVteHigh	High Tidal Volume
19	LimitVteLow	Low Tidal Volume

### 7.3 Session Settings

	Parameter ID	Description
1	SetAAC	Automatic Airway Compensation OFF/ON
2	SetAltitude	Altitude at which the device is being operated
3	SetCircComp	Circuit Compliance Compensation Constant
4	SetEttDia	Endotracheal Tube Diameter
5	SetEttLen	Endotracheal Tube Length
6	SetFiO2Monitoring	FiO2 monitoring OFF/ON
7	SetHumidifier	Active Humidifier OFF/ON
8	SetLanguage	Local Language for User Interface
9	SetLeakComp	Leak Compensation OFF/ON
10	SetMode	Breath delivery mode used
11	SetModellv	Independent Lung Ventilation configuration
12	SetNebulizerTime	Time the nebulizer is active
13	SetPanelLockEnable	Enable/disable front panel user input lockout availability
14	SetPatSize	Patient size setting
15	SetPatWt	Patient weight setting
16	SetPresBaro	Barometric pressure setting
17	SetIncrFiO2	Added percentage of oxygen concentration delivered to the patient for a specified time based on the FiO2 setting.
18	SetSensitivityLowVte	The number of consecutive breaths that is in violation of the low Vte threshold.
19	SetSpO2TargetHigh	Auto FiO2 Cmd is greater than or equal to the preset High Auto FiO2 Limit for a period of 60 seconds or more.
20	SetSpO2TargetLow	Auto FiO2 Cmd is less than or equal to the preset Low Auto FiO2 Limit for a period of 60 seconds or more.
21	SetVeLowOffEnable	Low Ve Limit Alert OFF setting is OFF/ON

### 7.4 Patient Information

	Parameter ID	Description
1	PatInfolD	Patient identification

### 7.5 Digital Monitors

	Parameter ID	Description
1	MntrAutoPEEP	Airway Auto PEEP
2	MntrAutoPEEPdelta	Delta Airway Auto PEEP
3	MntrAutoPEEPesoph	Esophageal Auto PEEP
4	MntrPresBaro	Barometric Pressure
5	MntrC20	Compliance Ratio
6	MntrCcw	Chestwall Compliance
7	MntrCdyn	Dynamic System Compliance
8	MntrCdynNorm	Dynamic System Compliance, Normalized
9	MntrClung	Lung Compliance

	Parameter ID	Description
10	MntrVentilationCO2	The patient's exhaled minute volume of CO2, calculated over the "VCO2 Average" interval.
11	MntrCstat	Static System Compliance
12	MntrCstatNorm	Static System Compliance, Normalized
13	MntrEndTidalCO2	Patient's peak expired CO2 level as measured and reported by the CO2 analyzer.
14	MntrFiO2	FiO2 – Fractional Inspiratory Oxygen Concentration
15	MntrFiO2Baseline	Average FiO2 required in order to maintain the patient in stable normoxemia over time.
16	MntrIE	Inspiratory:Expiratory Time Ratio
17	MntrLeak	$(V_{ti}-V_{te})/V_{ti} \times 100$
18	<i>MntrIveTotalNorm</i>	<i>Minute Volume, Normalized - <b>DEPRECATED</b><sup>2</sup></i> (See MntrVeTotalNorm)
19	MntrMIP	Maximum Inspiratory Pressure
20	MntrNcpapMeanFlow	Nasal CPAP mean inspiratory flow
21	MntrNcpapPres	Nasal CPAP pressure
22	MntrP100	Respiratory Drive
23	MntrPair	Air Supply Pressure
24	MntrPawDelta	Delta Airway Pressure
25	MntrPeep	Positive End Expiratory Pressure
26	MntrPefr	Peak Expiratory Flow Rate
27	MntrPesDelta	Delta Esophageal Pressure
28	MntrPifr	Peak Inspiratory Flow Rate
29	MntrPmean	Mean Airway Pressure
30	MntrPO2	Oxygen Supply Pressure
31	MntrPpeak	Peak Airway Pressure
32	MntrPplat	Plateau Airway Pressure
33	MntrPplatPtp	Transpulmonary Pressure at Airway Pressure Plateau
34	MntrPtpPEEP	Transpulmonary Pressure at Auto PEEP
35	MntrRate	Total Breath Rate
36	MntrRateMand	Mandatory Breath Rate
37	MntrRateSpon	Spontaneous Breath Rate
38	MntrRimp	Imposed Resistance
39	MntrRlung	Lung Resistance
40	MntrRpeak	Peak Expiratory Airway Resistance
41	MntrRrs	Respiratory System Resistance
42	MntrRSBIndex	Rapid Shallow Breathing Index
43	MntrSpO2	Patient's SpO2 as measured and reported by the pulse oximeter.

<sup>2</sup> "DEPRECATED" indicates that this term is no longer supported in new VOXP software releases.

	Parameter ID	Description
44	MntrSpO2PerfusionIndex	Percentage of pulsatile signal to non-pulsatile signal (pulse strength).
45	MntrSpO2PulseRate	Patient's pulse rate as measured and reported by the pulse oximeter.
46	MntrTe	Expiratory Time
47	MntrTi	Inspiratory Time
48	MntrVentilationAnatomicalDeadSpace	The patient's anatomical dead space, measured on each breath, and averaged over the "VCO2 Average" time interval.
49	MntrVentilationAnatomicalDeadSpaceVtRatio	The patient's airway dead space to tidal volume ratio, measured on each breath, and averaged over the "VCO2 Average" time interval.
50	MntrVdel	Machine Delivered Volume
51	MntrVeSpon	Spontaneous Minute Volume
52	MntrVeSponNorm	Spontaneous Minute Volume, Normalized
53	MntrVeTotal	Minute Volume
54	MntrVeTotalNorm	Minute Volume, Normalized
55	MntrVtCO2	The patient's exhaled tidal volume of CO2, calculated on a breath-by-breath basis.
56	MntrVte	Exhaled Tidal Volume
57	MntrVteMand	Mandatory Tidal Volume
58	MntrVteMandNorm	Mandatory Tidal Volume, Normalized
59	MntrVteNorm	Exhaled Tidal Volume, Normalized
60	MntrVteSpon	Spontaneous Tidal Volume
61	MntrVteSponNorm	Spontaneous Tidal Volume, Normalize
62	MntrVti	Inspired Tidal Volume
63	MntrVtiNorm	Inspired Tidal Volume, Normalized
64	MntrWobImposed	Work of Breathing, Imposed
65	MntrWobPatient	Work of Breathing, Patient
66	MntrWobVent	Work of Breathing, Ventilator

## 7.6 Trended Parameters

	Parameter ID	Description
1	TrendO2Index	Oxygen Index
2	TrendAutoPEEP	Airway Auto PEEP
3	TrendAutoPEEPdelta	Delta Airway Auto PEEP
4	TrendAutoPEEPesoph	Esophageal Auto PEEP
5	TrendC20	Compliance Ratio
6	TrendCcw	Chestwall Compliance
7	TrendCdyn	Dynamic System Compliance
8	TrendCdynNorm	Dynamic System Compliance, Normalized
9	TrendClung	Lung Compliance
10	TrendCPAPFlow	Continuous Positive Airway Pressure



	Parameter ID	Description
		Flow
11	TrendCstat	Static System Compliance
12	TrendCstatNorm	Static System Compliance, Normalized
13	TrendFiO2	FiO2 – Fractional Inspiratory Oxygen Concentration
14	TrendIE	Inspiratory:Expiratory Time Ratio
15	TrendLeak	$(V_{ti}-V_{te})/V_{ti} \times 100$
16	TrendIveTotalNorm	Minute Volume, Normalized
17	TrendMIP	Maximum Inspiratory Pressure
18	TrendNcpapMeanFlow	Nasal Continuous Positive Airway Pressure Mean Flow.
19	TrendNcpapPres	Nasal Continuous Positive Airway Pressure
20	TrendP100	Respiratory Drive
21	TrendPaCO2	Arterial CO2 set by the operator.
22	TrendPair	Air Supply Pressure
23	TrendPaO2	Arterial Oxygen set by the operator.
24	TrendPawDelta	Delta Airway Pressure
25	TrendPeep	Positive End Expiratory Pressure
26	TrendPefr	Peak Expiratory Flow Rate
27	TrendPesDelta	Delta Esophageal Pressure
28	TrendPFRatio	PaO2/FiO2 Ratio
29	TrendPifr	Peak Inspiratory Flow Rate
30	TrendPmean	Mean Airway Pressure
31	TrendPO2	Oxygen Supply Pressure
32	TrendPpeak	Peak Airway Pressure
33	TrendPplat	Plateau Airway Pressure
34	TrendPplatPtp	Transpulmonary Pressure at Airway Pressure Plateau
35	TrendPtpPEEP	Transpulmonary Pressure at Auto PEEP
36	TrendRate	Total Breath Rate
37	TrendRateMand	Mandatory Breath Rate
38	TrendRateSpon	Spontaneous Breath Rate
39	TrendRimp	Imposed Resistance
40	TrendRlung	Lung Resistance
41	TrendRpeak	Peak Expiratory Airway Resistance
42	TrendRrs	Respiratory System Resistance
43	TrendRSBIndex	Rapid Shallow Breathing Index
44	TrendTe	Expiratory Time
45	TrendTi	Inspiratory Time
46	TrendVdel	Machine Delivered Volume
47	TrendVentilationAlveolar	Alveolar Ventilation

	Parameter ID	Description
48	TrendVentilationAlveolarDeadSpace	Alveolar Dead Space
49	TrendVentilationPhysiologicalDeadSpace	Physiological Dead Space
50	TrendVentilationPhysiologicalDeadSpaceVtRatio	Physiological Dead Space / Tidal Volume Ratio
51	TrendVeSpon	Spontaneous Minute Volume
52	TrendVeSponNorm	Spontaneous Minute Volume, Normalized
53	TrendVeTotal	Minute Volume
54	TrendVte	Exhaled Tidal Volume
55	TrendVteMand	Mandatory Tidal Volume
56	TrendVteMandNorm	Mandatory Tidal Volume, Normalized
57	TrendVteNorm	Exhaled Tidal Volume, Normalized
58	TrendVteSpon	Spontaneous Tidal Volume
59	TrendVteSponNorm	Spontaneous Tidal Volume, Normalize
60	TrendVti	Inspired Tidal Volume
61	TrendVtiNorm	Inspired Tidal Volume, Normalized
62	TrendWobImposed	Work of Breathing, Imposed
63	TrendWobPatient	Work of Breathing, Patient
64	TrendWobVent	Work of Breathing, Ventilator

## 7.7 Alarms

### 7.7.1 Active Alarm Indicators

	String ID	Description
1	AlarmActive	Indicates whether any alarms are currently being annunciated
2	AlarmActivePriority	Indicates the highest priority of all active alarm conditions
3	AlarmApnea	Apnea Interval Exceeded
4	AlarmCO2CheckAirwayAdapter	CO2 device reported a CO2 Airway Adapter needs to be checked.
5	AlarmCheckEvents	Check event notifications on the ventilator
6	AlarmCircDisc	Patient Circuit Disconnect
7	AlarmClockBattLow	Clock Battery Low
8	AlarmCO2CommunicationError	CO2 device reported a communication error.
9	AlarmCO2OutOfRange	CO2 device reported a CO2 range error condition.
10	AlarmCO2SensorFault	CO2 device reported a fault condition with CO2 sensor.
11	AlarmCO2SensorOverTemp	CO2 device reported a fault condition due to temperature.
12	AlarmCO2ZeroRequired	CO2 device reported sensor requires to be initialized to zero.
13	AlarmDefaults	Ventilator set to factory defaults
14	AlarmDirtyFilter	Check for a dirty filter
15	AlarmEEPROMFault	EEPROM failure detected
16	AlarmFanFail	Fan failure

	String ID	Description
17	AlarmFiO2AutoLimitHigh	The High Auto FiO <sub>2</sub> Alarm” will alarm if Auto FiO <sub>2</sub> Cmd $\geq$ preset High Auto FiO <sub>2</sub> Limit for a period of 60 seconds or greater.
18	AlarmFiO2AutoLimitLow	The output Auto FiO <sub>2</sub> Cmd shall not be allowed to be less than the Low Auto FiO <sub>2</sub> Limit. The “Low Auto FiO <sub>2</sub> Limit” will alarm if Auto FiO <sub>2</sub> Cmd $\leq$ preset Low Auto FiO <sub>2</sub> Limit for a period of 60 seconds or greater.
19	AlarmFiO2BaselineLimitHigh	Alarm if the Baseline FiO <sub>2</sub> $\geq$ High Baseline FiO <sub>2</sub> Alarm setting.
20	AlarmFiO2Cal	FiO <sub>2</sub> Sensor needs calibration
21	AlarmFiO2High	High FiO <sub>2</sub>
22	AlarmFiO2Low	Low FiO <sub>2</sub>
23	AlarmFiO2Range	FiO <sub>2</sub> out of upper & lower limits
24	AlarmFiO2SensorFail	FiO <sub>2</sub> Sensor failure
25	AlarmFlowSensorDisc	Flow sensor disconnect
26	AlarmHwFault	A hardware failure has been detected
27	AlarmIlvSlaveDisc	ILV slave disconnect detected
28	AlarmInop	Vent Inoperable
29	AlarmEndTidalCO2High	Exceeded EtCO <sub>2</sub> High Limit
30	AlarmEndTidalCO2Invalid	Invalid EtCO <sub>2</sub> .
31	AlarmEndTidalCO2Low	Below EtCO <sub>2</sub> Low Limit
32	AlarmInvalidGasId	Invalid Gas ID
33	AlarmInvalidSN	Invalid Serial Number
34	AlarmLimitIE	I:E Ratio Limit
35	AlarmLimitTi	Inspiratory Time Limit
36	AlarmLimitVol	Volume Limit
37	AlarmLossAir	Loss of Air Supply
38	AlarmLossGas	Loss of All Gas Supply
39	AlarmLossHeliox	Loss of Heliox Supply
40	AlarmLossO2	Loss of O <sub>2</sub> Supply
41	AlarmMotorFault	The turbine motor has failed
42	AlarmNcpapHigh	Nasal CPAP high pressure
43	AlarmNcpapHighPresLimit	Nasal CPAP high pressure (time) limit
44	AlarmNcpapLow	Nasal CPAP low pressure
45	AlarmNoCalData	No calibration data is available
46	AlarmO2PressHigh	High O <sub>2</sub> Inlet Pressure
47	AlarmOcclusion	Patient Circuit Occlusion
48	AlarmOpenSV	Safety Valve Open
49	AlarmPatientDefaults	Configured with patient default settings
50	AlarmPeepHigh	High PEEP
51	AlarmPeepLow	Low PEEP
52	AlarmPpeakHigh	High Peak Airway Pressure
53	AlarmPpeakHighExt	Extended High Peak Airway Pressure
54	AlarmPpeakLow	Low Peak Airway Pressure
55	AlarmPulseRateHigh	Pulse rate is greater than the High Pulse Rate alarm setting.
56	AlarmPulseRateLow	Pulse rate is less than the Low Pulse Rate alarm setting.

	String ID	Description
57	AlarmPwrAcLoss	Loss of AC Power
58	AlarmPwrBattLow	Low Battery
59	AlarmPwrBattMed	Medium Battery
60	AlarmRateHigh	High Breath Rate
61	AlarmSilence	Alarm annunciation is temporarily suspended on the device
62	AlarmSpO2Failure	High priority alarm is issued when the MS-aa PCB reports a board failure (see Masimo CSD-1086 Rev C).
63	AlarmSpO2High	Monitored SpO2 is greater than the preset High SpO2 Alarm for more than the SpO2 Alarm Delay period.
64	AlarmSpO2Low	Monitored SpO2 is less than the preset Low SpO2 Alarm for more than the SpO2 Alarm Delay period.
65	AlarmSpO2NotConnected	High priority alarm is issued when the pulse oximeter is not connected to the ventilator (i.e. RS-232 not connected).
66	AlarmSpO2SensorDefective	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Defective.
67	AlarmSpO2SensorNotConnected	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Not Connected.
68	AlarmSpO2SensorUnrecognized	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Unrecognized.
69	AlarmSpO2SignalAmbientLight	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Ambient Light.
70	AlarmSpO2SignalAmbientLightExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Ambient Light for an extended period.
71	AlarmSpO2SignalInterference	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Interference.
72	AlarmSpO2SignalInterferenceExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Interference for an extended period.
73	AlarmSpO2SignalLowSIQ	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal is Low.
74	AlarmSpO2SignalLowSIQExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal is Low for an extended period.
75	AlarmSpO2SignalLowPerfusion	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Perfusion is Low.

	String ID	Description
76	AlarmSpO2SignalLowPerfusionExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Perfusion is Low for an extended period.
77	AlarmSpO2SignalPulseSearch	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Pulse Search is Low.
78	AlarmSpO2SignalPulseSearchExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Pulse Search is Low for an extended period.
79	AlarmSpO2SignalSensorOff	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Sensor is off.
80	AlarmSpO2SignalSensorOffExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Sensor is off for an extended period.
81	AlarmTest	Audible Alarm Test, High Priority
82	AlarmTransducerFault	Transducer failure detected
83	AlarmVeHigh	High Minute Volume
84	AlarmVeLow	Low Minute Volume
85	AlarmVteLow	Low Tidal Volume
86	AlarmVtHigh	High Tidal Volume

### 7.7.2 History Alarm Indicators

The following table of nomenclature terms represent alarm conditions that are no longer being actively asserted but still require operator confirmation.

	String ID	Description
1	AlarmHistApnea	Apnea Interval Exceeded (historical / not cleared)
2	AlarmHistCO2CheckAirwayAdapter	CO2 device reported a CO2 Airway Adapter needs to be checked (historical / not cleared).
3	AlarmHistCheckEvents	Check event notifications on the ventilator (historical / not cleared)
4	AlarmHistCircDisc	Patient Circuit Disconnect (historical / not cleared)
5	AlarmHistClockBattLow	Clock Battery Low (historical / not cleared)
6	AlarmHistCO2CommunicationError	CO2 device reported a communication error (historical / not cleared).
7	AlarmHistCO2OutOfRange	CO2 device reported a CO2 range error condition (historical / not cleared).
8	AlarmHistCO2SensorFault	CO2 device reported a fault condition with CO2 sensor (historical / not cleared).
9	AlarmHistCO2SensorOverTemp	CO2 device reported a fault condition due to temperature (historical / not cleared).
10	AlarmHistCO2ZeroRequired	CO2 device reported sensor requires to be initialized to zero (historical / not cleared).
11	AlarmHistDefaults	Ventilator set to factory defaults (historical / not cleared)

	String ID	Description
12	AlarmHistDirtyFilter	Check for a dirty filter (historical / not cleared)
13	AlarmHistEEPROMFault	EEPROM failure detected (historical / not cleared)
14	AlarmHistEndTidalCO2High	Exceeded EtCO2 High Limit (historical / not cleared)
15	AlarmHistEndTidalCO2Invalid	Invalid EtCO2 (historical / not cleared).
16	AlarmHistEndTidalCO2Low	Below EtCO2 Low Limit (historical / not cleared).
17	AlarmHistFanFail	Fan failure (historical / not cleared)
18	AlarmHistFiO2AutoLimitHigh	Auto FiO2 Cmd is greater than or equal to the preset High Auto FiO2 Limit for a period of 60 seconds or more (historical / not cleared).
19	AlarmHistFiO2AutoLimitLow	Auto FiO2 Cmd is less than or equal to the preset Low Auto FiO2 Limit for a period of 60 seconds or more (historical / not cleared).
20	AlarmHistFiO2BaselineLimitHigh	Baseline FiO2 is greater than or equal to the High Baseline FiO2 Alarm setting sensor (historical / not cleared).
21	AlarmHistFiO2Cal	FiO2 Sensor needs calibration (historical / not cleared)
22	AlarmHistFiO2High	High FiO2 (historical / not cleared)
23	AlarmHistFiO2Low	Low FiO2 (historical / not cleared)
24	AlarmHistFiO2Range	FiO2 out of upper & lower limits (historical / not cleared)
25	AlarmHistFiO2SensorFail	FiO2 Sensor failure (historical / not cleared)
26	AlarmHistFlowSensorDisc	Flow sensor disconnect (historical / not cleared)
27	AlarmHistHwFault	A hardware failure has been detected (historical / not cleared)
28	AlarmHistIlvSlaveDisc	ILV slave disconnect detected (historical / not cleared)
29	AlarmHistInop	Vent Inoperable (historical / not cleared)
30	AlarmHistInvalidGasId	Invalid Gas ID (historical / not cleared)
31	AlarmHistInvalidSN	Invalid Serial Number (historical / not cleared)
32	AlarmHistLimitIE	I:E Ratio Limit (historical / not cleared)
33	AlarmHistLimitTi	Inspiratory Time Limit (historical / not cleared)
34	AlarmHistLimitVol	Volume Limit (historical / not cleared)
35	AlarmHistLossAir	Loss of Air Supply (historical / not cleared)
36	AlarmHistLossGas	Loss of All Gas Supply (historical / not cleared)
37	AlarmHistLossHeliox	Loss of Heliox Supply (historical / not cleared)
38	AlarmHistLossO2	Loss of O2 Supply (historical / not cleared)
39	AlarmHistMotorFault	The turbine motor has failed (historical / not cleared)
40	AlarmHistNcpapHigh	Nasal CPAP high pressure
41	AlarmHistNcpapHighPresLimit	Nasal CPAP high pressure (time) limit
42	AlarmHistNcpapLow	Nasal CPAP low pressure
43	AlarmHistNoCalData	No calibration data is available (historical / not cleared)
44	AlarmHistO2PressHigh	High O2 Inlet Pressure (historical / not cleared)
45	AlarmHistOcclusion	Patient Circuit Occlusion (historical / not cleared)
46	AlarmHistOpenSV	Safety Valve Open (historical / not cleared)
47	AlarmHistPatientDefaults	Configured with patient default settings (historical / not cleared)
48	AlarmHistPeepHigh	High PEEP (historical / not cleared)
49	AlarmHistPeepLow	Low PEEP (historical / not cleared)

	String ID	Description
50	AlarmHistPpeakHigh	High Peak Airway Pressure (historical / not cleared)
51	AlarmHistPpeakHighExt	Extended High Peak Airway Pressure (historical / not cleared)
52	AlarmHistPpeakLow	Low Peak Airway Pressure (historical / not cleared)
53	AlarmHistPulseRateHigh	Pulse rate is greater than the High Pulse Rate alarm setting (historical / not cleared).
54	AlarmHistPulseRateLow	Pulse rate is less than the Low Pulse Rate alarm setting (historical / not cleared).
55	AlarmHistPwrAcLoss	Loss of AC Power (historical / not cleared)
56	AlarmHistPwrBattLow	Low Battery (historical / not cleared)
57	AlarmHistPwrBattMed	Medium Battery (historical / not cleared)
58	AlarmHistRateHigh	High Breath Rate (historical / not cleared)
59	AlarmHistSpO2Failure	High priority alarm is issued when the MS-aa PCB reports a board failure (see Masimo CSD-1086 Rev C) - (historical / not cleared).
60	AlarmHistSpO2High	Monitored SpO2 is greater than the preset High SpO2 Alarm for more than the SpO2 Alarm Delay period (historical / not cleared).
61	AlarmHistSpO2Low	Monitored SpO2 is less than the preset Low SpO2 Alarm for more than the SpO2 Alarm Delay period (historical / not cleared).
62	AlarmHistSpO2NotConnected	High priority alarm is issued when the pulse oximeter is not connected to the ventilator (i.e. RS-232 not connected) - (historical / not cleared).
63	AlarmHistSpO2SensorDefective	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Defective (historical / not cleared).
64	AlarmHistSpO2SensorNotConnected	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Not Connected (historical / not cleared).
65	AlarmHistSpO2SensorUnrecognized	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Unrecognized (historical / not cleared).
66	AlarmHistSpO2SignalAmbientLight	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Ambient Light (historical / not cleared).
67	AlarmHistSpO2SignalAmbientLightExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Ambient Light for an extended period (historical / not cleared).
68	AlarmHistSpO2SignalInterference	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Interference (historical / not cleared).
69	AlarmHistSpO2SignalInterferenceExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2

	String ID	Description
		Signal Interference for an extended period (historical / not cleared).
70	AlarmHistSpO2SignalLowPerfusion	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Perfusion is Low (historical / not cleared).
71	AlarmHistSpO2SignalLowPerfusionExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Perfusion is Low for an extended period (historical / not cleared).
72	AlarmHistSpO2SignalLowSIQ	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal is Low (historical / not cleared).
73	AlarmHistSpO2SignalLowSIQExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal is Low for an extended period (historical / not cleared).
74	AlarmHistSpO2SignalPulseSearch	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Pulse Search is Low (historical / not cleared).
75	AlarmHistSpO2SignalPulseSearchExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Pulse Search is Low for an extended period (historical / not cleared).
76	AlarmHistSpO2SignalSensorOff	Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Sensor is off (historical / not cleared).
77	AlarmHistSpO2SignalSensorOffExt	High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Sensor is off for an extended period (historical / not cleared).
78	AlarmHistTest	Audible Alarm Test, High Priority (historical / not cleared)
79	AlarmHistTransducerFault	Transducer failure detected (historical / not cleared)
80	AlarmHistVeHigh	High Minute Volume (historical / not cleared)
81	AlarmHistVeLow	Low Minute Volume (historical / not cleared)
82	AlarmHistVteLow	Low Tidal Volume (historical / not cleared)
83	AlarmHistVtHigh	High Tidal Volume (historical / not cleared)

## 7.8 Scalars

	Parameter ID	Description
1	WaveAnlg0	Analog Input Channel 0
2	WaveAnlg1	Analog Input Channel 1
3	WaveFexp	Expiratory Sensor Flow
4	WaveFinsp	Inspiratory Sensor Flow
5	WaveFlow	Airway Flow



	Parameter ID	Description
6	WaveMetric	Phase, Trigger, Type, etc.
7	WavePaw	Airway Pressure
8	WavePCO2	Measured CO2 level
9	WavePes	Esophageal Pressure
10	WavePinsp	Inspiratory Sensor Pressure
11	WavePtp	Transpulmonary Pressure
12	WavePtr	Tracheal Pressure
13	WaveSpO2Pleth	Raw IR signal over time from the MS-11 PCB. The range value is inverted and scaled from the raw data sent by the oximeter (oximeter +127 = 0%, oximeter – 128 = 100%).
14	WaveVt	Tidal Volume

## 7.9 System Information

	Parameter ID	Description
1	SysInfoConfig	Model Configuration
2	SysInfoDateTime	Current Date & Time Setting
3	SysInfoModel	Base Model
4	SysInfoOUI	Organizationally Unique Identifier
5	SysInfoSerial	Serial Number
6	SysInfoSwVer	Software Version
7	SysInfoSwVerBootLoader	Boot Loader Software Version Number
8	SysInfoTimeSession	Hours operated since Power ON
9	SysInfoTimeTotal	System Lifetime Hours Operated
10	SysInfoTurbineSerial	Turbine Serial Number

## 8 Constant Definitions

The following constant values are defined for this protocol:

VOXP Constant Definitions		
Constant	Value	Description
Parameter Not Available	+(MAX_NEG_INT)	Indicates that the associated parameter is either not currently available (e.g., due to the operational mode of the device) or is out-of-range and thus a valid value cannot be provided. For 32-bit integers, the value is 0x80000000; for 16-bit integers, 0x8000; and for 8-bit integers, 0x80
Parameter Off	+(MAX_NEG_INT + 1)	Indicates that a signed integer parameter is Off; typically used for settings that may be turned Off as opposed to On but with a zero value. For 32-bit integers, the value is 0x80000001; for 16-bit integers, 0x8001; and for 8-bit integers, 0x81

## 9 Appendix A – CRC Definition and Algorithm

To calculate the 16 bit CRC the message bits are considered to be the coefficients of a polynomial. This message polynomial is first multiplied by  $X^{16}$  and then divided by the generator polynomial ( $X^{16} + X^{12} + X^5 + 1$ ) using modulo two arithmetic. The remainder left after the division is the desired CRC. Since a message block in the Modem Protocol is 128 bytes or 1024 bits, the message polynomial will be of order  $X^{1023}$ . The hi order bit of the first byte of the message block is the coefficient of  $X^{1023}$  in the message polynomial. The lo order bit of the last byte of the message block is the coefficient of  $X^0$  in the message polynomial.

### 16-bit CRC Algorithm in C++:

```

unsigned short
usCalcCRC(unsigned char *pBuffer, unsigned short usLength) {
 unsigned short usCRC = 0;

 // Note: CRC-CCITT Polynomial
 const unsigned short usGeneratorPolynom = 0x1021;

 int iCarry;

 for(; usLength > 0; usLength--, pBuffer++)
 {
 usCRC ^= (unsigned short)*pBuffer << 8;
 for(int i = 0; i < 8; i++)
 {
 if(iCarry = ((usCRC & 0x8000) != 0))
 {
 usCRC ^= usGeneratorPolynom;
 }
 usCRC = (usCRC << 1) + iCarry;
 }
 }
 return usCRC ;
}

```

## 10 VIASYS VOXP Compliant Ventilators and Supported Data Set

As described in earlier sections, the VOX Protocol is a flexible and extensible communication protocol that supports any VIASYS ventilator conforming to the protocol. The external host system shall conform and be able to handle the specific communication supported on the various ventilators, including the ventilator supported data set using the common nomenclature.

### 10.1 Ventilator Configuration / Setup

This section describes the system level interface, controls, and other items required for an external host system using the VOXP communication protocol to communicate with the AVEA or VELA ventilator.

#### 10.1.1 System Configuration Requirements

Communication is established between a ventilator and an external system with a properly configured set of system level items from physical cables, adapters, communication parameters to application protocols. This section will describe the interface communication layers and the tasks that must be performed in order to establish proper communication to a ventilator from an external host system.

#### 10.1.2 Physical and Link Layer

The AVEA and VELA VOXP interface shall use the ventilator's external communication (or "MIB") port to communicate to an external host system. The physical interface is an RJ-45 connector with RS-232 signaling levels. Only the following pins are required:

External Host to AVEA MIB Cable			
External Host Label	Direction	MIB	
		RJ45 Pin	Label
RxD	⇐	5	dTxD
TxD	⇒	7	dRxD
Digital GND	↔	4	dGND

In order to make the connection to an external host system and AVEA or VELA ventilator, a 10BaseT STP cable, and RJ45 to DB-9 or DB-25 adapter shall be used to convert from the external host's serial cable to the standard connection used by the ventilator<sup>3</sup>. There is no configuration required of the external host system other than the proper cabling and serial data communication parameters:

Serial data communication parameters shall be configurable to the following:

BAUD Rate: 115200, 57600, 38400 (default), 19200, 9600  
 Data Bits: 8 (default) or 7  
 Stop Bits: 1  
 Parity: None (default), Odd, Even

<sup>3</sup> Defined in ISO/IEEE 11073-30200-2004.

**Note: the maximum number of waveforms that can be selected is 4. It is recommended to always select Wave Metric when another waveform is selected. When a waveform is selected, the baud rate should be 57600 or higher.**

### 10.1.3 Ventilator Configuration – Setup / Communication Protocol

This section defines several groups of identifiers in the nomenclature that are used in the VOXP protocol. However, the implementation specifics relating to any identifier in this nomenclature can be found in the target ventilator specification complying with the VOXP protocol. These ventilator specification details include, but are not limited to, data type and domain values.

### 10.1.4 Logical Connections, Link State, and Link Restarting

Although the protocol is not “technically” connection-oriented, it simulates a “connected” state to ensure the protocol setup information is not modified while there is an active connection with an external host.

The ventilator will indicate that a connection is active when it receives a valid reply message (i.e. link.ack()) from the external serial port. After the first 2-minutes of sending a link message (i.e. link.ping()) every 5-seconds, and not receiving a response from the host with a link.ack() message, the connection will be indicated as broken, enabling the communication settings or protocol selection to be modified.

The ventilator shall transmit a link.restart() message to the connected host system upon an acceptance of a change of any one of the following supported parameters:

- a) Language (English, Chinese, Spanish, etc.)
- b) Communication Parameter Change (i.e. Baud rate, data bit, parity, etc.)
- c) Patient Size (Neonate, Pediatric, or Adult)

## 10.2 VOXP Message, Syntax, and Sequence Description

Each VOXP compliant ventilator shall conform to the message format, syntax, and exchange sequence specified in the previous sections of this document. The protocol provides the flexibility to define the values for the “profile” message’s voxpVersion and profileVersion attributes, thus enabling the host system to support current and future versions of the VOXP protocol (e.g. new commands, data-sets, message types, etc.) on any VOXP compliant ventilator.

### 10.2.1 Epoch Data

The AVEA and VELA epoch wave data shall be updated every 500msec.

### 10.2.2 CRC Attribute

All data messages shall contain a cyclic redundancy check (CRC) on the data portion of the message.

For example:

```
<data class="monitor" crc="3A99"
msgID="0005">80008000800000067800005C580008000800080000015FF9D000A8000800
080000029001E000604218000023F0008000000248000800080000008000800008000800
000608000800002A9004502FCC25500000000013F02A43E5D0260170B8000800000000
000000002F52A7C80008000800000098</data>
```

### 10.2.3 AVEA Ventilator Specific Message, Syntax and Sequence

The table below defines the values for the “profile” message’s voxpVersion and profileVersion attributes for an AVEA ventilator.

The Profile message’s “model” attribute shall be set to Avea. For example: <profile model=”Avea”>.

For the TEXT data type, unless otherwise specified (in the profile element’s textEncoding attribute), the character encoding shall be UTF-16 (16-bit characters).

The VOXP “profile” message includes the attributes “voxpVersion” and “profileVersion” (see section 5.1). The following table correlates values for these attributes, namely major and minor specialization profile version numbers, identified by part number (PN) and revision of the respective specification.

AVEA Major/Minor Version Number Assignments					
PN / Revision	AVEA Profile <sup>4</sup>		VOXP <sup>5</sup> Version		Rationale for Version
	Major	Minor	Major	Minor	
91316/A					
91316/B	1	0	2	1	Initial release
91316/C	2	0	3	0	Updated for nasal CPAP semantics; harmonized with latest version of the base VOXP specification; synchronized format and content with the VELA VOXP specialization document.
ER-1908/A	2	1	3	2	

### 10.2.4 VELA Ventilator Specific Message, Syntax and Sequence

The table below defines the values for the “profile” message’s voxpVersion and profileVersion attributes for an VELA ventilator.

VELA Major/Minor Version Number Assignments					
PN / Revision	VELA Profile <sup>6</sup>		VOXP <sup>7</sup>		Rationale for Version
	Major	Minor	Major	Minor	
91415/A	1	0	3	1	Initial release
ER-1908/A	1	1	3	2	

<sup>4</sup> For example, the attribute in the profile message could be: profileVersion=”2.0”

<sup>5</sup> For example, the attribute in the profile message could be: voxpVersion=”3.2”

<sup>6</sup> For example, the attribute in the profile message would be: profileVersion=”1.0”

<sup>7</sup> For example, the attribute in the profile message would be: voxpVersion=”3.0”

The Profile message's "model" attribute shall be set to 'Vela (space) Model'.

For example: <profile model="Vela Comprehensive">  
 <profile model="Vela Plus">  
 <profile model="Vela Plus International">  
 <profile model="Vela Basic">

For the TEXT data type, unless otherwise specified (in the profile element's textEncoding attribute), the character encoding shall be UTF-32 (32-bit characters).

### 10.2.5 AVEA Ventilator Data type, and Bit Map Representation

The following identifiers shall conform to the respective data type, bit map definitions, and algorithms:

a) **SysInfoConfig**: A UWORD type that describes the options on the specified model.

Bit	Type	Identifier	Description
0	unsigned int	bCompressor	0=compressor option not installed; 1=compressor option installed.
1	unsigned int	bHeliOx	0=heliox option not installed; 1=heliox option installed.
2	unsigned int	bHeliOxNeb	0=heliox and nebulizer option installed; 1= heliox and nebulizer option not installed.
3	unsigned int	bPflex	0=pressure inflexion option not installed; 1=pressure inflexion option installed.
4	unsigned int	bNeo	0=neonate functions not available; 1=neonate functions available.
5	unsigned int	bPed	0=pediatric functions not available; 1=pediatric functions available.
6	unsigned int	bAdult	0=adult function not available; 1=adult functions available.
7	unsigned int	bManeuvers	0=maneuver options not available; 1=maneuver options available.
8	unsigned int	bNcpap	0=Nasal CPAP not available; 1=Nasal CPAP available
9-15	unused		Not used.

b) **WaveMetric**: A UWORD type with structure described below:

Bit	Type	Identifier	Description
0-2	enum	VpBreathPhase	0=NOS*_PHASE; 1=INSP_PHASE; 2=EXP_PHASE; 3=INSP_PAUSE_PHASE; 4=EXP_PAUSE_PHASE.
3-4	enum	VpBreathType	0=NOS*_BREATH TYPE; 1=SPONTANEOUS_BREATH_TYPE; 2=MANDATORY_BREATH_TYPE; 3=SPECIAL_BREATH TYPE.

Bit	Type	Identifier	Description
5-6	enum	VpBreathTrigger	0=PATIENT_TRIGGERED; 1=VENTILATOR_TRIGGER.
7	bool	VpSpO2SignalIQPresent	0=Signal IQ Not Included; 1=Signal IQ Included
8-9	unused		Not used.
10-15	unsigned int	VpSpO2SignalIQ	SpO2 Signal IQ Value (6 bits)

\*NOS = “Not otherwise specified”

### 10.2.6 VELA Ventilator Data type, and Bit Map Representation

The following identifiers shall conform to the respective data type, bit map definitions, and algorithms:

a) **SysInfoConfig**: A UWORD type that describes the options on the specified model.

Bit	Type	Identifier	Description
15	unsigned int	bLeakCompensation	0=Leak Compensation option not available; 1=Leak Compensation available.
14	unsigned int	bMIPNIF <sup>8</sup>	0=MIP/NIF option not available; 1= MIP/NIF option available.
13	unsigned int	bNPPV	1=NPPV option available; 0= NPPV option not available.
12	unsigned int	bAPRV	0=APRV option not available; 1=APRV option available.
11	unsigned int	bPRVC	0=PRVC not available; 1=PRVC available.
10	unsigned int	bVSYNC	0=VSync not available; 1= VSync available.
9	unsigned int	bSquareWaveform	0=Square Waveform not available; 1=Square Waveform available.
8	unsigned int	bAssuredVolumeVAPS	0=Assured Volume VAPS not available; 1=Assured Volume VAPS available.
7	unsigned int	bLoops	0=“Loop” breath displays are not available; 1=“Loop” breath displays are available
6	unsigned int	bTrends	0=Trended monitors are not available; 1=Trended monitors are available
5-0	unused		Not used.

b) **WaveMetric**: A UWORD type with structure described below:

Bit	Type	Identifier	Description
0-2	enum	VpBreathPhase	0=NOS* _PHASE; 1=INSP _PHASE; 2=EXP _PHASE; 3=INSP _PAUSE _PHASE; 4=EXP _PAUSE _PHASE.

<sup>8</sup> MIP/NIF = Maximum Inspiratory Pressure / Negative Inspiratory Force; related to a VELA Maneuver.

3-4	enum	VpBreathType	0=NOS*_BREATH_TYPE; 1=SPONTANEOUS_BREATH_TYPE; 2=MANDATORY_BREATH_TYPE.
5-6	enum	VpBreathTrigger	0=PATIENT_TRIGGERED; 1=VENTILATOR_TRIGGERED.
7	bool	VpSpO2SignalIQPresent	0=Signal IQ Not Included; 1=Signal IQ Included
8-9	unused		Not used.
10-15	unsigned int	VpSpO2SignalIQ	SpO2 Signal IQ Value (6 bits)

\*NOS = “Not otherwise specified”

### 10.3 VOXP Data Support

The AVEA and VELA graphical user interface (GUI) nomenclature shall support and be mapped to the VOXP nomenclature as presented in the class descriptions below. The ventilators will support 5-classes of data:

- a. Info These parameters represent general information that are typically unique to the system or may relate to the patient (e.g., patient identifier).
- b. Setting These parameters represent the operational control values that the caregiver has accepted for the ventilator.
- c. Monitor These parameters represent information that the ventilator/system is monitoring, and will be dynamic based on the specific parameter.
- d. Alarm These parameters indicate operational conditions of the ventilator that may require operator review and corrective action.
- e. Scalar These parameters are updated every 500ms and are typically represented graphically on a two-dimensional scale.

VOXP messages are queued for transmission according to a prioritization of the classes. The classes are prioritized as follows:

- 1) Alarms,
- 2) Scalars / Waves
- 3) Monitors
- 4) Settings
- 5) Info

Where the Alarm class has the highest priority, and the Info class has the lowest.

In the following sections, ‘ID’ is the nomenclature from the various tables in section 7 with the label of Parameter ID or String ID. ‘VELA GUI / Membrane’ or ‘AVEA GUI / Membrane’ refer to the nomenclature displayed on the respective ventilator touch screen or display membrane, thus providing a map to the corresponding VOXP nomenclature for the same item. ‘Description’ provides a brief definition for each specific item.



## 10.4 Info Class

### 10.4.1 AVEA Ventilator – Info Class

<b>10.4.1.1</b>	<b>ID:</b>	<b>SysInfoConfig</b>
	Description:	Describes the options that are installed/available on the AVEA ventilator (see section on Data type and Bit Map Representation for details.).
	Type:	UWORD
	AVEA GUI/Membrane:	none
<b>10.4.1.2</b>	<b>ID:</b>	<b>SysInfoModel</b>
	Description:	The base model configuration
	Label:	“Base Model”
	Type:	ENUM
	Enum value=label:	0=“Invalid”; 1=“Avea Comp”; 2=“Avea +”; 3= “Avea”
	AVEA GUI/Membrane:	none
<b>10.4.1.3</b>	<b>ID:</b>	<b>SysInfoOUI</b>
	Description:	The Organizationally Unique Identifier ('OUI') is a 64-bit value that uniquely identifies each VIASYS Respiratory Care device. The upper 24 bits is a unique value assigned to VIASYS by the IEEE Registration Authority (00-07-A1 hexadecimal). The lower 40 bits are used as a serialization number by VIASYS. Given that the AVEA OUI would consist of the following ASCII HEX bytes: 0007A100B9E92A00 the conversion to a 16-bit character would look like: 003000300030003700410031003000300042003900450 03900320041003000300000
	Type:	TEXT
	AVEA GUI/Membrane:	none
<b>10.4.1.4</b>	<b>ID:</b>	<b>SysInfoTimeTotal</b>
	Description:	The cumulative amount of time in 100 <sup>th</sup> of an hour (i.e. a value of 1=36 seconds) since the ventilator has been in operation.
	Type:	UINT
	Scale:	2
	Range (Adult/Ped/Neo):	0 – 9,999,999
	AVEA GUI/Membrane:	

- 10.4.1.5 ID: SysInfoSerial**  
 Description: AVEA serial number – consists of 3-alphabet characters and 5-digits (e.g. AEV01001).  
 Type: TEXT  
 AVEA GUI/Membrane: see back panel
- 10.4.1.6 ID: SysInfoSwVer**  
 Description: SysInfoSwVer describes the AVEA's software version release that is executing on the ventilator.  
 Type: TEXT  
 AVEA GUI/Membrane: On-power-up
- 10.4.1.7 ID: PatInfoID**  
 Description: This is an identifier corresponding to the patient identification.  
 Type: TEXT  
 Label: IDENTIFICATION  
 AVEA GUI/Membrane: Patient Identification

#### 10.4.2 VELA Ventilator – Info Class

- 10.4.2.1 ID: PatInfoID**  
 Description: This is an identifier corresponding to the patient identification.  
 Label: IDENTIFICATION  
 Type: TEXT  
 VELA GUI/Membrane: Patient Identification
- 10.4.2.2 ID: SysInfoConfig**  
 Description: This describes the options that are installed/available on the VELA ventilator (see Data type, and Bit Map Representation section).  
 Type: UWORD  
 VELA GUI/Membrane: none
- 10.4.2.3 ID: SysInfoDateTime**  
 Description: The ventilator's current date & time setting (string with a format of: YYYY-MM-DDTHH:MM:SS.SSS)  
 Type: TEXT  
 VELA GUI/Membrane: Date/Time Screen
- 10.4.2.4 ID: SysInfoModel**  
 Description: The base model configuration.  
 Label: "VENT MODEL"  
 Type: ENUM

Enum value=label: 0="Vela Comprehensive"; 1="Vela Plus"; 2="Vela Plus International"; 3="Vela Basic".

VELA GUI/Membrane: Model

#### 10.4.2.5 ID: **SysInfoOUI**

Description: The Organizationally Unique Identifier ('OUI') is a 64-bit value that uniquely identifies each VIASYS Respiratory Care device. The upper 24 bits is a unique value assigned to VIASYS by the IEEE Registration Authority (00-07-A1 hexadecimal). The lower 40 bits are used as a serialization number by VIASYS. Given that the VELA OUI would consist of the following ASCII HEX bytes: 0007A100B9E92A00 the conversion to a 16-bit character would look like:  
003000300030003700410031003000300042003900450  
03900320041003000300000.

Type: TEXT

VELA GUI/Membrane: EUI<sup>9</sup>

#### 10.4.2.6 ID: **SysInfoSerial**

Description: VELA serial number – consists of 3-alphabet characters and 5-digits (e.g. AET01230).

Type: TEXT

VELA GUI/Membrane: VENT S/N

#### 10.4.2.7 ID: **SysInfoSwVer**

Description: SysInfoSwVer describes the VELA's software version release that is executing on the ventilator.

Type: TEXT

VELA GUI/Membrane: MAIN S/W VERSION

#### 10.4.2.8 ID: **SysInfoSwVerBootLoader**

Description: VELA's boot loader software version release identifier.

Type: TEXT

VELA GUI/Membrane: BOOTLOADER VERSION

#### 10.4.2.9 ID: **SysInfoTimeTotal**

Description: The cumulative amount of time in 100<sup>th</sup> of an hour (i.e. a value of 1=36 seconds) since the ventilator has been in operation.

Type: UINT

Scale: 2

Range (Adult/Ped/Neo): 0 – 9,999,999

Label: VENT HOURS

<sup>9</sup> Select Extended Functions screen and then Version Info screen.

VELA GUI/Membrane: VENT HOURS

**10.4.2.10 ID: SysInfoTurbineSerial**  
 Description: Serial number for the ventilator's turbine subsystem.  
 Type: TEXT  
 VELA GUI/Membrane: TURB S/N

## 10.5 Setting Class

### 10.5.1 AVEA Ventilator – Setting Class

**10.5.1.1 ID: SetCO2MonitoringEnable**  
 Description: Enable or Disable CO2 Monitoring  
 Type: ENUM  
 Enum value = label: 0=Disabled; 1=Enabled  
 Label: CO2  
 AVEA GUI/Membrane: CO2

**10.5.1.2 ID: SetEndTidalCO2Average**  
 Description: Number of breaths EtCO2 is averaged over.  
 Type: WORD  
 Resolution: 1  
 Range: 1 - 8  
 Units: breath  
 Label: EtCO2 Avg.  
 AVEA GUI/Membrane: EtCO2 Avg.

**10.5.1.3 ID: SetFiO2**  
 Description: Percent of oxygen that the ventilator is set to administer.  
 Type: WORD  
 Resolution: 1  
 Range: 21 – 100  
 Units: %  
 Label: FiO2  
 AVEA GUI/Membrane: FiO2

**10.5.1.4 ID: SetFlowBias**  
 Description: Flow delivered to the breathing circuit during the expiratory phase.  
 Type: WORD  
 Scale: 1  
 Resolution: 1  
 Range: 4 - 50  
 Units: L/min  
 Label: Bias Flow  
 AVEA GUI/Membrane: Bias Flow

<b>10.5.1.5</b>	<b>ID:</b>	<b>SetFlowCycle</b>
	Description:	Percent of the peak inspiratory flow (Peak Flow), at which the inspiratory phase of Pressure Controlled (including TCPL) breaths are terminated.
	Type:	WORD
	Resolution:	5
	Range:	0 - 45
	Units:	%
	Label:	Flow Cycle
	AVEA GUI/Membrane:	Flow Cycle
<b>10.5.1.6</b>	<b>ID:</b>	<b>SetFlowCyclePsv</b>
	Description:	Percent of peak inspiratory flow at which the inspiratory phase of a PSV breath is terminated
	Type:	WORD
	Resolution:	5
	Range (Adult/Ped/Neo):	5 - 45
	Units:	%
	Label:	PSV Cycle
	AVEA GUI/Membrane:	Psv Cycle
<b>10.5.1.7</b>	<b>ID:</b>	<b>SetFlowDemand</b>
	Description:	Enabled/Disabled state for the Intra-Breath Demand Flow feature for Volume Controlled breaths.
	Type:	WORD
	Resolution:	1
	Range (Adult/Ped/Neo):	0 - 1
	Units:	0=Off 1=On
	Label:	Demand Flow
	AVEA GUI/Membrane:	Demand Flow
<b>10.5.1.8</b>	<b>ID:</b>	<b>SetFlowInsp</b>
	Description:	Flow targeted by the ventilator during the inspiratory phase of flow controlled breaths.
	Type:	WORD
	Scale:	1
	Resolution (Adult/Ped):	10
	Resolution (Neonate):	2
	Range (Adult):	30 - 1500
	Range (Pediatric):	10 - 750
	Range (Neonate):	4 - 300
	Units:	L/min
	Label:	Peak Flow
	AVEA GUI/Membrane:	Peak Flow

- 10.5.1.9 ID: SetPauseInsp**  
 Description: Time Inspiration is extended before exhalation occurs after a volume breath is delivered.  
 Type: WORD  
 Scale: 2  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 300  
 Units: sec  
 Label: Insp Pause  
 AVEA GUI/Membrane: Insp Pause
- 10.5.1.10 ID: SetPresHigh**  
 Description: Baseline target for the “high” phase of APRV/BiPhasic mode  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 90  
 Units: cmH2O  
 Label: Pres High  
 AVEA GUI/Membrane: Pres High
- 10.5.1.11 ID: SetPresInsp**  
 Description: Pressure target for mandatory Pressure Controlled breaths.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped): 0 - 90  
 Range (Neonate): 0 - 80  
 Units: cmH2O  
 Label: Insp Pres  
 AVEA GUI/Membrane: Insp Pres
- 10.5.1.12 ID: SetPresLow**  
 Description: Baseline target for the “low” phase of APRV/BiPhasic mode.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 45  
 Units: cmH2O  
 Label: Pres Low  
 AVEA GUI/Membrane: Pres Low
- 10.5.1.13 ID: SetPresNasalCPAP**  
 Description: Nasal CPAP level. Applies to neonates only.  
 Type: WORD  
 Resolution: 1  
 Range: 2 - 10  
 Units: cmH2O

Label: nCPAP  
 AVEA GUI/Membrane: nCPAP

- 10.5.1.14 ID: SetPresPeep**  
 Description: Positive End Expiratory Pressure. I.e., target pressure to be maintained at the end of exhalation.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 50  
 Units: cmH2O  
 Label: PEEP  
 AVEA GUI/Membrane: PEEP
- 10.5.1.15 ID: SetPresPsv**  
 Description: Pressure target for spontaneous (PSV or Pressure Support Ventilation) breaths.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped): 0 – 90  
 Range (Neonate): 0 – 80  
 Units: cmH2O  
 Label: PSV  
 AVEA GUI/Membrane: PSV
- 10.5.1.16 ID: SetRate**  
 Description: Mandatory breath rate setting.  
 Type: WORD  
 Resolution: 1  
 Range (Adult): 1 - 120  
 Range (Ped/Neo): 1 - 150  
 Units: bpm  
 Label: Rate  
 AVEA GUI/Membrane: Rate
- 10.5.1.17 ID: SetRiseInsp**  
 Description: Setting controls the slope of the pressure rise during a Pressure Controlled breath. (Fastest Rise =1 and Slowest Rise=9).  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 1 - 9  
 Label: Insp Rise  
 AVEA GUI/Membrane: Insp Rise

- 10.5.1.18 ID: SetRisePsv**  
 Description: Setting controls the slope of the pressure rise during Pressure Supported (PSV) breaths. (Fastest Rise =1 and Slowest Rise=9).  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 1 - 9  
 Label: PSV Rise  
 AVEA GUI/Membrane: PSV Rise
- 10.5.1.19 ID: SetRiseVsync**  
 Description: Setting controls the slope of the pressure rise during Vsync breaths (Vsync is a Pressure Controlled, Volume Targeted modification of Volume Control) (Fastest Rise =1 and Slowest Rise=9).  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 1 - 9  
 Label: Vsync Rise  
 AVEA GUI/Membrane: Vsync Rise
- 10.5.1.20 ID: SetTimeHigh**  
 Description: Time duration setting for the “high” phase of APRV/BiPhasic mode.  
 Type: WORD  
 Scale: 1  
 Resolution: 1  
 Range (Adult/Ped/Neo): 2 - 300  
 Units: sec  
 Label: Time High  
 AVEA GUI/Membrane: Time High
- 10.5.1.21 ID: SetTimeHighPsv**  
 Description: Enabled/Disabled state of PSV breaths is available during Time High by activating T High PSV (an advanced setting of Time High).  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 1  
 Units: 0=Off 1=On  
 Label: T High PSV  
 AVEA GUI/Membrane: T High PSV
- 10.5.1.22 ID: SetTimeHighSync**  
 Description: Percentage of Time High in APRV/BiPhasic mode allowed at the end of the Time High period for



synchronizing the high to low baseline transition with a patient exhalation.

Type: WORD  
 Resolution: 5  
 Range (Adult/Ped/Neo): 0 - 50  
 Units: %  
 Label: T High Sync  
 AVEA GUI/Membrane: T High Sync

#### 10.5.1.23 ID: **SetTimeInsp**

Description: Time duration setting for the inspiratory phase of Pressure Controlled (including TCPL) mandatory breaths.

Type: WORD  
 Scale: 2  
 Resolution: 1  
 Range (Adult/Ped): 20 - 500  
 Range (Neonate): 15 -300  
 Units: sec  
 Label: Insp Time  
 AVEA GUI/Membrane: Insp Time

#### 10.5.1.24 ID: **SetTimeLow**

Description: Time duration setting for the “low” phase of APRV/BiPhasic mode.

Type: WORD  
 Scale: 1  
 Resolution: 1  
 Range (Adult/Ped/Neo): 2 - 300  
 Units: sec  
 Label: Time Low  
 AVEA GUI/Membrane: Time Low

#### 10.5.1.25 ID: **SetTimeLowSync**

Description: Percentage of Time Low in APRV/BiPhasic mode allowed at the end of the Time Low period for synchronizing the low to high baseline transition with a patient inspiration.

Type: WORD  
 Resolution: 5  
 Range (Adult/Ped/Neo): 0 - 50  
 Units: %  
 Label: T Low Sync  
 AVEA GUI/Membrane: T Low Sync

- 10.5.1.26 ID: SetTmaxPsv**  
 Description: Maximum time duration allowed for the inspiratory phase of a Pressure Supported breath.  
 Type: WORD  
 Scale: 2  
 Resolution: 1  
 Range (Adult/Ped): 20 - 500  
 Range (Neonate): 15 - 300  
 Units: sec  
 Label: PSV Tmax  
 AVEA GUI/Membrane: PSV Tmax
- 10.5.1.27 ID: SetTrigFlow**  
 Description: Flow Trigger level (i.e. Net Flow level at which inspiration is triggered).  
 Type: WORD  
 Scale: 1  
 Resolution: 1  
 Range (Adult/Ped/Neo): 1 - 200  
 Units: L/min  
 Label: Flow Trig  
 AVEA GUI/Membrane: Flow Trig
- 10.5.1.28 ID: SetTrigPres**  
 Description: Pressure Trigger level. I.e., Pressure drop below PEEP at which inspiration is triggered.  
 Type: WORD  
 Scale: 1  
 Resolution: 1  
 Range (Adult/Ped/Neo): 1 - 200  
 Units: cmH2O  
 Label: Pres Trig  
 AVEA GUI/Membrane: Pres Trig
- 10.5.1.29 ID: SetVCO2Average**  
 Description: Patient's exhaled minute volume of CO2.  
 Type: WORD  
 Resolution: 3  
 Range (Adult/Ped/Neo): 3 - 12  
 Units: Minute  
 Label: VCO2 Avg.  
 AVEA GUI/Membrane: VCO2 Avg.
- 10.5.1.30 ID: SetVol**  
 Description: Volume targeted by the ventilator for delivery of Volume Controlled breaths.

Type: WORD  
 Scale (Adult): 4  
 Scale (Pediatric/Neo): 1  
 Resolution (Adult): 100  
 Resolution (Pediatric): 10  
 Resolution (Neonate): 1  
 Range (Adult): 1000 – 25000  
 Range (Pediatric): 250 - 5000  
 Range (Neonate): 20 - 3000  
 Units (Adult): L  
 Units (Pediatric/Neo): mL  
 Label: Volume  
 AVEA GUI/Membrane: Volume

#### 10.5.1.31 ID: **SetVolAssured**

Description: Minimum Volume desired from delivery of a Pressure Controlled breath.

Type: WORD  
 Scale (Adult): 4  
 Scale (Pediatric/Neo): 1  
 Resolution (Adult): 100  
 Resolution (Pediatric): 10  
 Resolution (Neonate): 1  
 Range (Adult): 0 - 25000  
 Range (Pediatric): 0 - 5000  
 Range (Neonate): 0 - 3000  
 Units (Adult): L  
 Units (Pediatric/Neo): mL  
 Label: Mach Vol  
 AVEA GUI/Membrane: Mach Vol

#### 10.5.1.32 ID: **SetVolLimit**

Description: Maximum Volume Limit for a Pressure Controlled breath.

Type: WORD  
 Scale (Adult): 4  
 Scale (Pediatric/Neo): 1  
 Resolution (Adult): 100  
 Resolution (Pediatric): 10  
 Resolution (Neonate): 1  
 Range (Adult): 1000 - 25000  
 Range (Pediatric): 250 – 7500  
 Range (Neonate): 20 – 3000  
 Units (Adult): L  
 Units (Pediatric/Neo): mL  
 Label: Vol Limit  
 AVEA GUI/Membrane: Vol Limit

- 10.5.1.33 ID: SetVolSigh**  
 Description: Enabled/Disabled state for the Sigh feature for Volume Controlled breaths. (Periodically delivers a proportionally higher volume than set.)  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 1  
 Units: 0=Off 1=On  
 Label: Sigh  
 AVEA GUI/Membrane: Sigh
- 10.5.1.34 ID: SetVolWave**  
 Description: Enumeration for the selected waveform shape (e.g., square or decelerating).  
 Type: ENUM  
 Enum value = label: 0=SQUARE; 1=DECELERATING  
 Label: Waveform  
 AVEA GUI/Membrane: Waveform
- 10.5.1.35 ID: SetVsync**  
 Description: Enabled/Disabled state for the Vsync feature for Volume Controlled breaths. (Vsync is a Pressure Controlled, Volume Targeted modification of Volume Control.)  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 1  
 Units: 0=Off 1=On  
 Label: Vsync  
 AVEA GUI/Membrane: Vsync
- 10.5.1.36 ID: LimitApnea**  
 Description: Time duration limit from the last onset of inspiration until the APNEA INTERVAL alarm is asserted.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 6 - 60  
 Units: sec  
 Label: Apnea Interval  
 AVEA GUI/Membrane: Apnea Interval
- 10.5.1.37 ID: LimitEndTidalCO2High**  
 Description: High End Tidal CO2 limit.  
 Type: WORD  
 Scale: 0 (mmHg); 1 (kPa)  
 Resolution: 1  
 Range (Adult/Ped/Neo): 6 – 150 (mmHg); 8 – 200 (kPa)

Units: mmHg; kPa  
 Label: High EtCO<sub>2</sub>  
 AVEA GUI/Membrane: High EtCO<sub>2</sub>

**10.5.1.38 ID: LimitEndTidalCO2Low**  
 Description: Low End Tidal CO<sub>2</sub> limit.  
 Type: WORD  
 Scale: 0 (mmHg); 1 (kPa)  
 Resolution: 1  
 Range (Adult/Ped/Neo): 1 – 145 (mmHg); 1 – 193 (kPa)  
 Units: mmHg; kPa  
 Label: Low EtCO<sub>2</sub>  
 AVEA GUI/Membrane: Low EtCO<sub>2</sub>

**10.5.1.39 ID: LimitFiO2AutoHigh**  
 Description: Upper bound for the High Auto FiO<sub>2</sub> Limit to trigger alarm.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 21 – 100  
 Units: %  
 Label: High Auto FiO<sub>2</sub>  
 AVEA GUI/Membrane: High Auto FiO<sub>2</sub>

**10.5.1.40 ID: LimitFiO2AutoLow**  
 Description: Lower bound for the Low Auto FiO<sub>2</sub> Limit to trigger alarm.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 21 – 100  
 Units: %  
 Label: Low Auto FiO<sub>2</sub>  
 AVEA GUI/Membrane: Low Auto FiO<sub>2</sub>

**10.5.1.41 ID: LimitFiO2BaselineHigh**  
 Description: Upper bound for the Baseline FiO<sub>2</sub> to trigger an alarm.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 21 – 100  
 Units: %  
 Label: High Base FiO<sub>2</sub>  
 AVEA GUI/Membrane: High Base FiO<sub>2</sub>

**10.5.1.42 ID: LimitPeepLow**  
 Description: Pressure limit for the LOW PEEP alarm. Alarm asserts when the airway pressure drops below the pressure limit.  
 Type: WORD

Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 60  
 Units: cmH2O  
 Label: Low PEEP  
 AVEA GUI/Membrane: Low PEEP

- 10.5.1.43 ID: LimitPeakHigh**  
 Description: Pressure limit for the HIGH  $P_{PEAK}$  alarm. Alarm is asserted when airway pressure exceeds the pressure limit.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped): 10 - 105  
 Range (Neonate): 10 - 85  
 Units: cmH2O  
 Label: High Ppeak  
 AVEA GUI/Membrane: High Ppeak
- 10.5.1.44 ID: LimitPeakLow**  
 Description: Pressure limit for the LOW  $P_{PEAK}$  alarm. Alarm is asserted if Peak Pressure ( $P_{PEAK}$ ) does not exceed the pressure limit by the end of the breath cycle.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped): 1 – 99  
 Range (Neo): 1 – 80  
 Units: cmH2O  
 Label: Low Ppeak  
 AVEA GUI/Membrane: Low Ppeak
- 10.5.1.45 ID: LimitPulseRateHigh**  
 Description: Upper bound of pulse rate to trigger an alarm  
 Type: WORD  
 Resolution: 5  
 Range (Adult/Ped/Neo): 30 – 240  
 Units: bpm  
 Label: High Pulse Limit  
 AVEA GUI/Membrane: High Pulse Limit
- 10.5.1.46 ID: LimitPulseRateLow**  
 Description: Lower bound of pulse rate to trigger an alarm  
 Type: WORD  
 Resolution: 5  
 Range (Adult/Ped/Neo): 25 – 235  
 Units: bpm  
 Label: Low Pulse Limit

## AVEA GUI/Membrane: Low Pulse Limit

- 10.5.1.47 ID: LimitRateHigh**  
 Description: Rate (breaths per minute) limit for the High Rate alarm. Alarm is asserted when monitored total breath rate exceeds the rate limit.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 1 - 200  
 Units: bpm  
 Label: High Rate  
 AVEA GUI/Membrane: High Rate
- 10.5.1.48 ID: LimitSpO2High**  
 Description: Upper bound SpO2 to trigger an alarm  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 3 – 99  
 Units: %  
 Label: High SpO2  
 AVEA GUI/Membrane: High SpO2
- 10.5.1.49 ID: LimitSpO2Low**  
 Description: Lower bound SpO2 to trigger an alarm  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 2 – 98  
 Units: %  
 Label: Low SpO2  
 AVEA GUI/Membrane: Low SpO2
- 10.5.1.50 ID: LimitVeHigh**  
 Description: Volume limit for the HIGH Ve (High Exhaled Minute Volume) alarm. Alarm is asserted when the monitored Minute Volume exceeds the volume limit.  
 Type: WORD  
 Scale: 2  
 Resolution (Adult): 100  
 Resolution (Pediatric): 10  
 Resolution (Neo): 1  
 Range (Adult): 0 – 7500  
 Range (Pediatric): 0 – 3000  
 Range (Neo): 0 – 500  
 Units: L  
 Label: High Ve  
 AVEA GUI/Membrane: High Ve

- 10.5.1.51 ID: LimitVeLow**  
 Description: Volume limit for the LOW Ve (Low Exhaled Minute Volume) alarm. Alarm is asserted when the monitored Minute Volume drops below the volume limit.  
 Type: WORD  
 Scale: 2  
 Resolution (Adult): 100  
 Resolution (Pediatric): 10  
 Resolution (Neo): 1  
 Range (Adult): 0 – 5000  
 Range (Pediatric): 0 – 3000  
 Range (Neonate): 0 - 500  
 Units: L  
 Label: Low Ve  
 AVEA GUI/Membrane: Low Ve
- 10.5.1.52 ID: LimitVteHigh**  
 Description: Volume limit for the HIGH Vte (High Exhaled Tidal Volume) alarm. Alarm is asserted when the monitored Tidal Volume exceeds the volume limit.  
 Type: WORD  
 Scale (Adult): 4  
 Scale (Pediatric/Neo): 1  
 Resolution (Adult): 100  
 Resolution (Pediatric): 10  
 Resolution (Neo): 1  
 Range (Adult): 1000 - 30000  
 Range (Pediatric): 250 - 10000  
 Range (Neonate): 20 - 3000  
 Units (Adult): L  
 Units (Pediatric/Neo): mL  
 Label: High Vte  
 AVEA GUI/Membrane: High Vte
- 10.5.1.53 ID: LimitVteLow**  
 Description: Volume limit for the LOW Vte (Low Exhaled Tidal Volume) alarm. Alarm is asserted when the monitored Tidal Volume is less than volume limit for the number of consecutive breaths set by the Low Vte Sensitivity setting.  
 Type: WORD  
 Scale (Adult): 4  
 Scale (Pediatric/Neo): 1  
 Resolution (Adult): 100  
 Resolution (Pediatric): 10  
 Resolution (Neo): 1



Range (Adult): 0 - 30000  
 Range (Pediatric): 0 - 10000  
 Range (Neonate): 0 - 3000  
 Units (Adult): L  
 Units (Pediatric/Neo): mL  
 Label: Low Vte  
 AVEA GUI/Membrane: Low Vte

**10.5.1.54 ID: SetAAC**  
 Description: Enabled/Disabled state for AAC (Artificial Airway Compensation) which augments pressures to compensate for ET Tube pressure drop.  
 Type: BOOL  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 1  
 Label: ACC On  
 AVEA GUI/Membrane: ACC (On/Off)

**10.5.1.55 ID: SetCircComp**  
 Description: Constant setting for Circuit Compliance Compensation. Defines volume stored in the patient circuit per unit pressure.  
 Type: WORD  
 Scale: 1  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 75  
 Units: mL/cmH2O  
 Label: Circ Comp  
 AVEA GUI/Membrane: Circ Comp

**10.5.1.56 ID: SetEttDia**  
 Description: Diameter of the patient Endotracheal Tube.  
 Type: WORD  
 Scale: 1  
 Resolution: 5  
 Range (Adult/Ped/Neo): 20 - 100  
 Units: mm  
 Label: Diameter  
 AVEA GUI/Membrane: Diameter

**10.5.1.57 ID: SetEttLen**  
 Description: Length of the patient Endotracheal Tube.  
 Type: WORD  
 Scale: 1  
 Resolution: 1  
 Range (Adult): 20 - 300

Range (Pediatric): 20 - 260  
 Range (Neonate): 20 - 150  
 Units: cm  
 Label: Length  
 AVEA GUI/Membrane: Length

**10.5.1.58 ID: SetHumidifier**  
 Description: Active/Passive state of the airway humidifier selection.  
 I.e. Active Humidifier selected or Passive Humidifier selected  
 Type: BOOL  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 1  
 Label: HUMIDIFIER  
 AVEA GUI/Membrane: Active (On/Off)

**10.5.1.59 ID: SetLanguage**  
 Description: Identifies the natural language of the user interface that has been selected.  
 Type: ENUM  
 Enum value = label: 0=English; 1=汉语; 2= Čeština; 3=Nederlands;  
 4=Français; 5=Deutsch; 6=Αγγλικά; 7=Magyar;  
 8=Italiano; 9=日本語; 10=Polski; 11=Português;  
 12=Русский; 13=Español; 14=Türkçe  
 Label: Language  
 AVEA GUI/Membrane: Language

**10.5.1.60 ID: SetLeakComp**  
 Description: Enabled/Disabled state for airway Leak Compensation.  
 When enabled, adds flow during exhalation to maintain PEEP.  
 Type: BOOL  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 1  
 Label: Leak Comp On  
 AVEA GUI/Membrane: Leak Comp (On/Off)

**10.5.1.61 ID: SetMode**  
 Description: Identifies the ventilation mode and breath-type that the ventilator is set to deliver.  
 Type: ENUM  
 Enum value = label: 1=APRV / BIPHASIC; 2=APRV / BIPHASIC; 3= PRVC SIMV;  
 4=PRVC A/C; 5=CPAP / PSV; 6=TCPL SIMV;  
 7=TCPL A/C; 8=CPAP / PSV; 9=PRESSURE SIMV;

10=PRESSURE A/C; 11=CPAP / PSV; 12=VOLUME  
SIMV; 13=VOLUME A/C; 14=Nasal CPAP; 15=SiPAP  
Label: MODE SELECT  
AVEA GUI/Membrane: Mode

- 10.5.1.62 ID: SetModellv**  
Description: Identifies the activity and role (i.e. OFF/Master/Slave) of the Independent Lung Ventilation feature.  
Type: ENUM  
Enum value = label: 0=Off; 1=Master; 2=Slave  
Label: ILV Mode  
AVEA GUI/Membrane: ILV Mode
- 10.5.1.63 ID: SetPatSize**  
Description: Identifies the size of the patient – neonate, pediatric, or adult.  
Type: ENUM  
Enum value = label: 0=Neo; 1=Ped; 2=Adult  
Label: Patient Size  
AVEA GUI/Membrane: Patient Size
- 10.5.1.64 ID: SetPatWt**  
Description: Weight of the patient.  
Type: WORD  
Scale (Adult/Ped/Neo): 2  
Resolution (Adult): 100  
Resolution (Pediatric): 10  
Resolution (Neonate): 1  
Range (Adult): 100 - 30000  
Range (Pediatric): 10 - 7500  
Range (Neonate): 10 - 1600  
Units: kg  
Label: Pt Weight  
AVEA GUI/Membrane: Pt Weight
- 10.5.1.65 ID: SetPresBaro**  
Description: Barometric Pressure. Absolute pressure of the ambient environment.  
Type: WORD  
Resolution: 1  
Range (Adult/Ped/Neo): 545 – 760  
Units: mmHg  
Label: Baro Pres  
AVEA GUI/Membrane: Baro Pres

- 10.5.1.66 ID: SetIncrFiO2**  
 Description: Percentage increase in delivered FiO2 when front membrane panel button “Increase O2” is activated.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 0 - 79  
 Units: %  
 Label: Increase FiO2 % Increment  
 AVEA GUI/Membrane: Increase FiO2 % Increment
- 10.5.1.67 ID: SetSensitivityLowVte**  
 Description: The number of consecutive breaths that are in violation of the Low Vte alarm limit that are required to assert the Low Vte Alarm.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 1 - 5  
 Label: Low Vte Alarm:  
 AVEA GUI/Membrane: Low Vte Alarm Sensitivity
- 10.5.1.68 ID: SetSpO2AlarmDelay**  
 Description: Determines the duration that a High or Low SpO2alarm or “SpO2 Invalid” alarm must persist before the alarm is activated.  
 Type: WORD  
 Range: 0 - 120  
 Resolution: 5  
 Units: seconds  
 Label: SpO2 Alrm Delay  
 AVEA GUI/Membrane: SpO2 Alrm Delay
- 10.5.1.69 ID: SetFiO2AutoControlEnable**  
 Description: Activates and de-activates the Automatic FiO2 Control system.  
 Type: ENUM  
 Enum value = label: 0=Disabled; 1=Enabled  
 Label: Auto FiO2  
 AVEA GUI/Membrane: Auto FiO2
- 10.5.1.70 ID: SetSpO2AveragingTime**  
 Description: Determines the time period for averaging the output of the SpO2 readings.  
 Type: WORD  
 Range: 2 - 16  
 Resolution: 2  
 Units: seconds

Label: Oximeter Averaging  
 AVEA GUI/Membrane: Oximeter Averaging

- 10.5.1.71 ID: SetSpO2Enable**  
 Description: Enable or Disables Pulse Oximeter related functions.  
 Type: ENUM  
 Enum value = label: 0=Disabled; 1=Enabled  
 Label: Oximeter Enable  
 AVEA GUI/Membrane: Oximeter Enable
- 10.5.1.72 ID: SetSpO2Mode**  
 Description: Configures the pulse oximeter module to a specified operating mode.  
 Type: ENUM  
 Enum value = label: 1=Maximum; 2=Normal; 3=APOD  
 Label: Oximeter Algorithm  
 AVEA GUI/Membrane: Oximeter Algorithm
- 10.5.1.73 ID: SetSpO2TargetHigh**  
 Description: Auto FiO2 Cmd is greater than or equal to the preset High Auto FiO2 Limit for a period of 60 seconds or more.  
 Type: WORD  
 Range: 82 – 100  
 Resolution: 1  
 Units: %  
 Label: High SpO2 Trgt  
 AVEA GUI/Membrane: High SpO2 Trgt
- 10.5.1.74 ID: SetSpO2TargetLow**  
 Description: Auto FiO2 Cmd is less than or equal to the preset Low Auto FiO2 Limit for a period of 60 seconds or more.  
 Type: WORD  
 Range: 80 – 98  
 Resolution: 1  
 Units: %  
 Label: Low SpO2 Trgt  
 AVEA GUI/Membrane: Low SpO2 Trgt

## 10.5.2 VELA Ventilator – Setting Class

- 10.5.2.1 ID: SetCO2MonitoringEnable**  
 Description: Enable or Disable CO2 Monitoring  
 Type: ENUM  
 Enum value = label: 0=Disabled; 1=Enabled  
 Label: CO2  
 VELA GUI/Membrane: CO2

- 10.5.2.2 ID: SetEndTidalCO2Average**  
 Description: Number of breaths EtCO2 is averaged over.  
 Type: WORD  
 Resolution: 1  
 Range (Adult/Ped/Neo): 1 - 8  
 Units: breath  
 Label: EtCO2 Avg.  
 VELA GUI/Membrane: EtCO2 Avg.
- 10.5.2.3 ID: SetFiO2**  
 Description: Percent of oxygen that the ventilator is set to administer.  
 Type: WORD  
 Resolution: 1  
 Range: 21 - 100  
 Units: %  
 Label: FiO2  
 VELA GUI/Membrane: O2
- 10.5.2.4 ID: SetFiO2IncreaseActive**  
 Description: Flow delivered to the breathing circuit during the expiratory phase.  
 Type: BOOL  
 Resolution: 1  
 Range: 0 - 1  
 Units: 0=Off 1=On  
 Label: 100% O2  
 VELA GUI/Membrane: 100% O2
- 10.5.2.5 ID: SetFlowBias**  
 Description: Flow delivered to the breathing circuit during the expiratory phase.  
 Type: WORD  
 Scale: 1  
 Resolution: 10  
 Range: 100 - 200  
 Units: L/min  
 Label: Bias Flow  
 VELA GUI/Membrane: Bias Flow
- 10.5.2.6 ID: SetFlowCycle**  
 Description: Percent of the peak inspiratory flow (Peak Flow), at which the inspiratory phase of Pressure Controlled (including TCPL) breaths are terminated.  
 Type: WORD  
 Resolution: 5  
 Range: 5 - 70

Units: %  
 Label: PC Flow Cycle  
 VELA GUI/Membrane: PC Flow Cycle

- 10.5.2.7 ID: SetFlowCyclePsv**  
 Description: Percent of peak inspiratory flow at which the inspiratory phase of a PSV breath is terminated.  
 Type: WORD  
 Resolution: 5  
 Range: 5 - 30  
 Units: %  
 Label: PSV Cycle  
 VELA GUI/Membrane: PSV Cycle
- 10.5.2.8 ID: SetFlowInsp**  
 Description: Flow targeted by the ventilator during the inspiratory phase of flow controlled breaths.  
 Type: WORD  
 Scale: 1  
 Resolution: 1  
 Range: 100 - 1400  
 Units: L/min  
 Label: Peak Flow  
 VELA GUI/Membrane: Peak Flow
- 10.5.2.9 ID: SetPanelLockActive**  
 Description: The front panel is locked for user input.  
 Type: BOOL  
 Resolution: 1  
 Range: 0 - 1  
 Units: 0=Off 1=On  
 Label: PANEL LOCK  
 VELA GUI/Membrane: PANEL LOCK
- 10.5.2.10 ID: SetPauseInsp**  
 Description: Time Inspiration is extended before exhalation occurs after a volume breath is delivered.  
 Type: WORD  
 Scale: 2  
 Resolution: 10  
 Range: 0 - 200  
 Units: Sec  
 Label: Insp Pause  
 VELA GUI/Membrane: Insp Pause

<b>10.5.2.11</b>	<b>ID:</b>	<b>SetPresHigh</b>
	Description:	Baseline target for the “high” phase of APRV/BiPhasic mode.
	Type:	WORD
	Resolution:	1
	Range:	0 - 60
	Units:	cmH2O
	Label:	Pres High
	VELA GUI/Membrane:	Pres High
<b>10.5.2.12</b>	<b>ID:</b>	<b>SetPresInsp</b>
	Description:	Pressure target for mandatory Pressure Controlled breaths.
	Type:	WORD
	Resolution:	1
	Range:	1 - 100
	Units:	cmH2O
	Label:	Insp Pres
	VELA GUI/Membrane:	Insp Pres
<b>10.5.2.13</b>	<b>ID:</b>	<b>SetPresInspNPPV</b>
	Description:	Pressure target for mandatory Pressure Controlled breaths, NPPV mode.
	Type:	WORD
	Resolution:	1
	Range:	1 - 40
	Units:	cmH2O
	Label:	NPPV P <sub>insp</sub>
	VELA GUI/Membrane:	NPPV P <sub>insp</sub>
<b>10.5.2.14</b>	<b>ID:</b>	<b>SetPresLow</b>
	Description:	Baseline target for the “low” phase of APRV/BiPhasic mode.
	Type:	WORD
	Resolution:	1
	Range:	0 - 45
	Units:	cmH2O
	Label:	Pres Low
	VELA GUI/Membrane:	Pres Low
<b>10.5.2.15</b>	<b>ID:</b>	<b>SetPresPeep</b>
	Description:	Positive End Expiratory Pressure. I.e., target pressure to be maintained at the end of exhalation.
	Type:	WORD
	Resolution:	1
	Range:	0 - 35



Units: cmH2O  
Label: PEEP  
VELA GUI/Membrane: PEEP

**10.5.2.16 ID: SetPresPsv**

Description: Pressure target for spontaneous (PSV or Pressure Support Ventilation) breaths.

Type: WORD  
Resolution: 1  
Range: 0 - 60  
Units: cmH2O  
Label: PSV  
VELA GUI/Membrane: PSV

**10.5.2.17 ID: SetPresPsvNPPV**

Description: Pressure target for spontaneous (PSV or Pressure Support Ventilation) breaths, NPPV mode.

Type: WORD  
Resolution: 1  
Range: 0 - 40  
Units: cmH2O  
Label: NPPV PSV  
VELA GUI/Membrane: NPPV PSV

**10.5.2.18 ID: SetRate**

Description: Mandatory breath rate setting.

Type: WORD  
Resolution: 1  
Range: 2 - 80  
Units: bpm  
Label: Rate  
VELA GUI/Membrane: Rate

**10.5.2.19 ID: SetTimeHigh**

Description: Time duration setting for the inspiratory phase of Pressure Controlled (including TCPL) mandatory breaths.

Type: WORD  
Scale: 1  
Resolution: 1  
Range: 3 - 300  
Units: Sec  
Label: Time High  
VELA GUI/Membrane: Time High

<b>10.5.2.20</b>	<b>ID:</b>	<b>SetTimeHighPsv</b>
	Description:	Enabled/Disabled state of PSV breaths is available during Time High by activating T High PSV (an advanced setting of Time High).
	Type:	BOOL
	Resolution:	1
	Range:	0 - 1
	Units:	0=Off 1=On
	Label:	T High PSV
	VELA GUI/Membrane:	T High PSV
<b>10.5.2.21</b>	<b>ID:</b>	<b>SetTimeHighSync</b>
	Description:	Percentage of Time High in APRV/BiPhasic mode allowed at the end of the Time High period for synchronizing the high to low baseline transition with a patient exhalation.
	Type:	WORD
	Resolution:	5
	Range:	0 - 50
	Units:	%
	Label:	T High Sync
	VELA GUI/Membrane:	T High Sync
<b>10.5.2.22</b>	<b>ID:</b>	<b>SetTimeInsp</b>
	Description:	Time duration setting for the inspiratory phase of Pressure Controlled (including TCPL) mandatory breaths.
	Type:	WORD
	Scale:	2
	Resolution:	10
	Range:	30 - 1000
	Units:	Sec
	Label:	Insp Time
	VELA GUI/Membrane:	Insp Time
<b>10.5.2.23</b>	<b>ID:</b>	<b>SetTimeLow</b>
	Description:	Time duration setting for the "low" phase of APRV/BiPhasic mode.
	Type:	WORD
	Scale:	1
	Resolution:	1
	Range:	3 - 300
	Units:	Sec
	Label:	Time Low
	VELA GUI/Membrane:	Time Low

- 10.5.2.24 ID: SetTimeLowSync**  
 Description: Percentage of Time Low in APRV/BiPhasic mode allowed at the end of the Time Low period for synchronizing the low to high baseline transition with a patient inspiration.  
 Type: WORD  
 Resolution: 5  
 Range: 0 - 50  
 Units: %  
 Label: T Low Sync  
 VELA GUI/Membrane: T Low Sync
- 10.5.2.25 ID: SetTmaxPsv**  
 Description: Maximum time duration allowed for the inspiratory phase of a Pressure Supported breath.  
 Type: WORD  
 Scale: 2  
 Resolution: 10  
 Range: 30 - 3000  
 Units: Sec  
 Label: PSV Tmax  
 VELA GUI/Membrane: PSV Tmax
- 10.5.2.26 ID: SetTrigFlow**  
 Description: Flow Trigger level. I.e., Net Flow level at which inspiration is triggered.  
 Type: WORD  
 Scale: 1  
 Resolution: 1  
 Range: 10 - 200  
 Units: L/min  
 Label: Flow Trig  
 VELA GUI/Membrane: Flow Trig
- 10.5.2.27 ID: SetVCO2Average**  
 Description: Patient's exhaled minute volume of CO<sub>2</sub>.  
 Type: WORD  
 Resolution: 3  
 Range (Adult/Ped/Neo): 3 - 12  
 Units: Minute  
 Label: VCO<sub>2</sub> Avg.  
 VELA GUI/Membrane: VCO<sub>2</sub> Avg.
- 10.5.2.28 ID: SetVol**  
 Description: Volume targeted by the ventilator for delivery of Volume Controlled breaths.

Type: WORD  
 Scale: 1  
 Resolution: 50  
 Range: 500 - 20000  
 Units: ml  
 Label: Volume  
 VELA GUI/Membrane: Volume

**10.5.2.29 ID: SetVolAssured**  
 Description: Minimum Volume desired from delivery of a Pressure Controlled breath.

Type: WORD  
 Scale: 1  
 Resolution: 50  
 Range: 500 - 20000  
 Units: ml  
 Label: Assured Vol  
 VELA GUI/Membrane: Assured Vol

**10.5.2.30 ID: SetVolLimit**  
 Description: Maximum Volume Limit for a Pressure Controlled breath.

Type: WORD  
 Scale: 1  
 Resolution: 50  
 Range: 500 - 25000  
 Units: ml  
 Label: Vol Limit  
 VELA GUI/Membrane: Vol Limit

**10.5.2.31 ID: SetVolSigh**  
 Description: Enabled/Disabled state for the Sigh feature for Volume Controlled breaths. (Periodically delivers a proportionally higher volume than set.)

Type: BOOL  
 Resolution: 1  
 Range: 0 - 1  
 Units: 0=Off 1=On  
 Label: Sigh  
 VELA GUI/Membrane: Sigh

**10.5.2.32 ID: SetVolWave**  
 Description: Enumeration for the selected waveform shape (e.g., square or decelerating).

Type: ENUM  
 Enum value=label: 0=SQUARE; 1=DECELERATING  
 Label: Waveform

VELA GUI/Membrane: Waveform

- 10.5.2.33 ID: SetVsync**  
 Description: Enabled/Disabled state for the Vsync feature for Volume Controlled breaths. (Vsync is a Pressure Controlled, Volume Targeted modification of Volume Control.)  
 Type: BOOL  
 Resolution: 1  
 Range: 0 - 1  
 Units: 0=Off 1=On  
 Label: Vsync  
 VELA GUI/Membrane: Vsync
- 10.5.2.34 ID: LimitApnea**  
 Description: Time duration limit from the last onset of inspiration until the APNEA INTERVAL alarm is asserted.  
 Type: WORD  
 Resolution: 1  
 Range: 10 - 60  
 Units: Sec  
 Label: Apnea Interval  
 VELA GUI/Membrane: Apnea Interval
- 10.5.2.35 ID: LimitEndTidalCO2High**  
 Description: High End Tidal CO2 limit.  
 Type: WORD  
 Scale: 0 (mmHg); 1 (kPa)  
 Resolution: 1  
 Range (Adult/Ped/Neo): 6 – 150 (mmHg); 8 - 200 (kPa)  
 Units: mmHg, kPa  
 Label: High EtCO2  
 VELA GUI/Membrane: High EtCO2
- 10.5.2.36 ID: LimitEndTidalCO2Low**  
 Description: Low End Tidal CO2 limit.  
 Type: WORD  
 Scale: 0 (mmHg); 1 (kPa)  
 Resolution: 1  
 Range (Adult/Ped/Neo): 1 – 145 (mmHg); 1 – 193 (kPa)  
 Units: mmHg; kPa  
 Label: Low EtCO2  
 VELA GUI/Membrane: Low EtCO2

<b>10.5.2.37</b>	<b>ID:</b>	<b>LimitPpeakHigh</b>
	Description:	Pressure limit for the HIGH P <sub>PEAK</sub> alarm. Alarm is asserted when airway pressure exceeds the pressure limit.
	Type:	WORD
	Resolution:	1
	Range:	5 - 120
	Units:	cmH2O
	Label:	High Ppeak
	VELA GUI/Membrane:	High Ppeak
<b>10.5.2.38</b>	<b>ID:</b>	<b>LimitPpeakLow</b>
	Description:	Pressure limit for the LOW P <sub>PEAK</sub> alarm. Alarm is asserted if Peak Pressure (P <sub>PEAK</sub> ) does not exceed the pressure limit by the end of the breath cycle.
	Type:	WORD
	Resolution:	1
	Range:	2 - 60
	Units:	cmH2O
	Label:	Low Ppeak
	VELA GUI/Membrane:	Low Ppeak
<b>10.5.2.39</b>	<b>ID:</b>	<b>LimitRateHigh</b>
	Description:	Rate (breaths per minute) limit for the High Rate alarm. Alarm is asserted when monitored total breath rate exceeds the rate limit.
	Type:	WORD
	Resolution:	1
	Range:	3 - 150
	Units:	bpm
	Label:	High Rate
	VELA GUI/Membrane:	High Rate
<b>10.5.2.40</b>	<b>ID:</b>	<b>LimitVeLow</b>
	Description:	Volume limit for the LOW Ve (Low Exhaled Minute Volume) alarm. Alarm is asserted when the monitored Minute Volume drops below the volume limit.
	Type:	WORD
	Scale:	2
	Resolution:	10
	Range:	10 - 9990
	Units:	L
	Label:	Low Ve
	VELA GUI/Membrane:	Low Ve

<b>10.5.2.41</b>	<b>ID:</b>	<b>SetAltitude</b>
	Description:	Altitude at which the ventilator will be operating.
	Type:	WORD
	Resolution:	30
	Range:	300 – 3060
	Units:	Meters
	Label:	Altitude
	VELA GUI/Membrane:	Altitude
<b>10.5.2.42</b>	<b>ID:</b>	<b>SetFiO2Monitoring</b>
	Description:	Enable / disable FiO2 monitoring.
	Type:	ENUM
	Enum value=label:	0=FiO2 Monitor Disabled; 1=FiO2 Monitor Enabled
	Label:	FiO2 Monitor Enabled
	VELA GUI/Membrane:	FiO2 Monitor (Enabled/Disabled)
<b>10.5.2.43</b>	<b>ID:</b>	<b>SetHumidifier</b>
	Description:	Active/Passive state of the airway humidifier selection. I.e. Active Humidifier selected or Passive Humidifier selected
	Type:	BOOL
	Resolution:	1
	Range:	0 - 1
	Units:	0=Off 1=On
	Label:	Humidifier
	VELA GUI/Membrane:	Humidifier
<b>10.5.2.44</b>	<b>ID:</b>	<b>SetLanguage</b>
	Description:	Identifies the natural language of the user interface that has been selected.
	Type:	ENUM
	Enum value=label:	0=Deutsch; 1=English; 2=中文; 3= Español; 4= Français; 5=Italiano; 6=日本語; 7=Polski; 8=Türkçe; 9=Русский; 10= Nederlands; 11=Magyar; 12=Čeština
	Label:	Language
	VELA GUI/Membrane:	Language
<b>10.5.2.45</b>	<b>ID:</b>	<b>SetLeakComp</b>
	Description:	Enabled/Disabled state for airway Leak Compensation. When enabled, adds flow during exhalation to maintain PEEP.
	Type:	BOOL
	Resolution:	1
	Range:	0 - 1
	Units:	0=Off 1=On

Label: Lk Comp  
VELA GUI/Membrane: Lk Comp

- 10.5.2.46 ID: SetMode**  
Description: Identifies the ventilation mode and breath-type that the ventilator is set to deliver.  
Type: ENUM  
Enum value=label: 1=NPPV A/C; 2=NPPV SIMV; 3=NPPVCPAP PSV-Volume; 4= NPPVCPAP PSV-Pressure; 5=APRV BiPhasic-Volume; 6=APRV BiPhasic-Pressure; 7=PRVC SIMV; 8=PRVC A/C; 9=CPAP PSV-Pressure; 10=PRESSURE SIMV; 11= PRESSURE A/C; 12=CPAP PSV-Volume; 13=VOLUME SIMV; 14=VOLUME A/C  
Label: MODE SELECT  
VELA GUI/Membrane: MODE SELECT
- 10.5.2.47 ID: SetNebulizerActive**  
Description: Enable / disable the nebulizer (button on the front panel)  
Type: BOOL  
Resolution: 1  
Range: 0 - 1  
Units: 0=Off 1=On  
Label: NEBULIZER  
VELA GUI/Membrane: NEBULIZER
- 10.5.2.48 ID: SetNebulizerTime**  
Description: Specify the number of minutes the nebulizer will be active.  
Type: WORD  
Resolution: 1  
Range: 1 - 60  
Units: Minutes  
Label: Neb Time  
VELA GUI/Membrane: Neb Time
- 10.5.2.49 ID: SetPanelLockEnable**  
Description: Enable / disable the ability to lock the front panel from user input.  
Type: ENUM  
Enum value = label: 0=Locks Disabled; 1=Locks Enabled  
Label: Locks Enabled  
VELA GUI/Membrane: Locks Enabled
- 10.5.2.50 ID: SetVeLowOffEnable**  
Description: Enable / disable the ability to turn OFF the Low Ve alarm limit



Type: ENUM  
 Enum value = label: 0=Low Min Vol Off Disabled; 1=Low Min Vol Off Enabled.  
 NOTE: There is a line feed character between the "Vol" and "Off" strings.  
 Label: Low Min Vol Off Enabled  
 VELA GUI/Membrane: Low Min Vol Off (Disabled / Enabled)

## 10.6 Monitor Class

### 10.6.1 AVEA Ventilator – Monitor Class

- 10.6.1.1 ID: MntrAutoPEEP**  
 Description: AutoPEEP. Pressure measured in the airway resulting from performance of an AutoPEEP Maneuver or Expiratory Hold  
 Type: WORD  
 Range: 0 - 50  
 Units: cmH2O  
 Label: AutoPEEP  
 AVEA GUI/Membrane: AutoPEEP
- 10.6.1.2 ID: MntrAutoPEEPdelta**  
 Description: Differential AutoPEEP. Pressure difference between baseline pressure (PEEP) and AutoPEEP at the time AutoPEEP was measured.  
 Type: WORD  
 Range: 0 - 50  
 Units: cmH2O  
 Label: dAutoPEEP  
 AVEA GUI/Membrane: dAutoPEEP
- 10.6.1.3 ID: MntrAutoPEEPesoph**  
 Description: Esophageal AutoPEEP. Pressure measured via an Esophageal Balloon and processed to calculate result.  
 Type: WORD  
 Range: 0 - 50  
 Units: cmH2O  
 Label: AutoPEEPes  
 AVEA GUI/Membrane: AutoPEEPes
- 10.6.1.4 ID: MntrC20**  
 Description: Compliance Ratio. The ratio of the dynamic compliance during the last 20% of inspiration ( $C_{20}$ ) to the total dynamic compliance (C).  
 Type: WORD  
 Scale: 2

Range: 0 - 500  
 Label: C20/C  
 AVEA GUI/Membrane: C20/C

- 10.6.1.5 ID: MntrCcw**  
 Description: Chestwall Compliance. The ratio of exhaled tidal volume to the Delta Esophageal Pressure ( $dP_{ES}$ ).  
 Type: WORD  
 Range: 0 – 300  
 Units: mL/cmH2O  
 Label: Ccw  
 AVEA GUI/Membrane: Ccw
- 10.6.1.6 ID: MntrCdyn**  
 Description: Dynamic Compliance. Volume capacity of patient per unit pressure measured while system is not at equilibrium.  
 Type: WORD  
 Scale (Neo): 2  
 Range (Adult/Ped): 0 – 300  
 Range (Neo): 0 – 30000  
 Units: mL/cmH2O  
 Label: Cdyn  
 AVEA GUI/Membrane: Cdyn
- 10.6.1.7 ID: MntrCdynNorm**  
 Description: Normalized Dynamic Compliance. Ratio of  $C_{DYN}$  to Patient Weight.  
 Type: WORD  
 Scale: 2  
 Range: 0 - 500  
 Units: mL/cmH2O/kg  
 Label: Cdyn/kg  
 AVEA GUI/Membrane: Cdyn/kg
- 10.6.1.8 ID: MntrClung**  
 Description: Lung Compliance. Volume capacity of patient's lungs per unit pressure.  
 Type: WORD  
 Range: 0 – 300  
 Units: mL/cmH2O  
 Label: Clung  
 AVEA GUI/Membrane: Clung

<b>10.6.1.9</b>	<b>ID:</b>	<b>MntrCstat</b>
	Description:	Static Compliance. A.k.a. Respiratory System Compliance ( $C_{RS}$ ). Volume capacity of patient per unit pressure measured while system is at equilibrium (i.e. Inspiratory Hold, Inspiratory Pause, Pressure Plateau, etc.).
	Type:	WORD
	Scale (Neo):	2
	Range (Adult/Ped):	0 – 300
	Range (Neo):	0 – 30000
	Units:	mL/cmH <sub>2</sub> O
	Label:	Cstat
	AVEA GUI/Membrane:	Cstat
<b>10.6.1.10</b>	<b>ID:</b>	<b>MntrCstatNorm</b>
	Description:	Normalized Static Compliance. Ratio of Cstat to patient Weight.
	Type:	WORD
	Scale:	2
	Range:	0 - 500
	Units:	mL/cmH <sub>2</sub> O/kg
	Label:	Cstat/kg
	AVEA GUI/Membrane:	Cstat/kg
<b>10.6.1.11</b>	<b>ID:</b>	<b>MntrEndTidalCO<sub>2</sub></b>
	Description:	Patient's peak expired CO <sub>2</sub> level as measured and reported by the CO <sub>2</sub> analyzer.
	Type:	WORD
	Scale:	1 (mmHg); 2 (kPa)
	Range:	0 – 1500 (mmHg); 0 – 2000 (kPa)
	Units:	mmHg; kPa
	Label:	EtCO <sub>2</sub>
	AVEA GUI/Membrane:	EtCO <sub>2</sub>
<b>10.6.1.12</b>	<b>ID:</b>	<b>MntrFiO<sub>2</sub></b>
	Description:	Fractional (percent) O <sub>2</sub> measured in inspiratory flow stream to patient.
	Type:	WORD
	Range:	0 – 106
	Units:	%
	Label:	FiO <sub>2</sub>
	AVEA GUI/Membrane:	FiO <sub>2</sub>
<b>10.6.1.13</b>	<b>ID:</b>	<b>MntrFiO<sub>2</sub>Baseline</b>
	Description:	Average FiO <sub>2</sub> required in order to maintain the patient in stable normoxemia over time.

Type: WORD  
 Range: 21 – 100  
 Units: %  
 Label: FiO2 Baseline  
 AVEA GUI/Membrane: FiO2 Baseline

**10.6.1.14 ID: MntrIE**  
 Description: Inspiratory/Expiratory Ratio. Ratio of time spent in the inspiratory phase of a breath cycle to the time spent in the expiratory phase.

Type: WORD  
 Scale: 1  
 Range: -999 - 999  
 Units:  
 Label: I:E  
 AVEA GUI/Membrane: I:E

**10.6.1.15 ID: MntrLeak**  
 Description: Leak. Expressed as the percentage difference between the inspired and exhaled tidal volumes.

Type: WORD  
 Range: 0 - 100  
 Units: %  
 Label: Leak  
 AVEA GUI/Membrane: Leak

**10.6.1.16 ID: MntrMIP**  
 Description: Maximum Inspiratory Pressure. A.k.a. Negative Inspiratory Force (NIF). The maximum negative airway pressure that is achieved by the patient during a MIP/P100 Maneuver or expiratory hold.

Type: WORD  
 Range: -60 - 120  
 Units: cmH2O  
 Label: MIP  
 AVEA GUI/Membrane: MIP

**10.6.1.17 ID: MntrNcpapMeanFlow**  
 Description: Nasal CPAP mean inspiratory flow

Type: WORD  
 Scale: 1  
 Range: 0 - 3000  
 Units: L/min  
 Label: CPAP Flow  
 AVEA GUI/Membrane: CPAP Flow

<b>10.6.1.18</b>	<b>ID:</b>	<b>MntrNcpapPres</b>
	Description:	Nasal CPAP pressure
	Type:	WORD
	Range:	0 - 120
	Units:	cmH2O
	Label:	nCPAP
	AVEA GUI/Membrane:	nCPAP
<b>10.6.1.19</b>	<b>ID:</b>	<b>MntrP100</b>
	Description:	Inspiratory Drive. The negative pressure measured 100ms after onset of a MIP/P100 Maneuver or expiratory hold.
	Type:	WORD
	Range:	-60 - 120
	Units:	cmH2O
	Label:	P100
	AVEA GUI/Membrane:	P100
<b>10.6.1.20</b>	<b>ID:</b>	<b>MntrPair</b>
	Description:	Air Supply Pressure. Pressure measured at the air inlet of the ventilator.
	Type:	WORD
	Range:	0 - 80
	Units:	psig
	Label:	Air Inlet
	AVEA GUI/Membrane:	Air Inlet
<b>10.6.1.21</b>	<b>ID:</b>	<b>MntrPawDelta</b>
	Description:	Differential Airway Pressure. The difference between peak airway pressure and baseline airway pressure throughout an entire breath cycle.
	Type:	WORD
	Range:	-120 - 120
	Units:	cmH2O
	Label:	dPaw
	AVEA GUI/Membrane:	dPaw
<b>10.6.1.22</b>	<b>ID:</b>	<b>MntrPeep</b>
	Description:	Positive End Expiratory Pressure. Airway pressure at the end of the expiratory phase of a breath cycle.
	Type:	WORD
	Range:	0 - 60
	Units:	cmH2O
	Label:	PEEP
	AVEA GUI/Membrane:	PEEP

<b>10.6.1.23</b>	<b>ID:</b>	<b>MntrPefr</b>
	Description:	Peak Expiratory Flow Rate. Highest flow rate measured during the expiratory phase of a breath cycle.
	Type:	WORD
	Scale (Neo):	1
	Range (Adult/Ped):	0 – 300
	Range (Neo):	0 – 3000
	Units:	L/min
	Label:	PEFR
	AVEA GUI/Membrane:	PEFR
<b>10.6.1.24</b>	<b>ID:</b>	<b>MntrPesDelta</b>
	Description:	Differential Esophageal Pressure. The pressure difference between peak esophageal pressure and baseline esophageal pressure throughout an entire breath cycle.
	Type:	WORD
	Range:	-120 - 120
	Units:	cmH2O
	Label:	dPes
	AVEA GUI/Membrane:	dPes
<b>10.6.1.25</b>	<b>ID:</b>	<b>MntrPifr</b>
	Description:	Peak Inspiratory Flow Rate. Highest flow rate measured during the inspiratory phase of a breath cycle.
	Type:	WORD
	Scale (Neo):	1
	Range (Adult/Ped):	0 – 300
	Range (Neo):	0 – 3000
	Units:	L/min
	Label:	PIFR
	AVEA GUI/Membrane:	PIFR
<b>10.6.1.26</b>	<b>ID:</b>	<b>MntrPmean</b>
	Description:	Mean Airway Pressure. Average airway pressure measured throughout an entire breath cycle.
	Type:	WORD
	Range:	0 - 120
	Units:	cmH2O
	Label:	Pmean
	AVEA GUI/Membrane:	Pmean
<b>10.6.1.27</b>	<b>ID:</b>	<b>MntrPO2</b>
	Description:	Oxygen Supply Pressure. Pressure measured at the oxygen inlet of the ventilator.
	Type:	WORD

Range: 0 - 80  
 Units: psig  
 Label: O2 Inlet  
 AVEA GUI/Membrane: O2 Inlet

**10.6.1.28 ID: MntrPpeak**  
 Description: Peak Airway Pressure. Highest airway pressure measured throughout an entire breath cycle.  
 Type: WORD  
 Range: 0 - 120  
 Units: cmH2O  
 Label: Ppeak  
 AVEA GUI/Membrane: Ppeak

**10.6.1.29 ID: MntrPplat**  
 Description: Plateau Pressure. Airway pressure measured when the pressure plateaus (i.e. pressure is neither increasing nor decreasing).  
 Type: WORD  
 Range: 0 - 120  
 Units: cmH2O  
 Label: Pplat  
 AVEA GUI/Membrane: Pplat

**10.6.1.30 ID: MntrPplatPtp**  
 Description: Transpulmonary Pressure, Plateau. Transpulmonary Pressure (Airway Pressure – Esophageal Pressure) measured at the time a Plateau Pressure is measured.  
 Type: WORD  
 Range: -60 - 120  
 Units: cmH2O  
 Label: Ptp Plat  
 AVEA GUI/Membrane: Ptp Plat

**10.6.1.31 ID: MntrPresBaro**  
 Description: Barometric Pressure  
 Type: WORD  
 Scale: 0 (mmHg); 1 (kPa)  
 Range: 545 – 760 (mmHg); 727 – 1013 (kPa)  
 Units: mmHg, kPa  
 Label: Pbaro  
 AVEA GUI/Membrane: Pbaro

**10.6.1.32 ID: MntrPtpPEEP**  
 Description: Transpulmonary Pressure, PEEP. Transpulmonary Pressure (Airway Pressure – Esophageal Pressure)

measured at the time an AutoPEEP measurement is taken.

Type: WORD  
 Range: -60 - 120  
 Units: cmH2O  
 Label: Ptp PEEP  
 AVEA GUI/Membrane: Ptp PEEP

**10.6.1.33 ID: MntrRate**  
 Description: Breath Rate. Number of inspirations per minute of all breath types.  
 Type: WORD  
 Range: 0 - 200  
 Units: bpm  
 Label: Rate  
 AVEA GUI/Membrane: Rate

**10.6.1.34 ID: MntrRateMand**  
 Description: Mandatory Breath Rate. Number of inspirations per minute of machine controlled breath types.  
 Type: WORD  
 Range: 0 - 200  
 Units: bpm  
 Label: Mand Rate  
 AVEA GUI/Membrane: Mand Rate

**10.6.1.35 ID: MntrRateSpon**  
 Description: Spontaneous breath rate. Number of inspirations per minute of patient controlled breath types.  
 Type: WORD  
 Range: 0 - 200  
 Units: bpm  
 Label: Spon Rate  
 AVEA GUI/Membrane: Spon Rate

**10.6.1.36 ID: MntrRimp**  
 Description: Imposed Resistance. Resistance imposed by the artificial airway situated between the circuit wye and the trachea.  
 Type: WORD  
 Scale: 1  
 Range: 0 - 1000  
 Units: cmH2O/L/Sec  
 Label: Rimp  
 AVEA GUI/Membrane: Rimp



<b>10.6.1.37</b>	<b>ID:</b>	<b>MntrRlung</b>
	Description:	Lung Resistance. Resistance attributed solely to structures within the lung.
	Type:	WORD
	Scale:	1
	Range:	0 - 1000
	Units:	cmH2O/L/Sec
	Label:	Rlung
	AVEA GUI/Membrane:	Rlung
<b>10.6.1.38</b>	<b>ID:</b>	<b>MntrRpeak</b>
	Description:	Peak Expiratory Resistance. Highest airway resistance measured throughout the expiratory phase of a breath cycle.
	Type:	WORD
	Scale:	1
	Range:	0 - 1000
	Units:	cmH2O/L/Sec
	Label:	Rpeak
	AVEA GUI/Membrane:	Rpeak
<b>10.6.1.39</b>	<b>ID:</b>	<b>MntrRrs</b>
	Description:	Respiratory System Resistance. Resistance measured during the inspiratory phase of a breath cycle.
	Type:	WORD
	Scale:	1
	Range:	0 - 1000
	Units:	cmH2O/L/Sec
	Label:	Rrs
	AVEA GUI/Membrane:	Rrs
<b>10.6.1.40</b>	<b>ID:</b>	<b>MntrRSBIndex</b>
	Description:	Rapid Shallow Breathing Index. Ratio of Rate to Tidal Volume
	Type:	WORD
	Range:	0 - 500
	Units:	b2/min/L
	Label:	f/Vt
	AVEA GUI/Membrane:	f/Vt
<b>10.6.1.41</b>	<b>ID:</b>	<b>MntrSpO2</b>
	Description:	Patient's SpO2 as measured and reported by the pulse oximeter.
	Type:	WORD
	Range:	1 – 100
	Resolution:	1

Units: %  
 Label: SpO2  
 AVEA GUI/Membrane: SpO2

- 10.6.1.42 ID: MntrSpO2PerfusionIndex**  
 Description: Percentage of pulsatile signal to non-pulsatile signal (pulse strength).  
 Type: WORD  
 Scale: 2  
 Range: 2 – 2000  
 Units: %  
 Label: P.I.  
 AVEA GUI/Membrane: P.I.
- 10.6.1.43 ID: MntrSpO2PulseRate**  
 Description: The patient's pulse rate as measured and reported by the pulse oximeter.  
 Type: WORD  
 Range: 25 – 240  
 Units: bpm  
 Label: Pulse Rate  
 AVEA GUI/Membrane: Pulse Rate
- 10.6.1.44 ID: MntrTe**  
 Description: Exhalation Time. Measured time duration of the expiratory phase of a breath cycle.  
 Type: WORD  
 Scale: 2  
 Range: 0 - 9999  
 Units: sec  
 Label: Te  
 AVEA GUI/Membrane: Te
- 10.6.1.45 ID: MntrTi**  
 Description: Inspiratory time. Measured time duration of the inspiratory phase of a breath cycle.  
 Type: WORD  
 Scale: 2  
 Range: 0 - 9999  
 Units: sec  
 Label: Ti  
 AVEA GUI/Membrane: Ti
- 10.6.1.46 ID: MntrVdel**  
 Description: Delivered Volume. Total volume of gas delivered to the patient circuit by the ventilator flow control valve.

Type: INT  
 Scale (Adult): 2  
 Scale (Pediatric): 0  
 Scale (Neo): 1  
 Range (Adult): 0 – 400  
 Range (Pediatric): 0 – 1999  
 Range (Neo): 0 – 9999  
 Units (Adult): L  
 Units (Pediatric/Neo): mL  
 Label: Vdel  
 AVEA GUI/Membrane: Vdel

- 10.6.1.47 ID: MntrVentilationAnatomicalDeadSpace**
- Description: The patient's anatomical dead space, measured on each breath, and averaged over the "VCO2 Average" time interval.
- Type: WORD  
 Scale: 1  
 Range: 0 – 9990  
 Units: ml  
 Label: Vd ana  
 AVEA GUI/Membrane: Vd ana
- 10.6.1.48 ID: MntrVentilationAnatomicalDeadSpaceVtRatio**
- Description: The patient's airway dead space to tidal volume ratio, measured on each breath, and averaged over the "VCO2 Average" time interval.
- Type: WORD  
 Range: 0 – 99  
 Units: %  
 Label: Vd ana / Vt  
 AVEA GUI/Membrane: Vd ana / Vt
- 10.6.1.49 ID: MntrVentilationCO2**
- Description: The patient's exhaled minute volume of CO2, calculated over the "VCO2 Average" interval.
- Type: WORD  
 Scale: 1  
 Range: 0 – 9990  
 Units: ml/min  
 Label: VCO2  
 AVEA GUI/Membrane: VCO2

<b>10.6.1.50</b>	<b>ID:</b>	<b>MntrVeSpon</b>
	Description:	Spontaneous minute volume. Volume of gas exhaled over the period of one minute from patient controlled breath types.
	Type:	WORD
	Scale (Adult/Ped):	1
	Scale (Neo):	2
	Range (Adult/Ped):	0 - 999
	Range (Neo):	0 - 9990
	Units:	L
	Label:	Spon Ve
	AVEA GUI/Membrane:	Spon Ve
<b>10.6.1.51</b>	<b>ID:</b>	<b>MntrVeSponNorm</b>
	Description:	Normalized Spontaneous Minute Volume. Ratio of Spon Ve to Patient Weight.
	Type:	WORD
	Range:	0 - 999
	Units:	mL/kg
	Label:	Spon Ve/kg
	AVEA GUI/Membrane:	Spon Ve/kg
<b>10.6.1.52</b>	<b>ID:</b>	<b>MntrVeTotal</b>
	Description:	Total Minute Volume. Volume of gas exhaled over the period of one minute from all breath types.
	Type:	WORD
	Scale (Adult/Ped):	1
	Scale (Neo):	2
	Range (Adult/Ped):	0 - 999
	Range (Neo):	0 - 9990
	Units:	L
	Label:	Total Ve
	AVEA GUI/Membrane:	Total Ve
<b>10.6.1.53</b>	<b>ID:</b>	<b>MntrVeTotalNorm</b>
	Description:	Normalized Total Minute Volume. Ratio of Total Ve to Patient Weight.
	Type:	WORD
	Range:	0 - 999
	Units:	mL/kg
	Label:	Total Ve/kg
	AVEA GUI/Membrane:	Total Ve/kg
<b>10.6.1.54</b>	<b>ID:</b>	<b>MntrVtCO2</b>
	Description:	Physiological Dead Space / Tidal Volume Ratio.
	Type:	WORD

Scale: 1  
 Range: 0 – 2990  
 Units: ml  
 Label: VtCO2  
 AVEA GUI/Membrane: VtCO2

- 10.6.1.55 ID: MntrVte**  
 Description: Exhaled Tidal Volume. Volume of gas exhaled in one breath. cycle.  
 Type: INT  
 Scale (Adult): 2  
 Scale (Pediatric): 0  
 Scale (Neo): 1  
 Range (Adult): 0 – 400  
 Range (Pediatric): 0 – 1999  
 Range (Neo): 0 – 9999  
 Units (Adult): L  
 Units (Pediatric/Neo): mL  
 Label: Vte  
 AVEA GUI/Membrane: Vte
- 10.6.1.56 ID: MntrVteMand**  
 Description: Mandatory Exhaled Tidal Volume. Volume of gas exhaled in one machine controlled breath type cycle.  
 Type: INT  
 Scale (Adult): 2  
 Scale (Pediatric): 0  
 Scale (Neo): 1  
 Range (Adult): 0 – 400  
 Range (Pediatric): 0 – 1999  
 Range (Neo): 0 – 9999  
 Units (Adult): L  
 Units (Pediatric/Neo): mL  
 Label: Mand Vte  
 AVEA GUI/Membrane: Mand Vte
- 10.6.1.57 ID: MntrVteMandNorm**  
 Description: Normalized Mandatory Exhaled Tidal Volume. Ratio of Mand Vte to Patient Weight.  
 Type: WORD  
 Scale: 1  
 Range: 0 - 300  
 Units: mL/kg  
 Label: Mand Vte/kg  
 AVEA GUI/Membrane: Mand Vte/kg

<b>10.6.1.58</b>	<b>ID:</b>	<b>MntrVteNorm</b>
	Description:	Normalized Exhaled Tidal Volume. Ratio of Vte to Patient Weight.
	Type:	WORD
	Scale:	1
	Range:	0 - 300
	Units:	mL/kg
	Label:	Vte/kg
	AVEA GUI/Membrane:	Vte/kg
<b>10.6.1.59</b>	<b>ID:</b>	<b>MntrVteSpon</b>
	Description:	Spontaneous Exhaled Tidal Volume. Volume of gas exhaled in one patient controlled breath type cycle.
	Type:	INT
	Scale (Adult):	2
	Scale (Pediatric):	0
	Scale (Neo):	1
	Range (Adult):	0 – 400
	Range (Pediatric):	0 – 1999
	Range (Neo):	0 – 9999
	Units (Adult):	L
	Units (Pediatric/Neo):	mL
	Label:	Spon Vte
	AVEA GUI/Membrane:	Spon Vte
<b>10.6.1.60</b>	<b>ID:</b>	<b>MntrVteSponNorm</b>
	Description:	Normalized Spontaneous Exhaled Tidal Volume. Ratio of Spon Vte to Patient Weight.
	Type:	WORD
	Scale:	1
	Range:	0 - 300
	Units:	mL/kg
	Label:	Spon Vte/kg
	AVEA GUI/Membrane:	Spon Vte/kg
<b>10.6.1.61</b>	<b>ID:</b>	<b>MntrVti</b>
	Description:	Inspired Tidal Volume. Volume of gas inspired in one breath cycle.
	Type:	INT
	Scale (Adult):	2
	Scale (Pediatric):	0
	Scale (Neo):	1
	Range (Adult):	0 – 400
	Range (Pediatric):	0 – 1999
	Range (Neo):	0 – 9999
	Units (Adult):	L

Units (Pediatric/Neo): mL  
 Label: Vti  
 AVEA GUI/Membrane: Vti

- 10.6.1.62 ID: MntrVtiNorm**  
 Description: Normalized Inspired Tidal Volume. Ratio of Vti to Patient Weight.  
 Type: WORD  
 Scale: 1  
 Range: 0 - 300  
 Units: mL/kg  
 Label: Vti/kg  
 AVEA GUI/Membrane: Vti/kg
- 10.6.1.63 ID: MntrWoblImposed**  
 Description: Imposed Work of Breathing. The additional energy imposed by the breathing apparatus during a patient's spontaneous breath.  
 Type: WORD  
 Scale: 2  
 Range: 0 - 2000  
 Units: joules/L  
 Label: WOBi  
 AVEA GUI/Membrane: WOBi
- 10.6.1.64 ID: MntrWobPatient**  
 Description: Patient Work of Breathing. Energy expended by the patient (lung and chest wall) by spontaneous breathing.  
 Type: WORD  
 Scale: 2  
 Range: 0 - 2000  
 Units: joules/L  
 Label: WOBp  
 AVEA GUI/Membrane: WOBp
- 10.6.1.65 ID: MntrWobVent**  
 Description: Ventilator Work of Breathing. Energy required to deliver the tidal volume to the patient.  
 Type: WORD  
 Scale: 2  
 Range: 0 - 2000  
 Units: joules/L  
 Label: WOBv  
 AVEA GUI/Membrane: WOBv

## 10.6.2 VELA Ventilator – Monitor Class

<b>10.6.2.1</b>	<b>ID:</b>	<b>MntrEndTidalCO2</b>
	Description:	Patient's peak expired CO2 level as measured and reported by the CO2 analyzer.
	Type:	WORD
	Scale:	1 (mmHg); 2 (kPa)
	Range:	0 – 1500 (mmHg); 0 – 2000 (kPa)
	Units:	mmHg; kPa
	Label:	EtCO2
	VELA GUI/Membrane:	EtCO2
<b>10.6.2.2</b>	<b>ID:</b>	<b>MntrFiO2</b>
	Description:	Fractional (percent) O <sub>2</sub> (measured in inspiratory flow stream) delivered to the patient.
	Type:	WORD
	Range:	0 - 100
	Units:	%
	Label:	FiO2
	VELA GUI/Membrane:	FiO2
<b>10.6.2.3</b>	<b>ID:</b>	<b>MntrIE</b>
	Description:	Inspiratory/Expiratory Ratio. Ratio of time spent in the inspiratory phase of a breath cycle to the time spent in the expiratory phase.
	Type:	WORD
	Scale:	1
	Range:	-999 - 999
	Label:	I:E
	VELA GUI/Membrane:	I:E
<b>10.6.2.4</b>	<b>ID:</b>	<b>MntrPeep</b>
	Description:	Positive End Expiratory Pressure. Airway pressure at the end of the expiratory phase of a breath cycle.
	Type:	WORD
	Range:	0 - 99
	Units:	cmH2O
	Label:	PEEP
	VELA GUI/Membrane:	PEEP
<b>10.6.2.5</b>	<b>ID:</b>	<b>MntrPmean</b>
	Description:	Mean Airway Pressure. Average airway pressure measured throughout an entire breath cycle.
	Type:	WORD
	Range:	0 - 99
	Units:	cmH2O



	Label:	Pmean
	VELA GUI/Membrane:	Pmean
<b>10.6.2.6</b>	<b>ID:</b>	<b>MntrPO2</b>
	Description:	Oxygen Supply Pressure. Pressure measured at the oxygen inlet of the ventilator.
	Type:	WORD
	Range:	0 – 85
	Units:	psig
	Label:	O2 Inlet
	VELA GUI/Membrane:	O2 Inlet
<b>10.6.2.7</b>	<b>ID:</b>	<b>MntrPpeak</b>
	Description:	Peak Airway Pressure. Highest airway pressure measured throughout an entire breath cycle.
	Type:	WORD
	Range:	0 - 140
	Units:	cmH2O
	Label:	Ppeak
	VELA GUI/Membrane:	Ppeak
<b>10.6.2.8</b>	<b>ID:</b>	<b>MntrPresBaro</b>
	Description:	Barometric Pressure
	Type:	WORD
	Scale:	0 (mmHg); 1 (kPa)
	Range:	545 – 760 (mmHg); 727 – 1013 (kPa)
	Units:	mmHg, kPa
	Label:	Pbaro
	VELA GUI/Membrane:	Pbaro
<b>10.6.2.9</b>	<b>ID:</b>	<b>MntrRate</b>
	Description:	Breath Rate. Number of inspirations per minute of all breath types.
	Type:	WORD
	Range:	0 - 250
	Units:	bpm
	Label:	Rate
	VELA GUI/Membrane:	Rate
<b>10.6.2.10</b>	<b>ID:</b>	<b>MntrRateSpon</b>
	Description:	Spontaneous breath rate. Number of inspirations per minute of patient controlled breath types.
	Type:	WORD
	Range:	0 - 250
	Units:	bpm
	Label:	Spon Rate

VELA GUI/Membrane: Spon Rate

- 10.6.2.11 ID: MntrTe**  
 Description: Exhalation Time. Measured time duration of the expiratory phase of a breath cycle.  
 Type: WORD  
 Scale: 2  
 Range: 0 - 999  
 Units: Sec  
 Label: Te  
 VELA GUI/Membrane: Te
- 10.6.2.12 ID: MntrTi**  
 Description: Inspiratory time. Measured time duration of the inspiratory phase of a breath cycle.  
 Type: WORD  
 Scale: 2  
 Range: 0 - 999  
 Units: Sec  
 Label: Ti  
 VELA GUI/Membrane: Ti
- 10.6.2.13 ID: MntrVeSpon**  
 Description: Spontaneous minute volume. Volume of gas exhaled over the period of one minute from patient controlled breath types.  
 Type: WORD  
 Scale: 1  
 Range: 0 - 999  
 Units: L  
 Label: Spon Ve  
 VELA GUI/Membrane: Spon Ve
- 10.6.2.14 ID: MntrVeTotal**  
 Description: Total Minute Volume. Volume of gas exhaled over the period of one minute from all breath types.  
 Type: WORD  
 Scale: 1  
 Range: 0 - 999  
 Units: L  
 Label: Ve  
 VELA GUI/Membrane: Ve
- 10.6.2.15 ID: MntrVte**  
 Description: Exhaled Tidal Volume. Volume of gas exhaled in one breath cycle.

Type: INT  
 Range: 0 - 4000  
 Units: ml  
 Label: Vte  
 VELA GUI/Membrane: Vte

**10.6.2.16 ID: MntrVteMand**

Description: Mandatory Exhaled Tidal Volume. Volume of gas exhaled in one breath cycle for machine controlled breath types.

Type: INT  
 Range: 0 - 4000  
 Units: ml  
 Label: Mand Vt  
 VELA GUI/Membrane: Mand Vt

**10.6.2.17 ID: MntrVteSpon**

Description: Spontaneous Exhaled Tidal Volume. Volume of gas exhaled in one breath cycle for patient controlled breath types.

Type: INT  
 Range: 0 - 4000  
 Units: ml  
 Label: Spon Vt  
 VELA GUI/Membrane: Spon Vt

**10.6.2.18 ID: MntrVti**

Description: Inspired Tidal Volume. Volume of gas inspired in one breath cycle.

Type: INT  
 Range: 0 - 4000  
 Units: ml  
 Label: Vti  
 VELA GUI/Membrane: Vti

## **10.7 Alarm Class**

### **10.7.1 AVEA Ventilator – Alarm Class**

**10.7.1.1 ID: AlarmActive**

Description: Active/Inactive state of alarms in general. Active if any alarm is asserted.

Type: BOOL  
 Label:  
 AVEA GUI/Membrane:

<b>10.7.1.2</b>	<b>ID:</b>	<b>AlarmActivePriority</b>
	Description:	Highest priority across all active alarm conditions.
	Type:	ENUM
	Enum value = label:	1=HIGH; 2=MED; 3=LOW; 4=VISUAL ALERT
	Label:	
	AVEA GUI/Membrane:	
<b>10.7.1.3</b>	<b>ID:</b>	<b>AlarmApnea</b>
	Description:	Ventilator does not detect a breath initiation within the Apnea Interval time of the previous breath initiation.
	Type:	BOOL
	Level:	HIGH
	Label:	APNEA INTERVAL
	AVEA GUI/Membrane:	APNEA INTERVAL
<b>10.7.1.4</b>	<b>ID:</b>	<b>AlarmCO2CheckAirwayAdapter</b>
	Description:	CO2 device reported the CO2 Airway Adapter needs to be checked.
	Type:	BOOL
	Level:	MED
	Label:	CO2 Check Adapter
	AVEA GUI/Membrane:	CO2 Check Adapter
<b>10.7.1.5</b>	<b>ID:</b>	<b>AlarmCircDisc</b>
	Description:	Patient circuit disconnected from the ventilator or patient.
	Type:	BOOL
	Level:	HIGH
	Label:	CIRCUIT DISCONNECT
	AVEA GUI/Membrane:	CIRCUIT DISCONNECT
<b>10.7.1.6</b>	<b>ID:</b>	<b>AlarmCO2CommunicationError</b>
	Description:	CO2 device reported a communication error.
	Type:	BOOL
	Level:	MED
	Label:	CO2 Comms Error
	AVEA GUI/Membrane:	CO2 Comms Error
<b>10.7.1.7</b>	<b>ID:</b>	<b>AlarmCO2OutOfRange</b>
	Description:	CO2 device reported a CO2 range error condition.
	Type:	BOOL
	Level:	MED
	Label:	CO2 Out Of Range
	AVEA GUI/Membrane:	CO2 Out Of Range
<b>10.7.1.8</b>	<b>ID:</b>	<b>AlarmCO2SensorFault</b>
	Description:	CO2 device reported a fault condition with CO2 sensor.

	Type:	BOOL
	Level:	MED
	Label:	CO2 Sensor Fault
	AVEA GUI/Membrane:	CO2 Sensor Fault
<b>10.7.1.9</b>	<b>ID:</b>	<b>AlarmCO2SensorOverTemp</b>
	Description:	CO2 device reported a fault condition due to temperature.
	Type:	BOOL
	Level:	MED
	Label:	CO2 Sensor Temp
	AVEA GUI/Membrane:	CO2 Sensor Temp
<b>10.7.1.10</b>	<b>ID:</b>	<b>AlarmCO2ZeroRequired</b>
	Description:	CO2 device reported sensor requires to be initialized to zero.
	Type:	BOOL
	Level:	MED
	Label:	CO2 Zero Req'd
	AVEA GUI/Membrane:	CO2 Zero Req'd
<b>10.7.1.11</b>	<b>ID:</b>	<b>AlarmFanFail</b>
	Description:	Internal cooling/enclosure ventilation fan has failed.
	Type:	BOOL
	Level:	LOW
	Label:	FAN FAILURE
	AVEA GUI/Membrane:	FAN FAILURE
<b>10.7.1.12</b>	<b>ID:</b>	<b>AlarmFiO2AutoLimitHigh</b>
	Description:	Measured FiO <sub>2</sub> exceeded the high oxygen concentration limit.
	Type:	BOOL
	Level:	HIGH
	Label:	High Auto FiO <sub>2</sub>
	AVEA GUI/Membrane:	High Auto FiO <sub>2</sub>
<b>10.7.1.13</b>	<b>ID:</b>	<b>AlarmFiO2AutoLimitLow</b>
	Description:	The output Auto FiO <sub>2</sub> Cmd shall not be allowed to be less than the Low Auto FiO <sub>2</sub> Limit. The "Low Auto FiO <sub>2</sub> Limit" will alarm if Auto FiO <sub>2</sub> Cmd $\leq$ preset Low Auto FiO <sub>2</sub> Limit for a period of 60 seconds or greater.
	Type:	BOOL
	Level:	HIGH
	Label:	Low Auto FiO <sub>2</sub>
	AVEA GUI/Membrane:	Low Auto FiO <sub>2</sub>

<b>10.7.1.14</b>	<b>ID:</b>	<b>AlarmFiO2BaselineLimitHigh</b>
	Description:	Alarm if the Baseline FiO <sub>2</sub> $\geq$ High Baseline FiO <sub>2</sub> Alarm setting.
	Type:	BOOL
	Level:	LOW
	Label:	High Base FiO2
	AVEA GUI/Membrane:	High Base FiO2
<b>10.7.1.15</b>	<b>ID:</b>	<b>AlarmFiO2High</b>
	Description:	Measured FiO <sub>2</sub> exceeded the high oxygen concentration limit.
	Type:	BOOL
	Level:	HIGH
	Label:	HIGH FiO2
	AVEA GUI/Membrane:	HIGH FiO2
<b>10.7.1.16</b>	<b>ID:</b>	<b>AlarmFiO2Low</b>
	Description:	Measured FiO <sub>2</sub> dropped below the low oxygen concentration limit.
	Type:	BOOL
	Level:	HIGH
	Label:	LOW FiO2
	AVEA GUI/Membrane:	LOW FiO2
<b>10.7.1.17</b>	<b>ID:</b>	<b>AlarmIlvSlaveDisc</b>
	Description:	Lost detection of slave ventilator during Independent Lung Ventilation.
	Type:	BOOL
	Level:	HIGH
	Label:	ILV DISCONNECT
	AVEA GUI/Membrane:	ILV DISCONNECT
<b>10.7.1.18</b>	<b>ID:</b>	<b>AlarmInop</b>
	Description:	Indicates ventilator is unable to ventilate the patient. Most commonly due to detection of an unrecoverable internal problem or failure.
	Type:	BOOL
	Level:	HIGH
	Label:	VENT INOP
	AVEA GUI/Membrane:	VENT INOP
<b>10.7.1.19</b>	<b>ID:</b>	<b>AlarmEndTidalCO2Invalid</b>
	Description:	Invalid EtCO <sub>2</sub> .
	Type:	BOOL
	Level:	MED
	Label:	Invalid EtCO2

AVEA GUI/Membrane: Invalid EtCO2

**10.7.1.20 ID: AlarmEndTidalCO2High**

Description: Exceeded EtCO2 High Limit  
 Type: BOOL  
 Level: LOW  
 Label: High EtCO2  
 AVEA GUI/Membrane: High EtCO2

**10.7.1.21 ID: AlarmEndTidalCO2Low**

Description: Below EtCO2 Low Limit  
 Type: BOOL  
 Level: LOW  
 Label: Low EtCO2  
 AVEA GUI/Membrane: Low EtCO2

**10.7.1.22 ID: AlarmInvalidGasId**

Description: Missing or defective Gas Type Identifier plug on ventilator rear panel.  
 Type: BOOL  
 Level: MED  
 Label: INVALID GAS ID  
 AVEA GUI/Membrane: INVALID GAS ID

**10.7.1.23 ID: AlarmLimitIE**

Description: Maximum I:E Ratio limit exceeded.  
 Type: BOOL  
 Level: LOW  
 Label: I:E LIMIT  
 AVEA GUI/Membrane: I:E LIMIT

**10.7.1.24 ID: AlarmLimitTi**

Description: Maximum inspiratory time limit exceeded.  
 Type: BOOL  
 Level: LOW  
 Label: MAX INSP TIME  
 AVEA GUI/Membrane: MAX INSP TIME

**10.7.1.25 ID: AlarmLimitVol**

Description: Pressure Control Volume Limit was invoked. The pressure controlled breath was terminated when the volume limit was exceeded without achieving the set pressure.  
 Type: BOOL  
 Level: VISUAL ALERT  
 Label: VOLUME LIMIT

AVEA GUI/Membrane: VOLUME LIMIT

- 10.7.1.26 ID: AlarmLossAir**  
 Description: Ventilator cannot detect a source of air supply. Not connected, insufficient pressure, no functional internal compressor or the compressor output is insufficient to meet demand.  
 Type: BOOL  
 Level: HIGH  
 Label: LOSS OF AIR  
 AVEA GUI/Membrane: LOSS OF AIR
- 10.7.1.27 ID: AlarmLossGas**  
 Description: Ventilator cannot detect any source of gas supply including Air, Oxygen, Heliox, or Internal Compressor. Not connected or insufficient pressure.  
 Type: BOOL  
 Level: HIGH  
 Label: LOSS OF GAS  
 AVEA GUI/Membrane: LOSS OF GAS
- 10.7.1.28 ID: AlarmLossHeliox**  
 Description: Ventilator cannot detect a source of Heliox gas supply. Not connected or insufficient pressure.  
 Type: BOOL  
 Level: HIGH  
 Label: LOSS OF HELIOX  
 AVEA GUI/Membrane: LOSS OF HELIOX
- 10.7.1.29 ID: AlarmLossO2**  
 Description: Ventilator cannot detect a source of oxygen gas supply. Not connected or insufficient pressure.  
 Type: BOOL  
 Level: HIGH  
 Label: LOSS OF O2  
 AVEA GUI/Membrane: LOSS OF O2
- 10.7.1.30 ID: AlarmNcpapHigh**  
 Description: Nasal CPAP high pressure  
 Type: BOOL  
 Level: HIGH  
 Label: HIGH nCPAP PRES  
 AVEA GUI/Membrane: HIGH nCPAP PRES
- 10.7.1.31 ID: AlarmNcpapHighPresLimit**  
 Description: Nasal CPAP high pressure (time) limit



	Type:	BOOL
	Level:	HIGH
	Label:	nCPAP PRES LIMIT
	AVEA GUI/Membrane:	nCPAP PRES LIMIT
<b>10.7.1.32</b>	<b>ID:</b>	<b>AlarmNcpapLow</b>
	Description:	Nasal CPAP low pressure
	Type:	BOOL
	Level:	HIGH
	Label:	LOW nCPAP PRES
	AVEA GUI/Membrane:	LOW nCPAP PRES
<b>10.7.1.33</b>	<b>ID:</b>	<b>AlarmOcclusion</b>
	Description:	Differential pressure across the Inspiratory/Expiratory systems of the ventilator indicates excessive resistance or occlusion.
	Type:	BOOL
	Level:	HIGH
	Label:	CIRCUIT OCCLUSION
	AVEA GUI/Membrane:	CIRCUIT OCCLUSION
<b>10.7.1.34</b>	<b>ID:</b>	<b>AlarmOpenSV</b>
	Description:	Asserts whenever system conditions cause the Safety Valve to open.
	Type:	BOOL
	Level:	HIGH
	Label:	SAFETY VALVE
	AVEA GUI/Membrane:	SAFETY VALVE
<b>10.7.1.35</b>	<b>ID:</b>	<b>AlarmPeepLow</b>
	Description:	Airway pressure has dropped below the low baseline pressure limit.
	Type:	BOOL
	Level:	HIGH
	Label:	LOW PEEP
	AVEA GUI/Membrane:	LOW PEEP
<b>10.7.1.36</b>	<b>ID:</b>	<b>AlarmPpeakHigh</b>
	Description:	Airway pressure has exceeded the high peak pressure limit.
	Type:	BOOL
	Level:	HIGH
	Label:	HIGH Ppeak
	AVEA GUI/Membrane:	HIGH Ppeak

- 10.7.1.37 ID: AlarmPpeakHighExt**  
Description: High P<sub>PEAK</sub> or Occlusion persistent for greater than 5 seconds.  
Type: BOOL  
Level: HIGH  
Label: EXT HIGH Ppeak  
AVEA GUI/Membrane: EXT HIGH Ppeak
- 10.7.1.38 ID: AlarmPpeakLow**  
Description: Airway pressure did not exceed the low peak pressure limit for the previous breath cycle.  
Type: BOOL  
Level: HIGH  
Label: LOW Ppeak  
AVEA GUI/Membrane: LOW Ppeak
- 10.7.1.39 ID: AlarmPulseRateHigh**  
Description: Pulse rate is greater than the High Pulse Rate alarm setting.  
Type: BOOL  
Level: MED  
Label: High Pulse  
AVEA GUI/Membrane: High Pulse
- 10.7.1.40 ID: AlarmPulseRateLow**  
Description: Pulse rate is less than the Low Pulse Rate alarm setting.  
Type: BOOL  
Level: MED  
Label: Low Pulse  
AVEA GUI/Membrane: Low Pulse
- 10.7.1.41 ID: AlarmPwrAcLoss**  
Description: AC power has been removed from the ventilator (i.e. power cord disconnect or loss of supply power).  
Type: BOOL  
Level: HIGH  
Label: LOSS OF A/C  
AVEA GUI/Membrane: LOSS OF A/C
- 10.7.1.42 ID: AlarmPwrBattLow**  
Description: Batteries have been depleted to a level that indicates two minutes or less of safe operation.  
Type: BOOL  
Level: HIGH  
Label: LOW BATTERY  
AVEA GUI/Membrane: LOW BATTERY

- 10.7.1.43 ID: AlarmRateHigh**  
Description: Respiratory Rate exceeded the high inspirations per minute limit.  
Type: BOOL  
Level: MED  
Label: HIGH RATE  
AVEA GUI/Membrane: HIGH RATE
- 10.7.1.44 ID: AlarmSilence**  
Description: Active/Inactive state of the capability to locally silence the audible ventilator alarms.  
Type: BOOL  
Label: ALARM SILENCE  
AVEA GUI/Membrane: ALARM SILENCE
- 10.7.1.45 ID: AlarmSpO2Failure**  
Description: High priority alarm is issued when the MS-11 PCB reports a board failure (see Masimo CSD-1086 Rev C).  
Type: BOOL  
Level: HIGH  
Label: Oximeter Failure  
AVEA GUI/Membrane: Oximeter Failure
- 10.7.1.46 ID: AlarmSpO2High**  
Description: Monitored SpO2 is greater than the preset High SpO2 Alarm for more than the SpO2 Alarm Delay period.  
Type: BOOL  
Level: HIGH  
Label: High SpO2  
AVEA GUI/Membrane: High SpO2
- 10.7.1.47 ID: AlarmSpO2Low**  
Description: Monitored SpO2 is less than the preset Low SpO2 Alarm for more than the SpO2 Alarm Delay period.  
Type: BOOL  
Level: HIGH  
Label: Low SpO2  
AVEA GUI/Membrane: Low SpO2
- 10.7.1.48 ID: AlarmSpO2NotConnected**  
Description: High priority alarm is issued when the pulse oximeter is not connected to the ventilator (i.e. RS-232 not connected).  
Type: BOOL  
Level: HIGH  
Label: Oximeter Disconnect

AVEA GUI/Membrane: Oximeter Disconnect

- 10.7.1.49 ID: AlarmSpO2SensorDefective**  
 Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Defective.  
 Type: BOOL  
 Level: HIGH  
 Label: Defective Sensor  
 AVEA GUI/Membrane: Defective Sensor
- 10.7.1.50 ID: AlarmSpO2SensorNotConnected**  
 Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Not Connected.  
 Type: BOOL  
 Level: HIGH  
 Label: No Sensor Connected  
 AVEA GUI/Membrane: No Sensor Connected
- 10.7.1.51 ID: AlarmSpO2SensorUnrecognized**  
 Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Unrecognized.  
 Type: BOOL  
 Level: HIGH  
 Label: Unrecognized Sensor  
 AVEA GUI/Membrane: Unrecognized Sensor
- 10.7.1.52 ID: AlarmSpO2SignalAmbientLight**  
 Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Ambient Light.  
 Type: BOOL  
 Level: LOW  
 Label: Ambient Light  
 AVEA GUI/Membrane: Ambient Light
- 10.7.1.53 ID: AlarmSpO2SignalAmbientLightExt**  
 Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Ambient Light for an extended period.  
 Type: BOOL  
 Level: HIGH  
 Label: Ambient Light  
 AVEA GUI/Membrane: Ambient Light

- 10.7.1.54 ID: AlarmSpO2SignalInterference**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Interference.  
Type: BOOL  
Level: LOW  
Label: Interference  
AVEA GUI/Membrane: Interference
- 10.7.1.55 ID: AlarmSpO2SignalInterferenceExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Interference for an extended period.  
Type: BOOL  
Level: HIGH  
Label: Interference  
AVEA GUI/Membrane: Interference
- 10.7.1.56 ID: AlarmSpO2SignalLowPerfusion**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Perfusion is Low.  
Type: BOOL  
Level: LOW  
Label: Low Perfusion  
AVEA GUI/Membrane: Low Perfusion
- 10.7.1.57 ID: AlarmSpO2SignalLowPerfusionExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Perfusion is Low for an extended period.  
Type: BOOL  
Level: HIGH  
Label: Low Perfusion  
AVEA GUI/Membrane: Low Perfusion
- 10.7.1.58 ID: AlarmSpO2SignalLowSIQ**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal is Low.  
Type: BOOL  
Level: LOW  
Label: Low SIQ  
AVEA GUI/Membrane: Low SIQ

- 10.7.1.59 ID: AlarmSpO2SignalLowSIQExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal is Low for an extended period.  
Type: BOOL  
Level: HIGH  
Label: Low SIQ  
AVEA GUI/Membrane: Low SIQ
- 10.7.1.60 ID: AlarmSpO2SignalPulseSearch**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Pulse Search is Low.  
Type: BOOL  
Level: LOW  
Label: Pulse Search  
AVEA GUI/Membrane: Pulse Search
- 10.7.1.61 ID: AlarmSpO2SignalPulseSearchExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Pulse Search is Low for an extended period.  
Type: BOOL  
Level: HIGH  
Label: Pulse Search  
AVEA GUI/Membrane: Pulse Search
- 10.7.1.62 ID: AlarmSpO2SignalSensorOff**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Sensor is off.  
Type: BOOL  
Level: LOW  
Label: Sensor Off  
AVEA GUI/Membrane: Sensor Off
- 10.7.1.63 ID: AlarmSpO2SignalSensorOffExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Sensor is off for an extended period.  
Type: BOOL  
Level: HIGH  
Label: Sensor Off  
AVEA GUI/Membrane: Sensor Off

<b>10.7.1.64</b>	<b>ID:</b>	<b>AlarmTest</b>
	Description:	Operator demand to test audible and visual indications of the Alarm system (e.g. during Alarm Loudness Increase/Decrease).
	Type:	BOOL
	Level:	HIGH
	Label:	ALARM TEST
	AVEA GUI/Membrane:	ALARM TEST
<b>10.7.1.65</b>	<b>ID:</b>	<b>AlarmVeHigh</b>
	Description:	Minute Volume exceeded the high minute volume limit.
	Type:	BOOL
	Level:	MED
	Label:	HIGH Ve
	AVEA GUI/Membrane:	HIGH Ve
<b>10.7.1.66</b>	<b>ID:</b>	<b>AlarmVeLow</b>
	Description:	Exhaled tidal volume was less than the low tidal volume limit for the number of consecutive breaths set by the Low Vte Sensitivity setting.
	Type:	BOOL
	Level:	HIGH
	Label:	LOW Ve
	AVEA GUI/Membrane:	LOW Ve
<b>10.7.1.67</b>	<b>ID:</b>	<b>AlarmVteLow</b>
	Description:	Exhaled tidal volume was less than the low tidal volume limit for the number of consecutive breaths set by the Low Vte Sensitivity setting.
	Type:	BOOL
	Level:	HIGH
	Label:	LOW Vte
	AVEA GUI/Membrane:	LOW Vte
<b>10.7.1.68</b>	<b>ID:</b>	<b>AlarmVtHigh</b>
	Description:	Exhaled tidal volume was higher than the high tidal volume limit for the previous breath cycle.
	Type:	BOOL
	Level:	LOW
	Label:	HIGH Vte
	AVEA GUI/Membrane:	HIGH Vte
<b>10.7.1.69</b>	<b>ID:</b>	<b>AlarmHistApnea</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmApnea for description.

Type: BOOL  
Level: HIGH  
Label: APNEA INTERVAL  
AVEA GUI/Membrane: APNEA INTERVAL

- 10.7.1.70 ID: AlarmHist CO2CheckAirwayAdapter**  
Description: CO2 device reported a CO2 Airway Adapter needs to be checked (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Check Adapter  
AVEA GUI/Membrane: CO2 Check Adapter
- 10.7.1.71 ID: AlarmHistCO2CommunicationError**  
Description: CO2 device reported a communication error (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Comms Error  
AVEA GUI/Membrane: CO2 Comms Error
- 10.7.1.72 ID: AlarmHistCO2OutOfRange**  
Description: CO2 device reported a CO2 range error condition (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Out Of Range  
AVEA GUI/Membrane: CO2 Out Of Range
- 10.7.1.73 ID: AlarmHistCO2SensorFault**  
Description: CO2 device reported a fault condition with CO2 sensor (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Sensor Fault  
AVEA GUI/Membrane: CO2 Sensor Fault
- 10.7.1.74 ID: AlarmHistCO2SensorOverTemp**  
Description: CO2 device reported a fault condition due to temperature (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Sensor Temp  
AVEA GUI/Membrane: CO2 Sensor Temp



- 10.7.1.75 ID: AlarmHistCO2ZeroRequired**  
Description: CO2 device reported sensor requires to be initialized to zero (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Zero Req'd  
AVEA GUI/Membrane: CO2 Zero Req'd
- 10.7.1.76 ID: AlarmHistEndTidalCO2High**  
Description: Exceeded EtCO2 High Limit (historical / not cleared).  
Type: BOOL  
Level: LOW  
Label: High EtCO2  
AVEA GUI/Membrane: High EtCO2
- 10.7.1.77 ID: AlarmHistEndTidalCO2Low**  
Description: Below EtCO2 Low Limit(historical / not cleared).  
Type: BOOL  
Level: LOW  
Label: Low EtCO2  
AVEA GUI/Membrane: Low EtCO2
- 10.7.1.78 ID: AlarmHistCircDisc**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmCircDisc description above.  
Type: BOOL  
Level: HIGH  
Label: CIRCUIT DISCONNECT  
AVEA GUI/Membrane: CIRCUIT DISCONNECT
- 10.7.1.79 ID: AlarmHistEndTidalCO2Invalid**  
Description: Invalid EtCO2 (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: Invalid EtCO2  
AVEA GUI/Membrane: Invalid EtCO2
- 10.7.1.80 ID: AlarmHistFanFail**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmFanFail description above.  
Type: BOOL  
Level: LOW  
Label: FAN FAILURE  
AVEA GUI/Membrane: FAN FAILURE

<b>10.7.1.81</b>	<b>ID:</b>	<b>AlarmHistFiO2AutoLimitHigh</b>
	Description:	Auto FiO2 Cmd is greater than or equal to the preset High Auto FiO2 Limit for a period of 60 seconds or more (historical / not cleared).
	Type:	BOOL
	Level:	HIGH
	Label:	High Auto FiO2
	AVEA GUI/Membrane:	High Auto FiO2
<b>10.7.1.82</b>	<b>ID:</b>	<b>AlarmHistFiO2AutoLimitLow</b>
	Description:	Auto FiO2 Cmd is less than or equal to the preset Low Auto FiO2 Limit for a period of 60 seconds or more (historical / not cleared).
	Type:	BOOL
	Level:	HIGH
	Label:	Low Auto FiO2
	AVEA GUI/Membrane:	Low Auto FiO2
<b>10.7.1.83</b>	<b>ID:</b>	<b>AlarmHistFiO2BaselineLimitHigh</b>
	Description:	Baseline FiO2 is greater than or equal to the High Baseline FiO2 Alarm setting sensor (historical / not cleared).
	Type:	BOOL
	Level:	LOW
	Label:	High Base FiO2
	AVEA GUI/Membrane:	High Base FiO2
<b>10.7.1.84</b>	<b>ID:</b>	<b>AlarmHistFiO2High</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmFiO2High description above.
	Type:	BOOL
	Level:	HIGH
	Label:	HIGH FiO2
	AVEA GUI/Membrane:	HIGH FiO2
<b>10.7.1.85</b>	<b>ID:</b>	<b>AlarmHistFiO2Low</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmFiO2Low description above.
	Type:	BOOL
	Level:	HIGH
	Label:	LOW FiO2
	AVEA GUI/Membrane:	LOW FiO2

<b>10.7.1.86</b>	<b>ID:</b>	<b>AlarmHistIlvSlaveDisc</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmIlvsSlaveDisc description above.
	Type:	BOOL
	Level:	HIGH
	Label:	ILV DISCONNECT
	AVEA GUI/Membrane:	ILV DISCONNECT
<b>10.7.1.87</b>	<b>ID:</b>	<b>AlarmHistInop</b>
	Description:	Indication that alarm asserted in the past, is no longer active, and has not been reset. See AlarmInop description above.
	Type:	BOOL
	Level:	HIGH
	Label:	VENT INOP
	AVEA GUI/Membrane:	VENT INOP
<b>10.7.1.88</b>	<b>ID:</b>	<b>AlarmHistInvalidGasId</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmInvalidGasId description above.
	Type:	BOOL
	Level:	MED
	Label:	INVALID GAS ID
	AVEA GUI/Membrane:	INVALID GAS ID
<b>10.7.1.89</b>	<b>ID:</b>	<b>AlarmHistLimitIE</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmLimitIE description above.
	Type:	BOOL
	Level:	LOW
	Label:	I:E LIMIT
	AVEA GUI/Membrane:	I:E LIMIT
<b>10.7.1.90</b>	<b>ID:</b>	<b>AlarmHistLimitTi</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmLimitTi description above.
	Type:	BOOL
	Level:	LOW
	Label:	MAX INSP TIME
	AVEA GUI/Membrane:	MAX INSP TIME

<b>10.7.1.91</b>	<b>ID:</b>	<b>AlarmHistLimitVol</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmLimitVol description above.
	Type:	BOOL
	Level:	VISUAL ALERT
	Label:	VOL LIMIT
	AVEA GUI/Membrane:	VOL LIMIT
<b>10.7.1.92</b>	<b>ID:</b>	<b>AlarmHistLossAir</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmLossAir description above.
	Type:	BOOL
	Level:	HIGH
	Label:	LOSS OF AIR
	AVEA GUI/Membrane:	LOSS OF AIR
<b>10.7.1.93</b>	<b>ID:</b>	<b>AlarmHistLossGas</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmLossGas description above.
	Type:	BOOL
	Level:	HIGH
	Label:	LOSS OF GAS
	AVEA GUI/Membrane:	LOSS OF GAS
<b>10.7.1.94</b>	<b>ID:</b>	<b>AlarmHistLossHeliox</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmLossHeliox description above.
	Type:	BOOL
	Level:	HIGH
	Label:	LOSS OF HELIOX
	AVEA GUI/Membrane:	LOSS OF HELIOX
<b>10.7.1.95</b>	<b>ID:</b>	<b>AlarmHistLossO2</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmLossO2 description above.
	Type:	BOOL
	Level:	HIGH
	Label:	LOSS OF O2
	AVEA GUI/Membrane:	LOSS OF O2

<b>10.7.1.96</b>	<b>ID:</b>	<b>AlarmHistNcpapHigh</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmNcpapHigh description above.
	Type:	BOOL
	Level:	HIGH
	Label:	HIGH nCPAP PRES
	AVEA GUI/Membrane:	HIGH nCPAP PRES
<b>10.7.1.97</b>	<b>ID:</b>	<b>AlarmHistNcpapHighPresLimit</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmNcpapHighPresLimit description above.
	Type:	BOOL
	Level:	HIGH
	Label:	nCPAP PRES LIMIT
	AVEA GUI/Membrane:	nCPAP PRES LIMIT
<b>10.7.1.98</b>	<b>ID:</b>	<b>AlarmHistNcpapLow</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmNcpapLow description above.
	Type:	BOOL
	Level:	HIGH
	Label:	LOW nCPAP PRES
	AVEA GUI/Membrane:	LOW nCPAP PRES
<b>10.7.1.99</b>	<b>ID:</b>	<b>AlarmHistOcclusion</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmOcclusion description above.
	Type:	BOOL
	Level:	HIGH
	Label:	CIRCUIT OCCLUSION
	AVEA GUI/Membrane:	CIRCUIT OCCLUSION
<b>10.7.1.100</b>	<b>ID:</b>	<b>AlarmHistOpenSV</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmOpenSV description above.
	Type:	BOOL
	Level:	HIGH
	Label:	SAFETY VALVE
	AVEA GUI/Membrane:	SAFETY VALVE

- 10.7.1.101 ID: AlarmHistPeepLow**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPeepLow description above.  
Type: BOOL  
Level: HIGH  
Label: LOW PEEP  
AVEA GUI/Membrane: LOW PEEP
- 10.7.1.102 ID: AlarmHistPpeakHigh**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPpeakHigh description above.  
Type: BOOL  
Level: HIGH  
Label: HIGH Ppeak  
AVEA GUI/Membrane: HIGH Ppeak
- 10.7.1.103 ID: AlarmHistPpeakHighExt**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPpeakHighExt description above.  
Type: BOOL  
Level: HIGH  
Label: EXT HIGH Ppeak  
AVEA GUI/Membrane: EXT HIGH Ppeak
- 10.7.1.104 ID: AlarmHistPpeakLow**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPpeakLow description above.  
Type: BOOL  
Level: HIGH  
Label: LOW Ppeak  
AVEA GUI/Membrane: LOW Ppeak
- 10.7.1.105 ID: AlarmHistPulseRateHigh**  
Description: Pulse rate is greater than the High Pulse Rate alarm setting (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: High Pulse  
AVEA GUI/Membrane: High Pulse

- 10.7.1.106 ID: AlarmHistPulseRateLow**  
Description: Pulse rate is less than the Low Pulse Rate alarm setting (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: Low Pulse  
AVEA GUI/Membrane: Low Pulse
- 10.7.1.107 ID: AlarmHistPwrAcLoss**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPwrAcLoss description above.  
Type: BOOL  
Level: HIGH  
Label: LOSS OF A/C  
AVEA GUI/Membrane: LOSS OF A/C
- 10.7.1.108 ID: AlarmHistPwrBattLow**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPwrBattLow description above.  
Type: BOOL  
Level: HIGH  
Label: LOW BATTERY  
AVEA GUI/Membrane: LOW BATTERY
- 10.7.1.109 ID: AlarmHistRateHigh**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmRateHigh description above.  
Type: BOOL  
Level: MED  
Label: HIGH RATE  
AVEA GUI/Membrane: HIGH RATE
- 10.7.1.110 ID: AlarmHistSpO2Failure**  
Description: High priority alarm is issued when the MS-11 PCB reports a board failure (see Masimo CSD-1086 Rev C) - (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Oximeter Failure  
AVEA GUI/Membrane: Oximeter Failure

- 10.7.1.111 ID: AlarmHistSpO2High**  
Description: Monitored SpO2 is greater than the preset High SpO2 Alarm for more than the SpO2 Alarm Delay period (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: High SpO2  
AVEA GUI/Membrane: High SpO2
- 10.7.1.112 ID: AlarmHistSpO2Low**  
Description: Monitored SpO2 is less than the preset Low SpO2 Alarm for more than the SpO2 Alarm Delay period (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Low SpO2  
AVEA GUI/Membrane: Low SpO2
- 10.7.1.113 ID: AlarmHistSpO2NotConnected**  
Description: High priority alarm is issued when the pulse oximeter is not connected to the ventilator (i.e. RS-232 not connected) - (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Oximeter Disconnect  
AVEA GUI/Membrane: Oximeter Disconnect
- 10.7.1.114 ID: AlarmHistSpO2SensorDefective**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Defective (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Defective Sensor  
AVEA GUI/Membrane: Defective Sensor
- 10.7.1.115 ID: AlarmHistSpO2SensorNotConnected**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Not Connected (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: No Sensor Connected  
AVEA GUI/Membrane: No Sensor Connected



- 10.7.1.116 ID: AlarmHistSpO2SensorUnrecognized**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Sensor is Unrecognized (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Unrecognized Sensor  
AVEA GUI/Membrane: Unrecognized Sensor
- 10.7.1.117 ID: AlarmHistSpO2SignalAmbientLight**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Ambient Light (historical / not cleared).  
Type: BOOL  
Level: LOW  
Label: Ambient Light  
AVEA GUI/Membrane: Ambient Light
- 10.7.1.118 ID: AlarmHistSpO2SignalAmbientLightExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Ambient Light for an extended period (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Ambient Light  
AVEA GUI/Membrane: Ambient Light
- 10.7.1.119 ID: AlarmHistSpO2SignalInterference**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Interference (historical / not cleared).  
Type: BOOL  
Level: LOW  
Label: Interference  
AVEA GUI/Membrane: Interference
- 10.7.1.120 ID: AlarmHistSpO2SignalInterferenceExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Interference for an extended period (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Interference  
AVEA GUI/Membrane: Interference

- 10.7.1.121 ID: AlarmHistSpO2SignalLowPerfusion**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Perfusion is Low (historical / not cleared).  
Type: BOOL  
Level: LOW  
Label: Low Perfusion  
AVEA GUI/Membrane: Low Perfusion
- 10.7.1.122 ID: AlarmHistSpO2SignalLowPerfusionExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Perfusion is Low for an extended period (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Low Perfusion  
AVEA GUI/Membrane: Low Perfusion
- 10.7.1.123 ID: AlarmHistSpO2SignalLowSIQ**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal is Low (historical / not cleared).  
Type: BOOL  
Level: LOW  
Label: Low SIQ  
AVEA GUI/Membrane: Low SIQ
- 10.7.1.124 ID: AlarmHistSpO2SignalLowSIQExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal is Low for an extended period (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Low SIQ  
AVEA GUI/Membrane: Low SIQ
- 10.7.1.125 ID: AlarmHistSpO2SignalPulseSearch**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Pulse Search is Low (historical / not cleared).  
Type: BOOL  
Level: LOW  
Label: Pulse Search  
AVEA GUI/Membrane: Pulse Search

- 10.7.1.126 ID: AlarmHistSpO2SignalPulseSearchExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Pulse Search is Low for an extended period (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Pulse Search  
AVEA GUI/Membrane: Pulse Search
- 10.7.1.127 ID: AlarmHistSpO2SignalSensorOff**  
Description: Low priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Sensor is off (historical / not cleared).  
Type: BOOL  
Level: LOW  
Label: Sensor Off  
AVEA GUI/Membrane: Sensor Off
- 10.7.1.128 ID: AlarmHistSpO2SignalSensorOffExt**  
Description: High priority alarm is issued when the MS-11 PCB (see Masimo CSD-1086 Rev C) detects the SpO2 Signal Sensor is off for an extended period (historical / not cleared).  
Type: BOOL  
Level: HIGH  
Label: Sensor Off  
AVEA GUI/Membrane: Sensor Off
- 10.7.1.129 ID: AlarmHistTest**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmTest description above.  
Type: BOOL  
Level: HIGH  
Label: ALARM TEST  
AVEA GUI/Membrane: ALARM TEST
- 10.7.1.130 ID: AlarmHistVeHigh**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmVeHigh description above.  
Type: BOOL  
Level: MED  
Label: HIGH Ve  
AVEA GUI/Membrane: HIGH Ve

- 10.7.1.131 ID: AlarmHistVeLow**  
 Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmVeLow description above.  
 Type: BOOL  
 Level: HIGH  
 Label: LOW Ve  
 AVEA GUI/Membrane: LOW Ve
- 10.7.1.132 ID: AlarmHistVteLow**  
 Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmVteLow description above.  
 Type: BOOL  
 Level: HIGH  
 Label: LOW Vte  
 AVEA GUI/Membrane: LOW Vte
- 10.7.1.133 ID: AlarmHistVtHigh**  
 Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmVtHigh description above.  
 Type: BOOL  
 Level: LOW  
 Label: HIGH Vte  
 AVEA GUI/Membrane: HIGH Vte

## 10.7.2 VELA Ventilator – Alarm Class

- 10.7.2.1 ID: AlarmActive**  
 Description: Active/Inactive state of alarms in general. Active if any alarm is asserted.  
 Type: BOOL  
 VELA GUI/Membrane:
- 10.7.2.2 ID: AlarmActivePriority**  
 Description: Highest priority across all active alarm conditions.  
 Type: ENUM  
 Enum value = label: 1=HIGH; 2=MED; 3=LOW; 4=ALERT  
 Label:  
 VELA GUI/Membrane:
- 10.7.2.3 ID: AlarmApnea**  
 Description: Ventilator does not detect a breath initiation within the Apnea Interval time of the previous breath initiation.  
 Type: BOOL

	Level:	HIGH
	Label:	APNEA INTERVAL
	VELA GUI/Membrane:	APNEA INTERVAL
<b>10.7.2.4</b>	<b>ID:</b>	<b>AlarmCO2CheckAirwayAdapter</b>
	Description:	CO2 device reported a CO2 Airway Adapter needs to be checked.
	Type:	BOOL
	Level:	MED
	Label:	CO2 Check Adapter
	VELA GUI/Membrane:	CO2 Check Adapter
<b>10.7.2.5</b>	<b>ID:</b>	<b>AlarmCheckEvents</b>
	Description:	Indicates that an unacceptable number of anomalous software detected events have occurred.
	Type:	BOOL
	Level:	MED
	Label:	CHECK EVENTS
	VELA GUI/Membrane:	CHECK EVENTS
<b>10.7.2.6</b>	<b>ID:</b>	<b>AlarmCircDisc</b>
	Description:	Indicates that the patient circuit is disconnected from the ventilator or patient.
	Type:	BOOL
	Level:	HIGH
	Label:	CIRCUIT FAULT
	VELA GUI/Membrane:	CIRCUIT FAULT
<b>10.7.2.7</b>	<b>ID:</b>	<b>AlarmClockBattLow</b>
	Description:	The battery for the real-time clock is low and should be replaced.
	Type:	BOOL
	Level:	MED
	Label:	LOW CLOCK BATTERY
	VELA GUI/Membrane:	LOW CLOCK BATTERY
<b>10.7.2.8</b>	<b>ID:</b>	<b>AlarmCO2CommunicationError</b>
	Description:	CO2 device reported a communication error.
	Type:	BOOL
	Level:	MED
	Label:	CO2 Comms Error
	VELA GUI/Membrane:	CO2 Comms Error
<b>10.7.2.9</b>	<b>ID:</b>	<b>AlarmCO2OutOfRange</b>
	Description:	CO2 device reported a CO2 range error condition.
	Type:	BOOL

Level: MED  
Label: CO2 Out Of Range  
VELA GUI/Membrane: CO2 Out Of Range

**10.7.2.10 ID: AlarmCO2SensorFault**  
Description: CO2 device reported a fault condition with CO2 sensor.  
Type: BOOL  
Level: MED  
Label: CO2 Sensor Fault  
VELA GUI/Membrane: CO2 Sensor Fault

**10.7.2.11 ID: AlarmCO2SensorOverTemp**  
Description: CO2 device reported a fault condition due to temperature.  
Type: BOOL  
Level: MED  
Label: CO2 Sensor Temp  
VELA GUI/Membrane: CO2 Sensor Temp

**10.7.2.12 ID: AlarmCO2ZeroRequired**  
Description: CO2 device reported sensor requires to be initialized to zero.  
Type: BOOL  
Level: MED  
Label: CO2 Zero Req'd  
VELA GUI/Membrane: CO2 Zero Req'd

**10.7.2.13 ID: AlarmDefaults**  
Description: A condition has occurred which has caused the ventilator to use default patient settings.  
Type: BOOL  
Level: MED  
Label: DEFAULTS  
VELA GUI/Membrane: DEFAULTS

**10.7.2.14 ID: AlarmDirtyFilter**  
Description: The inlet air filter may need to be cleaned.  
Type: BOOL  
Level: ALERT  
Label: CHECK FILTER  
VELA GUI/Membrane: CHECK FILTER

**10.7.2.15 ID: AlarmEEPROMFault**  
Description: The EEPROM can no longer be written.  
Type: BOOL

	Level:	ALERT
	Label:	EEPROM FAULT
	VELA GUI/Membrane:	EEPROM FAULT
<b>10.7.2.16</b>	<b>ID:</b>	<b>AlarmEndTidalCO2High</b>
	Description:	Exceeded EtCO2 High Limit
	Type:	BOOL
	Level:	LOW
	Label:	High EtCO2
	VELA GUI/Membrane:	High EtCO2
<b>10.7.2.17</b>	<b>ID:</b>	<b>AlarmEndTidalCO2Invalid</b>
	Description:	Invalid EtCO2.
	Type:	BOOL
	Level:	MED
	Label:	Invalid EtCO2
	VELA GUI/Membrane:	Invalid EtCO2
<b>10.7.2.18</b>	<b>ID:</b>	<b>AlarmEndTidalCO2Low</b>
	Description:	Below EtCO2 Low Limit
	Type:	BOOL
	Level:	LOW
	Label:	Low EtCO2
	VELA GUI/Membrane:	Low EtCO2
<b>10.7.2.19</b>	<b>ID:</b>	<b>AlarmFanFail</b>
	Description:	Internal cooling/enclosure ventilation fan has failed.
	Type:	BOOL
	Level:	MED
	Label:	FAN FAILURE
	VELA GUI/Membrane:	FAN FAILURE
<b>10.7.2.20</b>	<b>ID:</b>	<b>AlarmFiO2Range</b>
	Description:	The O <sub>2</sub> reading is outside the range of the delivered setting on the front panel.
	Type:	BOOL
	Level:	HIGH
	Label:	%O2 RANGE ERROR
	VELA GUI/Membrane:	%O2 RANGE ERROR
<b>10.7.2.21</b>	<b>ID:</b>	<b>AlarmFlowSensorDisc</b>
	Description:	Tidal volume accuracy has been degraded to the lowest accuracy supported by the ventilator.
	Type:	BOOL
	Level:	ALERT
	Label:	FLOW SENSOR DISC

VELA GUI/Membrane: FLOW SENSOR DISC

- 10.7.2.22 ID: AlarmHwFault**  
Description: A general equipment failure condition has been detected  
Type: BOOL  
Level: HIGH  
Label: H/W FAULT  
VELA GUI/Membrane: H/W FAULT
- 10.7.2.23 ID: AlarmInop**  
Description: Indicates ventilator is unable to ventilate the patient.  
Most likely due to detection of an unrecoverable internal problem or failure.  
Type: BOOL  
Level: HIGH  
Label: VENT INOP  
VELA GUI/Membrane: VENT INOP
- 10.7.2.24 ID: AlarmInvalidSN**  
Description: The ventilator's internal serial number setting is not correct.  
Type: BOOL  
Level: LOW  
Label: INVALID SERIAL NUMBER  
VELA GUI/Membrane: INVALID SERIAL NUMBER
- 10.7.2.25 ID: AlarmLossO2**  
Description: Ventilator cannot detect a source of oxygen gas supply.  
Not connected or insufficient pressure.  
Type: BOOL  
Level: HIGH  
Label: O2 Inlet LOW  
VELA GUI/Membrane: O2 Inlet LOW
- 10.7.2.26 ID: AlarmMotorFault**  
Description: Turbine motor has indicated a failure condition.  
Type: BOOL  
Level: HIGH  
Label: Motor Fault  
VELA GUI/Membrane: MOTOR FAULT
- 10.7.2.27 ID: AlarmNoCalData**  
Description: The transducer calibration data contained in the EEPROM is invalid.  
Type: BOOL  
Level: LOW



	Label:	NO CAL DATA
	VELA GUI/Membrane:	NO CAL DATA
<b>10.7.2.28</b>	<b>ID:</b>	<b>AlarmFiO2Cal</b>
	Description:	The O <sub>2</sub> reading is outside the range of the delivered setting on the front panel.
	Type:	BOOL
	Level:	MED
	Label:	CHECK O2 CAL
	VELA GUI/Membrane:	CHECK O2 CAL
<b>10.7.2.29</b>	<b>ID:</b>	<b>AlarmO2PressHigh</b>
	Description:	Indicates that the pressure measured at the regulated oxygen outlet is too high.
	Type:	BOOL
	Level:	MED
	Label:	O2 INLET HIGH
	VELA GUI/Membrane:	O2 INLET HIGH
<b>10.7.2.30</b>	<b>ID:</b>	<b>AlarmFiO2SensorFail</b>
	Description:	A failure has been detected in the FiO <sub>2</sub> monitor sensor.
	Type:	BOOL
	Level:	MED
	Label:	O2 SENSOR FAILURE
	VELA GUI/Membrane:	O2 SENSOR FAILURE
<b>10.7.2.31</b>	<b>ID:</b>	<b>AlarmPatientDefaults</b>
	Description:	A condition has occurred which has caused the ventilator to use default patient settings.
	Type:	BOOL
	Level:	ALERT
	Label:	PATIENT DEFAULTS
	VELA GUI/Membrane:	PATIENT DEFAULTS
<b>10.7.2.32</b>	<b>ID:</b>	<b>AlarmPeepHigh</b>
	Description:	Indicates that the circuit pressure is not returning to the set PEEP level during exhalation.
	Type:	BOOL
	Level:	HIGH
	Label:	HIGH PEEP
	VELA GUI/Membrane:	HIGH PEEP
<b>10.7.2.33</b>	<b>ID:</b>	<b>AlarmPpeakHigh</b>
	Description:	Airway pressure has exceeded the high peak pressure limit.
	Type:	BOOL

Level: HIGH  
Label: HIGH PIP  
VELA GUI/Membrane: HIGH PIP

- 10.7.2.34 ID: AlarmPpeakHighExt**  
Description: High P<sub>PEAK</sub> or occlusion persistent for greater than 5 seconds.  
Type: BOOL  
Level: HIGH  
Label: HIGH PIP, SUST.  
VELA GUI/Membrane: HIGH PIP, SUST.
- 10.7.2.35 ID: AlarmPpeakLow**  
Description: Airway pressure did not exceed the low peak pressure limit for the previous breath cycle.  
Type: BOOL  
Level: HIGH  
Label: LOW PIP  
VELA GUI/Membrane: LOW PIP
- 10.7.2.36 ID: AlarmPwrAcLoss**  
Description: AC power has been removed from the ventilator (i.e. power cord disconnect or loss of supply power).  
Type: BOOL  
Level: MED  
Label: ON BATTERY POWER  
VELA GUI/Membrane: ON BATTERY POWER
- 10.7.2.37 ID: AlarmPwrBattLow**  
Description: Batteries have been depleted to a level that indicates two minutes or less of safe operation.  
Type: BOOL  
Level: HIGH  
Label: LOW BATTERY  
VELA GUI/Membrane: LOW BATTERY
- 10.7.2.38 ID: AlarmPwrBattMed**  
Description: The power available in the main battery has dropped below the medium power threshold.  
Type: BOOL  
Level: MED  
Label: MED BATTERY  
VELA GUI/Membrane: MED BATTERY

- 10.7.2.39 ID: AlarmRateHigh**  
Description: Respiratory Rate exceeded the high inspirations per minute limit.  
Type: BOOL  
Level: MED  
Label: HIGH RATE  
VELA GUI/Membrane: HIGH RATE
- 10.7.2.40 ID: AlarmTransducerFault**  
Description: Invalid range detected during the zero-pressure test readings.  
Type: BOOL  
Level: MED  
Label: XDCR FAULT  
VELA GUI/Membrane: XDCR FAULT
- 10.7.2.41 ID: AlarmVeLow**  
Description: Minute Volume dropped below the low minute volume limit.  
Type: BOOL  
Level: HIGH  
Label: LOW Ve  
VELA GUI/Membrane: LOW Ve
- 10.7.2.42 ID: AlarmHistApnea**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmApnea description above.  
Type: BOOL  
Level: HIGH  
Label: APNEA INTERVAL  
VELA GUI/Membrane: APNEA INTERVAL
- 10.7.2.43 ID: AlarmHistCheckEvents**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmCheckEvents description above.  
Type: BOOL  
Level: MED  
Label: CHECK EVENTS  
VELA GUI/Membrane: CHECK EVENTS
- 10.7.2.44 ID: AlarmHistCircDisc**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmCircDisc description above.

Type: BOOL  
Level: HIGH  
Label: CIRCUIT FAULT  
VELA GUI/Membrane: CIRCUIT FAULT

**10.7.2.45 ID: AlarmHistClockBattLow**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmClockBattLow description above.  
Type: BOOL  
Level: MED  
Label: LOW CLOCK BATTERY  
VELA GUI/Membrane: LOW CLOCK BATTERY

**10.7.2.46 ID: AlarmHistCO2CheckAirwayAdapter**  
Description: CO2 device reported a CO2 Airway Adapter needs to be checked (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Check Adapter  
VELA GUI/Membrane: CO2 Check Adapter

**10.7.2.47 ID: AlarmHistCO2CommunicationError**  
Description: CO2 device reported a communication error (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Comms Error  
VELA GUI/Membrane: CO2 Comms Error

**10.7.2.48 ID: AlarmHistCO2OutOfRange**  
Description: CO2 device reported a CO2 range error condition (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Out Of Range  
VELA GUI/Membrane: CO2 Out Of Range

**10.7.2.49 ID: AlarmHistCO2SensorFault**  
Description: CO2 device reported a fault condition with CO2 sensor (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Sensor Fault  
VELA GUI/Membrane: CO2 Sensor Fault

- 10.7.2.50 ID: AlarmHistCO2SensorOverTemp**  
Description: CO2 device reported a fault condition due to temperature (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Sensor Temp  
VELA GUI/Membrane: CO2 Sensor Temp
- 10.7.2.51 ID: AlarmHistCO2ZeroRequired**  
Description: CO2 device reported sensor requires to be initialized to zero (historical / not cleared).  
Type: BOOL  
Level: MED  
Label: CO2 Zero Req'd  
VELA GUI/Membrane: CO2 Zero Req'd
- 10.7.2.52 ID: AlarmHistDefaults**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmDefaults description above.  
Type: BOOL  
Level: MED  
Label: DEFAULTS  
VELA GUI/Membrane: DEFAULTS
- 10.7.2.53 ID: AlarmHistDirtyFilter**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmDirtyFilter description above.  
Type: BOOL  
Level: ALERT  
Label: CHECK FILTER  
VELA GUI/Membrane: CHECK FILTER
- 10.7.2.54 ID: AlarmHistEEPROMFault**  
Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmEEPROMFault description above.  
Type: BOOL  
Level: ALERT  
Label: EEPROM FAULT  
VELA GUI/Membrane: EEPROM FAULT
- 10.7.2.55 ID: AlarmHistEndTidalCO2High**  
Description: Exceeded EtCO2 High Limit (historical / not cleared).  
Type: BOOL

Level: LOW  
 Label: High EtCO2  
 VELA GUI/Membrane: High EtCO2

- 10.7.2.56 ID: AlarmHistEndTidalCO2Invalid**  
 Description: Invalid EtCO2 (historical / not cleared).  
 Type: BOOL  
 Level: MED  
 Label: Invalid EtCO2  
 VELA GUI/Membrane: Invalid EtCO2
- 10.7.2.57 ID: AlarmHistEndTidalCO2Low**  
 Description: Below EtCO2 Low Limit (historical / not cleared).  
 Type: BOOL  
 Level: LOW  
 Label: Low EtCO2  
 VELA GUI/Membrane: Low EtCO2
- 10.7.2.58 ID: AlarmHistFanFail**  
 Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmFanFail description above.  
 Type: BOOL  
 Level: MED  
 Label: FAN FAILURE  
 VELA GUI/Membrane: FAN FAILURE
- 10.7.2.59 ID: AlarmHistFiO2Range**  
 Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmFiO2Range description above.  
 Type: BOOL  
 Level: HIGH  
 Label: %O2 RANGE ERROR  
 VELA GUI/Membrane: %O2 RANGE ERROR
- 10.7.2.60 ID: AlarmHistFlowSensorDisc**  
 Description: Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmFlowSensorDisc description above.  
 Type: BOOL  
 Level: ALERT  
 Label: FLOW SENSOR DISC  
 VELA GUI/Membrane: FLOW SENSOR DISC

<b>10.7.2.61</b>	<b>ID:</b>	<b>AlarmHistHwFault</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmHwFault description above.
	Type:	BOOL
	Level:	HIGH
	Label:	H/W FAULT
	VELA GUI/Membrane:	H/W FAULT
<b>10.7.2.62</b>	<b>ID:</b>	<b>AlarmHistInop</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmInop description above.
	Type:	BOOL
	Level:	HIGH
	Label:	VENT INOP
	VELA GUI/Membrane:	VENT INOP
<b>10.7.2.63</b>	<b>ID:</b>	<b>AlarmHistInvalidSN</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmInvalidSN description above.
	Type:	BOOL
	Level:	LOW
	Label:	INVALID SERIAL NUMBER
	VELA GUI/Membrane:	INVALID SERIAL NUMBER
<b>10.7.2.64</b>	<b>ID:</b>	<b>AlarmHistLossO2</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmLossO2 description above.
	Type:	BOOL
	Level:	HIGH
	Label:	O2 Inlet LOW
	VELA GUI/Membrane:	O2 Inlet LOW
<b>10.7.2.65</b>	<b>ID:</b>	<b>AlarmHistMotorFault</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmMotorFault description above.
	Type:	BOOL
	Level:	HIGH
	Label:	Motor Fault
	VELA GUI/Membrane:	MOTOR FAULT

<b>10.7.2.66</b>	<b>ID:</b>	<b>AlarmHistNoCalData</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmNoCalData description above.
	Type:	BOOL
	Level:	LOW
	Label:	NO CAL DATA
	VELA GUI/Membrane:	NO CAL DATA
<b>10.7.2.67</b>	<b>ID:</b>	<b>AlarmHistFiO2Cal</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmFiO2Cal description above.
	Type:	BOOL
	Level:	MED
	Label:	CHECK O2 CAL
	VELA GUI/Membrane:	CHECK O2 CAL
<b>10.7.2.68</b>	<b>ID:</b>	<b>AlarmHistO2PressHigh</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmO2PressHigh description above.
	Type:	BOOL
	Level:	MED
	Label:	O2 INLET HIGH
	VELA GUI/Membrane:	O2 INLET HIGH
<b>10.7.2.69</b>	<b>ID:</b>	<b>AlarmHistFiO2SensorFail</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmFiO2SensorFail description above.
	Type:	BOOL
	Level:	MED
	Label:	O2 SENSOR FAILURE
	VELA GUI/Membrane:	O2 SENSOR FAILURE
<b>10.7.2.70</b>	<b>ID:</b>	<b>AlarmHistPatientDefaults</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPatientDefaults description above.
	Type:	BOOL
	Level:	ALERT
	Label:	PATIENT DEFAULTS
	VELA GUI/Membrane:	PATIENT DEFAULTS



<b>10.7.2.71</b>	<b>ID:</b>	<b>AlarmHistPeepHigh</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPeepHigh description above.
	Type:	BOOL
	Level:	HIGH
	Label:	HIGH PEEP
	VELA GUI/Membrane:	HIGH PEEP
<b>10.7.2.72</b>	<b>ID:</b>	<b>AlarmHistPpeakHigh</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPpeakHigh description above.
	Type:	BOOL
	Level:	HIGH
	Label:	HIGH PIP
	VELA GUI/Membrane:	HIGH PIP
<b>10.7.2.73</b>	<b>ID:</b>	<b>AlarmHistPpeakHighExt</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPpeakHighExt description above.
	Type:	BOOL
	Level:	HIGH
	Label:	HIGH PIP, SUST.
	VELA GUI/Membrane:	HIGH PIP, SUST.
<b>10.7.2.74</b>	<b>ID:</b>	<b>AlarmHistPpeakLow</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPpeakLow description above.
	Type:	BOOL
	Level:	HIGH
	Label:	LOW PIP
	VELA GUI/Membrane:	LOW PIP
<b>10.7.2.75</b>	<b>ID:</b>	<b>AlarmHistPwrAcLoss</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPwrAcLoss description above.
	Type:	BOOL
	Level:	MED
	Label:	ON BATTERY POWER
	VELA GUI/Membrane:	ON BATTERY POWER

<b>10.7.2.76</b>	<b>ID:</b>	<b>AlarmHistPwrBattLow</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPwrBattLow description above.
	Type:	BOOL
	Level:	HIGH
	Label:	LOW BATTERY
	VELA GUI/Membrane:	LOW BATTERY
<b>10.7.2.77</b>	<b>ID:</b>	<b>AlarmHistPwrBattMed</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmPwrBattMed description above.
	Type:	BOOL
	Level:	MED
	Label:	MED BATTERY
	VELA GUI/Membrane:	MED BATTERY
<b>10.7.2.78</b>	<b>ID:</b>	<b>AlarmHistRateHigh</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmRateHigh description above.
	Type:	BOOL
	Level:	MED
	Label:	HIGH RATE
	VELA GUI/Membrane:	HIGH RATE
<b>10.7.2.79</b>	<b>ID:</b>	<b>AlarmHistTransducerFault</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmTransducerFault description above.
	Type:	BOOL
	Level:	MED
	Label:	XDCR FAULT
	VELA GUI/Membrane:	XDCR FAULT
<b>10.7.2.80</b>	<b>ID:</b>	<b>AlarmHistVeLow</b>
	Description:	Indication that alarm asserted in the past, is no longer active and has not been reset. See AlarmVeLow description above.
	Type:	BOOL
	Level:	HIGH
	Label:	LOW Ve
	VELA GUI/Membrane:	LOW Ve

<b>10.7.2.81</b>	<b>ID:</b>	<b>AlarmSilence</b>
	Description:	Active/Inactive state of the capability to locally silence the audible ventilator alarms.
	Type:	BOOL
	Label:	ALARM SILENCE
	VELA GUI/Membrane:	ALARM SILENCE

## 10.8 Scalar Class

### 10.8.1 AVEA Ventilator – Scalar Class

<b>10.8.1.1</b>	<b>ID:</b>	<b>WaveAnlg0</b>
	Description:	Echo of Analog Input Channel 0
	Type:	WORD
	Scale:	3
	Range:	0 – 12000
	Epoch:	500
	Size:	50
	Label:	Analog 0
	AVEA GUI/Membrane:	Analog 0

<b>10.8.1.2</b>	<b>ID:</b>	<b>WaveAnlg1</b>
	Description:	Echo of Analog Input Channel 1
	Type:	WORD
	Scale:	3
	Range:	0 – 12000
	Epoch:	500
	Size:	50
	Label:	Analog 1
	AVEA GUI/Membrane:	Analog 1

<b>10.8.1.3</b>	<b>ID:</b>	<b>WaveFexp</b>
	Description:	Expiratory Flow
	Type:	WORD
	Scale:	2
	Range:	-30000 - 30000
	Epoch:	500
	Size:	50
	Units:	L/min
	Label:	Fexp
	AVEA GUI/Membrane:	Fexp

<b>10.8.1.4</b>	<b>ID:</b>	<b>WaveFinsp</b>
	Description:	Inspiratory Flow
	Type:	WORD

Scale: 2  
Range: -30000 - 30000  
Epoch: 500  
Size: 50  
Units: L/min  
Label: Finsp  
AVEA GUI/Membrane: Finsp

**10.8.1.5 ID: WaveFlow**  
Description: Airway Net Flow  
Type: WORD  
Scale: 2  
Range: -30000 - 30000  
Epoch: 500  
Size: 50  
Units: L/min  
Label: Flow  
AVEA GUI/Membrane: Flow

**10.8.1.6 ID: WaveMetric**  
Description: Metrics with respect to waveforms - refer to section 10.2.5 for bit field definitions  
Type: UWORD  
Epoch: 500  
Size: 50  
Units:  
Label:  
AVEA GUI/Membrane:

**10.8.1.7 ID: WavePaw**  
Description: Airway Pressure  
Type: WORD  
Scale: 2  
Range: -6000 - 12000  
Epoch: 500  
Size: 50  
Units: cmH2O  
Label: Paw  
AVEA GUI/Membrane: Paw

**10.8.1.8 ID: WavePes**  
Description: Esophageal Pressure  
Type: WORD  
Scale: 2  
Range: -6000 - 12000  
Epoch: 500

Size: 50  
 Units: cmH2O  
 Label: Pes  
 AVEA GUI/Membrane: Pes

**10.8.1.9 ID: WavePinsp**  
 Description: Airway Pressure at Machine Outlet  
 Type: WORD  
 Scale: 2  
 Range: -6000 - 12000  
 Epoch: 500  
 Size: 50  
 Units: cmH2O  
 Label: Pinsp  
 AVEA GUI/Membrane: Pinsp

**10.8.1.10 ID: WavePtp**  
 Description: Transpulmonary Pressure  
 Type: WORD  
 Scale: 2  
 Range: -6000 - 12000  
 Epoch: 500  
 Size: 50  
 Units: cmH2O  
 Label: Ptp  
 AVEA GUI/Membrane: Ptp

**10.8.1.11 ID: WavePtr**  
 Description: Tracheal Pressure  
 Type: WORD  
 Scale: 2  
 Range: -6000 - 12000  
 Epoch: 500  
 Size: 50  
 Units: cmH2O  
 Label: Ptr  
 AVEA GUI/Membrane: Ptr

**10.8.1.12 ID: WaveSpO2Pleth**  
 Description: Raw IR signal from the MS-11 PCB. The range value is inverted and scaled from the raw data sent by the oximeter (oximeter +127 = 0%, oximeter – 128 = 100%).  
 Type: WORD  
 Range: 0 – 100  
 Resolution: 1  
 Epoch: 100

Size: 10  
 Units: %  
 Label: Pleth  
 AVEA GUI/Membrane: Pleth

**10.8.1.13 ID: WaveVt**  
 Description: Airway Tidal Volume  
 Type: WORD  
 Scale: 1  
 Range: -10000 - 30000  
 Epoch: 500  
 Size: 50  
 Units: mL  
 Label: Vt  
 AVEA GUI/Membrane: Vt

## 10.8.2 VELA Ventilator – Scalar Class

**10.8.2.1 ID: WaveFlow**  
 Description: Airway Net Flow  
 Type: WORD  
 Scale: 2  
 Range: -30000 - 30000  
 Epoch: 500  
 Size: 50  
 Units: L/min  
 Label: V (lpm)  
 VELA GUI/Membrane: Flow

**10.8.2.2 ID: WaveMetric**  
 Description: Metrics with respect to waveforms - refer to section 10.2. for bit field definitions.  
 Type: UWORD  
 Range: 0 - 127  
 Epoch: 500  
 Size: 50  
 Units:  
 Label:  
 VELA GUI/Membrane:

**10.8.2.3 ID: WavePaw**  
 Description: Airway Pressure  
 Type: WORD  
 Scale: 2  
 Range: -6000 - 12000

Epoch: 500  
 Size: 50  
 Units: cmH2O  
 Label: Paw  
 VELA GUI/Membrane: Paw

**10.8.2.4 ID: WaveVt**  
 Description: Airway Tidal Volume  
 Type: WORD  
 Scale: 1  
 Range: -7000 - 21000  
 Epoch: 500  
 Size: 50  
 Units: ml  
 Label: Vt(ml)  
 VELA GUI/Membrane: Vt(ml)

## 11 Foreign Language Support

The VOX Protocol shall support all languages selectable from the graphical user interface of both the AVEA and VELA ventilators.

## 12 VOXP Profile, Config, Data, and Link Messages

Below are AVEA VOXP messages examples, which have been provided for reference.

### 12.1 Profile Message – AVEA Adult

```
<profile model="Avea Comp" profileVersion="2.0" voxpVersion="3.0" textEncoding="UTF-16"
msgID="0003">
 <unit class="setting" ID="SetFiO2" type="WORD" resolution="0001" range="0015:0064"
 units="0025" label="00460069004F0032"/>
 <unit class="setting" ID="SetFlowBias" type="WORD" scale="E+1" resolution="0001"
 range="0004:0032" units="004C002F006D0069006E"
 label="004200690061007300200046006C006F0077"/>
 <unit class="setting" ID="SetFlowCycle" type="WORD" resolution="0005"
 range="0000:002D" units="0025"
 label="0046006C006F00770020004300790063006C0065"/>
 <unit class="setting" ID="SetFlowCyclePsv" type="WORD" resolution="0005"
 range="0005:002D" units="0025"
 label="0050005300560020004300790063006C0065"/>
 <unit class="setting" ID="SetFlowDemand" type="WORD" resolution="0001"
 range="0000:0001" units="0030003D004F0066006600200031003D004F006E"
 label="00440065006D0061006E006400200046006C006F0077"/>
 <unit class="setting" ID="SetFlowInsp" type="WORD" scale="E+1" resolution="0001"
 range="0003:0096" units="004C002F006D0069006E"
 label="005000650061006B00200046006C006F0077"/>
```

```

<unit class="setting" ID="SetPauseInsp" type="WORD" scale="E+2" resolution="0001"
 range="0000:012C" units="007300650063"
 label="0049006E00730070002000500061007500730065"/>
<unit class="setting" ID="SetPresHigh" type="WORD" resolution="0001"
 range="0000:005A" units="0063006D00480032004F"
 label="005000720065007300200048006900670068"/>
<unit class="setting" ID="SetPresInsp" type="WORD" resolution="0001"
 range="0000:005A" units="0063006D00480032004F"
 label="0049006E0073007000200050007200650073"/>
<unit class="setting" ID="SetPresLow" type="WORD" resolution="0001"
 range="0000:002D" units="0063006D00480032004F"
 label="00500072006500730020004C006F0077"/>
<unit class="setting" ID="SetPresNasalCPAP" type="WORD" resolution="0001"
 range="0002:000A" units="0063006D00480032004F"
 label="006E0043005000410050"/>
<unit class="setting" ID="SetPresPeep" type="WORD" resolution="0001"
 range="0000:0032" units="0063006D00480032004F" label="0050004500450050"/>
<unit class="setting" ID="SetPresPsv" type="WORD" resolution="0001"
 range="0000:005A" units="0063006D00480032004F" label="005000530056"/>
<unit class="setting" ID="SetRate" type="WORD" resolution="0001" range="0001:0078"
 units="00620070006D" label="0052006100740065"/>
<unit class="setting" ID="SetRiseInsp" type="WORD" resolution="0001"
 range="0001:0009" units="" label="0049006E0073007000200052006900730065"/>
<unit class="setting" ID="SetRisePsv" type="WORD" resolution="0001" range="0001:0009"
 units="" label="00500053005600200052006900730065"/>
<unit class="setting" ID="SetRiseVsync" type="WORD" resolution="0001"
 range="0001:0009" units=""
 label="005600730079006E006300200052006900730065"/>
<unit class="setting" ID="SetTimeHigh" type="WORD" scale="E+1" resolution="0001"
 range="0002:012C" units="007300650063"
 label="00540069006D006500200048006900670068"/>
<unit class="setting" ID="SetTimeHighPsv" type="WORD" resolution="0001"
 range="0000:0001" units="0030003D004F0066006600200031003D004F006E"
 label="0054002000480069006700680020005000530056"/>
<unit class="setting" ID="SetTimeHighSync" type="WORD" resolution="0001"
 range="0000:0032" units="0025"
 label="005400200048006900670068002000530079006E0063"/>
<unit class="setting" ID="SetTimeInsp" type="WORD" scale="E+2" resolution="0001"
 range="0014:01F4" units="007300650063"
 label="0049006E00730070002000540069006D0065"/>
<unit class="setting" ID="SetTimeLow" type="WORD" scale="E+1" resolution="0001"
 range="0002:012C" units="007300650063"
 label="00540069006D00650020004C006F0077"/>
<unit class="setting" ID="SetTimeLowSync" type="WORD" resolution="0001"
 range="0000:0032" units="0025"
 label="00540020004C006F0077002000530079006E0063"/>

```



```

<unit class="setting" ID="SetTmaxPsv" type="WORD" scale="E+2" resolution="0001"
 range="0014:01F4" units="007300650063"
 label="00500053005600200054006D00610078"/>
<unit class="setting" ID="SetTrigFlow" type="WORD" scale="E+1" resolution="0001"
 range="0001:00C8" units="004C002F006D0069006E"
 label="0046006C006F007700200054007200690067"/>
<unit class="setting" ID="SetTrigPres" type="WORD" scale="E+1" resolution="0001"
 range="0001:00C8" units="0063006D00480032004F"
 label="005000720065007300200054007200690067"/>
<unit class="setting" ID="SetVol" type="WORD" scale="E+4" resolution="0001"
 range="000A:00FA" units="004C" label="0056006F006C0075006D0065"/>
<unit class="setting" ID="SetVolAssured" type="WORD" scale="E+4" resolution="0001"
 range="0000:00FA" units="004C" label="004D00610063006800200056006F006C"/>
<unit class="setting" ID="SetVolLimit" type="WORD" scale="E+4" resolution="0001"
 range="000A:00FA" units="004C"
 label="0056006F006C0020004C0069006D00690074"/>
<unit class="setting" ID="SetVolSigh" type="WORD" resolution="0001" range="0000:0001"
 units="0030003D004F0066006600200031003D004F006E"
 label="0053006900670068"/>
<unit class="setting" ID="SetVolWave" type="ENUM"
 label="00570061007600650066006F0072006D">
 <enum value="0000" label="005300510055004100520045"/>
 <enum value="0001"
 label="0044004500430045004C00450052004100540049004E0047"/>
</unit>
<unit class="setting" ID="SetVsync" type="WORD" resolution="0001" range="0000:0001"
 units="0030003D004F0066006600200031003D004F006E"
 label="005600730079006E0063"/>
<unit class="setting" ID="LimitApnea" type="WORD" resolution="0001" range="0006:003C"
 units="007300650063"
 label="00410070006E0065006100200049006E00740065007200760061006C"/>
<unit class="setting" ID="LimitPeepLow" type="WORD" resolution="0001"
 range="0000:003C" units="0063006D00480032004F"
 label="004C006F007700200050004500450050"/>
<unit class="setting" ID="LimitPpeakHigh" type="WORD" resolution="0001"
 range="000A:0069" units="0063006D00480032004F"
 label="004800690067006800200050007000650061006B"/>
<unit class="setting" ID="LimitPpeakLow" type="WORD" resolution="0001"
 range="0001:0063" units="0063006D00480032004F"
 label="004C006F007700200050007000650061006B"/>
<unit class="setting" ID="LimitRateHigh" type="WORD" resolution="0001"
 range="0001:00C8" units="00620070006D"
 label="004800690067006800200052006100740065"/>
<unit class="setting" ID="LimitVeHigh" type="WORD" scale="E+2" resolution="0001"
 range="0000:004B" units="004C" label="0048006900670068002000560065"/>

```

```

<unit class="setting" ID="LimitVeLow" type="WORD" scale="E+2" resolution="0001"
 range="0000:0032" units="004C" label="004C006F0077002000560065"/>
<unit class="setting" ID="LimitVteHigh" type="WORD" scale="E+4" resolution="0001"
 range="000A:012C" units="004C" label="00480069006700680020005600740065"/>
<unit class="setting" ID="LimitVteLow" type="WORD" scale="E+4" resolution="0001"
 range="0000:012C" units="004C" label="004C006F00770020005600740065"/>
<unit class="setting" ID="SetAAC" type="BOOL" resolution="1" range="0:1" units=""
 label="004100410043000A004F006E"/>
<unit class="setting" ID="SetCircComp" type="WORD" scale="E+1" resolution="0001"
 range="0000:004B" units="006D004C002F0063006D00480032004F"
 label="004300690072006300200043006F006D0070"/>
<unit class="setting" ID="SetEttDia" type="WORD" scale="E+1" resolution="0001"
 range="0014:0064" units="006D006D"
 label="004400690061006D0065007400650072"/>
<unit class="setting" ID="SetEttLen" type="WORD" scale="E+1" resolution="0001"
 range="0014:012C" units="0063006D" label="004C0065006E006700740068"/>
<unit class="setting" ID="SetHumidifier" type="BOOL" resolution="1" range="0:1" units=""
 label="00480055004D0049004400490046004900450052"/>
<unit class="setting" ID="SetLanguage" type="ENUM"
 label="004C0061006E00670075006100670065003A">
 <enum value="0000" label="0045006E0067006C006900730068"/>
 <enum value="0001" label="004600720061006E00E7006100690073"/>
 <enum value="0002" label="0044006500750074007300630068"/>
 <enum value="0003" label="004900740061006C00690061006E006F"/>
 <enum value="0004" label="0050006F0072007400750067007500EA0073"/>
 <enum value="0005" label="004500730070006100F1006F006C"/>
 <enum value="0006" label="6C498BED"/>
 <enum value="0007" label="004E0065006400650072006C0061006E00640073"/>
</unit>
<unit class="setting" ID="SetLeakComp" type="BOOL" resolution="1" range="0:1" units=""
 label="004C00650061006B00200043006F006D0070000A004F006E"/>
<unit class="setting" ID="SetMode" type="ENUM"
 label="004D004F00440045002000530045004C004500430054">
 <enum value="0001"
 label="00410050005200560020002F002000420049005000480041005300490043"/>
 <enum value="0002"
 label="00410050005200560020002F002000420049005000480041005300490043"/>
 <enum value="0003" label="0050005200560043002000530049004D0056"/>
 <enum value="0004" label="005000520056004300200041002F0043"/>
 <enum value="0005" label="00430050004100500020002F0020005000530056"/>
 <enum value="0006" label="005400430050004C002000530049004D0056"/>
 <enum value="0007" label="005400430050004C00200041002F0043"/>
 <enum value="0008" label="00430050004100500020002F0020005000530056"/>

```

```

 <enum value="0009"
 label="00500052004500530053005500520045002000530049004D0056"/>
 <enum value="000A"
 label="0050005200450053005300550052004500200041002F0043"/>
 <enum value="000B" label="00430050004100500020002F0020005000530056"/>
 <enum value="000C"
 label="0056004F004C0055004D0045002000530049004D0056"/>
 <enum value="000D" label="0056004F004C0055004D004500200041002F0043"/>
 <enum value="000E" label="004E006100730061006C00200043005000410050"/>
 </unit>
 <unit class="setting" ID="SetModellv" type="ENUM"
 label="0049004C00560020004D006F00640065003A">
 <enum value="0000" label="004F00660066"/>
 <enum value="0001" label="004D00610073007400650072"/>
 <enum value="0002" label="0053006C006100760065"/>
 </unit>
 <unit class="setting" ID="SetPatSize" type="ENUM"
 label="00500041005400490045004E0054002000530049005A0045002000530045004C004500430054">
 <enum value="0000" label="004E0065006F"/>
 <enum value="0001" label="005000650064"/>
 <enum value="0002" label="004100640075006C0074"/>
 </unit>
 <unit class="setting" ID="SetPatWt" type="WORD" scale="E+2" resolution="0001"
 range="0001:012C" units="006B0067"
 label="005000740020005700650069006700680074"/>
 <unit class="setting" ID="SetPresBaro" type="WORD" resolution="0001"
 range="0221:02F8" units="006D006D00480067"
 label="004200610072006F00200050007200650073"/>
 <unit class="setting" ID="SetIncrFiO2" type="WORD" resolution="0001"
 range="0000:004F" units="0025"
 label="0049006E006300720065006100730065002000460069004F0032003A"/>
 <unit class="setting" ID="SetSensitivityLowVte" type="WORD" resolution="0001"
 range="0001:0005" units=""
 label="004C006F0077002000560074006500200041006C00610072006D003A"/>
 <unit class="monitor" ID="MntrAutoPEEP" type="WORD" range="0000:0032"
 units="0063006D00480032004F" label="004100750074006F0050004500450050"/>
 <unit class="monitor" ID="MntrAutoPEEPdelta" type="WORD" range="0000:0032"
 units="0063006D00480032004F"
 label="0064004100750074006F0050004500450050"/>
 <unit class="monitor" ID="MntrAutoPEEPesoph" type="WORD" range="0000:0032"
 units="0063006D00480032004F"
 label="004100750074006F005000450045005000650073"/>
 <unit class="monitor" ID="MntrC20" type="WORD" scale="E+2" range="0000:01F4"
 units="" label="004300320030002F0043"/>

```

```

<unit class="monitor" ID="MntrCcw" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="004300630077"/>
<unit class="monitor" ID="MntrCdyn" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="004300640079006E"/>
<unit class="monitor" ID="MntrCdynNorm" type="WORD" scale="E+2" range="0000:01F4"
 units="006D004C002F0063006D00480032004F002F006B0067"
 label="004300640079006E002F006B0067"/>
<unit class="monitor" ID="MntrClung" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="0043006C0075006E0067"/>
<unit class="monitor" ID="MntrCstat" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="00430073007400610074"/>
<unit class="monitor" ID="MntrCstatNorm" type="WORD" scale="E+2" range="0000:01F4"
 units="006D004C002F0063006D00480032004F002F006B0067"
 label="00430073007400610074002F006B0067"/>
<unit class="monitor" ID="MntrFiO2" type="WORD" range="0000:0064" units="0025"
 label="00460069004F0032"/>
<unit class="monitor" ID="MntrIE" type="WORD" scale="E+1" range="FC19:03E7" units=""
 label="0049003A0045"/>
<unit class="monitor" ID="MntrLeak" type="WORD" range="0000:0064" units="0025"
 label="004C00650061006B"/>
<unit class="monitor" ID="MntrMIP" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="004D00490050"/>
<unit class="monitor" ID="MntrNcpapMeanFlow" type="WORD" scale="E+1"
 range="0000:0BB8" units="004C002F006D0069006E"
 label="004300500041005000200046006C006F0077"/>
<unit class="monitor" ID="MntrNcpapPres" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="006E0043005000410050"/>
<unit class="monitor" ID="MntrP100" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="0050003100300030"/>
<unit class="monitor" ID="MntrPair" type="WORD" range="0000:0050"
 units="0070007300690067" label="00410069007200200049006E006C00650074"/>
<unit class="monitor" ID="MntrPawDelta" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0064005000610077"/>
<unit class="monitor" ID="MntrPeep" type="WORD" range="0000:0032"
 units="0063006D00480032004F" label="0050004500450050"/>
<unit class="monitor" ID="MntrPefr" type="WORD" scale="E+1" range="0000:0BB8"
 units="004C002F006D0069006E" label="0050004500460052"/>
<unit class="monitor" ID="MntrPesDelta" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0064005000650073"/>
<unit class="monitor" ID="MntrPifr" type="WORD" scale="E+1" range="0000:0BB8"
 units="004C002F006D0069006E" label="0050004900460052"/>
<unit class="monitor" ID="MntrPmean" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0050006D00650061006E"/>
<unit class="monitor" ID="MntrPO2" type="WORD" range="0000:0050"
 units="0070007300690067" label="004F003200200049006E006C00650074"/>

```

```
<unit class="monitor" ID="MntrPpeak" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0050007000650061006B"/>
<unit class="monitor" ID="MntrPplat" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="00500070006C00610074"/>
<unit class="monitor" ID="MntrPplatPtp" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="00500074007000200050006C00610074"/>
<unit class="monitor" ID="MntrPtpPEEP" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="00500074007000200050004500450050"/>
<unit class="monitor" ID="MntrRate" type="WORD" range="0000:00C8"
 units="00620070006D" label="0052006100740065"/>
<unit class="monitor" ID="MntrRateMand" type="WORD" range="0000:00C8"
 units="00620070006D" label="004D0061006E006400200052006100740065"/>
<unit class="monitor" ID="MntrRateSpon" type="WORD" range="0000:00C8"
 units="00620070006D" label="00530070006F006E00200052006100740065"/>
<unit class="monitor" ID="MntrRimp" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="00520069006D0070"/>
<unit class="monitor" ID="MntrRlung" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="0052006C0075006E0067"/>
<unit class="monitor" ID="MntrRpeak" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="0052007000650061006B"/>
<unit class="monitor" ID="MntrRrs" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="005200720073"/>
<unit class="monitor" ID="MntrRSBIndex" type="WORD" range="0000:01F4"
 units="00620032002F006D0069006E002F004C" label="0066002F00560074"/>
<unit class="monitor" ID="MntrTe" type="WORD" scale="E+2" range="0000:270F"
 units="007300650063" label="00540065"/>
<unit class="monitor" ID="MntrTi" type="WORD" scale="E+2" range="0000:270F"
 units="007300650063" label="00540069"/>
<unit class="monitor" ID="MntrVdel" type="INT" scale="E+8" range="00000000:3B8B87C0"
 units="004C" label="005600640065006C"/>
<unit class="monitor" ID="MntrVeSpon" type="WORD" scale="E+2" range="0000:2706"
 units="004C" label="00530070006F006E002000560065"/>
<unit class="monitor" ID="MntrVeSponNorm" type="WORD" range="0000:03E7"
 units="006D004C002F006B0067"
 label="00530070006F006E002000560065002F006B0067"/>
<unit class="monitor" ID="MntrVeTotal" type="WORD" scale="E+2" range="0000:2706"
 units="004C" label="0054006F00740061006C002000560065"/>
<unit class="monitor" ID="MntrVeTotalNorm" type="WORD" range="0000:03E7"
 units="006D004C002F006B0067"
 label="0054006F00740061006C002000560065002F006B0067"/>
<unit class="monitor" ID="MntrVte" type="INT" scale="E+8" range="00000000:3B8B87C0"
 units="004C" label="005600740065"/>
```

```

<unit class="monitor" ID="MntrVteMand" type="INT" scale="E+8"
 range="00000000:3B8B87C0" units="004C"
 label="004D0061006E00640020005600740065"/>
<unit class="monitor" ID="MntrVteMandNorm" type="WORD" scale="E+2"
 range="0000:0BB8" units="006D004C002F006B0067"
 label="004D0061006E00640020005600740065002F006B0067"/>
<unit class="monitor" ID="MntrVteNorm" type="WORD" scale="E+2" range="0000:0BB8"
 units="006D004C002F006B0067" label="005600740065002F006B0067"/>
<unit class="monitor" ID="MntrVteSpon" type="INT" scale="E+8"
 range="00000000:3B8B87C0" units="004C"
 label="00530070006F006E0020005600740065"/>
<unit class="monitor" ID="MntrVteSponNorm" type="WORD" scale="E+2"
 range="0000:0BB8" units="006D004C002F006B0067"
 label="00530070006F006E0020005600740065002F006B0067"/>
<unit class="monitor" ID="MntrVti" type="INT" scale="E+8" range="00000000:3B8B87C0"
 units="004C" label="005600740069"/>
<unit class="monitor" ID="MntrVtiNorm" type="WORD" scale="E+2" range="0000:0BB8"
 units="006D004C002F006B0067" label="005600740069002F006B0067"/>
<unit class="monitor" ID="MntrWobImposed" type="WORD" scale="E+2"
 range="0000:07D0" units="006A006F0075006C00650073002F004C"
 label="0057004F00420069"/>
<unit class="monitor" ID="MntrWobPatient" type="WORD" scale="E+2" range="0000:07D0"
 units="006A006F0075006C00650073002F004C" label="0057004F00420070"/>
<unit class="monitor" ID="MntrWobVent" type="WORD" scale="E+2" range="0000:07D0"
 units="006A006F0075006C00650073002F004C" label="0057004F00420076"/>
<unit class="alarm" ID="AlarmActive" type="BOOL" label=""/>
<unit class="alarm" ID="AlarmActivePriority" type="ENUM" label="">
 <enum value="0001" label="0048004900470048"/>
 <enum value="0002" label="004D00450044"/>
 <enum value="0003" label="004C004F0057"/>
 <enum value="0004" label="0041004C004500520054"/>
</unit>
<unit class="alarm" ID="AlarmApnea" type="BOOL" level="HIGH"
 label="00410050004E0045004100200049004E00540045005200560041004C"/>
<unit class="alarm" ID="AlarmCircDisc" type="BOOL" level="HIGH"
 label="004300490052004300550049005400200044004900530043004F004E004E0045
 00430054"/>
<unit class="alarm" ID="AlarmFanFail" type="BOOL" level="LOW"
 label="00460041004E0020004600410049004C005500520045"/>
<unit class="alarm" ID="AlarmFiO2High" type="BOOL" level="HIGH"
 label="0048004900470048002000460069004F0032"/>
<unit class="alarm" ID="AlarmFiO2Low" type="BOOL" level="HIGH"
 label="004C004F0057002000460069004F0032"/>
<unit class="alarm" ID="AlarmIlvSlaveDisc" type="BOOL" level="HIGH"
 label="0049004C005600200044004900530043004F004E004E004500430054"/>

```

```

<unit class="alarm" ID="AlarmInop" type="BOOL" level="HIGH"
 label="00560045004E005400200049004E004F0050"/>
<unit class="alarm" ID="AlarmInvalidGasId" type="BOOL" level="MED"
 label="0049004E00560041004C004900440020004700410053002000490044"/>
<unit class="alarm" ID="AlarmLimitIE" type="BOOL" level="LOW"
 label="0049003A00450020004C0049004D00490054"/>
<unit class="alarm" ID="AlarmLimitTi" type="BOOL" level="LOW"
 label="004D0041005800200049004E00530050002000540049004D0045"/>
<unit class="alarm" ID="AlarmLimitVol" type="BOOL" level="LOW"
 label="0056004F004C0055004D00450020004C0049004D00490054"/>
<unit class="alarm" ID="AlarmLossAir" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004100490052"/>
<unit class="alarm" ID="AlarmLossGas" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004700410053"/>
<unit class="alarm" ID="AlarmLossHeliox" type="BOOL" level="HIGH"
 label="004C004F005300530020004F0046002000480045004C0049004F0058"/>
<unit class="alarm" ID="AlarmLossO2" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004F0032"/>
<unit class="alarm" ID="AlarmNcpapHigh" type="BOOL" level="HIGH"
 label="00480049004700480020006E004300500041005000200050005200450053"/>
<unit class="alarm" ID="AlarmNcpapHighPresLimit" type="BOOL" level="HIGH"
 label="006E0043005000410050002000500052004500530020004C0049004D0049005
 4"/>
<unit class="alarm" ID="AlarmNcpapLow" type="BOOL" level="HIGH"
 label="004C004F00570020006E004300500041005000200050005200450053"/>
<unit class="alarm" ID="AlarmOcclusion" type="BOOL" level="HIGH"
 label="00430049005200430055004900540020004F00430043004C005500530049004F
 004E"/>
<unit class="alarm" ID="AlarmOpenSV" type="BOOL" level="HIGH"
 label="005300410046004500540059002000560041004C00560045"/>
<unit class="alarm" ID="AlarmPeepLow" type="BOOL" level="HIGH"
 label="004C004F005700200050004500450050"/>
<unit class="alarm" ID="AlarmPpeakHigh" type="BOOL" level="HIGH"
 label="004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmPpeakHighExt" type="BOOL" level="HIGH"
 label="0045005800540020004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmPpeakLow" type="BOOL" level="HIGH"
 label="004C004F005700200050007000650061006B"/>
<unit class="alarm" ID="AlarmPwrAcLoss" type="BOOL" level="HIGH"
 label="004C004F005300530020004F004600200041002F0043"/>
<unit class="alarm" ID="AlarmPwrBattLow" type="BOOL" level="HIGH"
 label="004C004F005700200042004100540054004500520059"/>
<unit class="alarm" ID="AlarmRateHigh" type="BOOL" level="MED"
 label="004800490047004800200052004100540045"/>
<unit class="alarm" ID="AlarmSilence" type="BOOL" label=""/>

```

```

<unit class="alarm" ID="AlarmTest" type="BOOL" level="HIGH"
 label="0041004C00410052004D00200054004500530054"/>
<unit class="alarm" ID="AlarmVeHigh" type="BOOL" level="MED"
 label="0048004900470048002000560065"/>
<unit class="alarm" ID="AlarmVeLow" type="BOOL" level="HIGH"
 label="004C004F0057002000560065"/>
<unit class="alarm" ID="AlarmVteLow" type="BOOL" level="HIGH"
 label="004C004F00570020005600740065"/>
<unit class="alarm" ID="AlarmVtHigh" type="BOOL" level="LOW"
 label="00480049004700480020005600740065"/>
<unit class="alarm" ID="AlarmHistApnea" type="BOOL" level="HIGH"
 label="00410050004E0045004100200049004E00540045005200560041004C"/>
<unit class="alarm" ID="AlarmHistCircDisc" type="BOOL" level="HIGH"
 label="004300490052004300550049005400200044004900530043004F004E004E0045
00430054"/>
<unit class="alarm" ID="AlarmHistFanFail" type="BOOL" level="LOW"
 label="00460041004E0020004600410049004C005500520045"/>
<unit class="alarm" ID="AlarmHistFiO2High" type="BOOL" level="HIGH"
 label="0048004900470048002000460069004F0032"/>
<unit class="alarm" ID="AlarmHistFiO2Low" type="BOOL" level="HIGH"
 label="004C004F0057002000460069004F0032"/>
<unit class="alarm" ID="AlarmHistIlvSlaveDisc" type="BOOL" level="HIGH"
 label="0049004C005600200044004900530043004F004E004E004500430054"/>
<unit class="alarm" ID="AlarmHistInop" type="BOOL" level="HIGH"
 label="00560045004E005400200049004E004F0050"/>
<unit class="alarm" ID="AlarmHistInvalidGasId" type="BOOL" level="MED"
 label="0049004E00560041004C004900440020004700410053002000490044"/>
<unit class="alarm" ID="AlarmHistLimitIE" type="BOOL" level="LOW"
 label="0049003A00450020004C0049004D00490054"/>
<unit class="alarm" ID="AlarmHistLimitTi" type="BOOL" level="LOW"
 label="004D0041005800200049004E00530050002000540049004D0045"/>
<unit class="alarm" ID="AlarmHistLimitVol" type="BOOL" level="LOW"
 label="0056004F004C0055004D00450020004C0049004D00490054"/>
<unit class="alarm" ID="AlarmHistLossAir" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004100490052"/>
<unit class="alarm" ID="AlarmHistLossGas" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004700410053"/>
<unit class="alarm" ID="AlarmHistLossHeliox" type="BOOL" level="HIGH"
 label="004C004F005300530020004F0046002000480045004C0049004F0058"/>
<unit class="alarm" ID="AlarmHistLossO2" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004F0032"/>
<unit class="alarm" ID="AlarmHistNcpapHigh" type="BOOL" level="HIGH"
 label="00480049004700480020006E004300500041005000200050005200450053"/>
<unit class="alarm" ID="AlarmHistNcpapHighPresLimit" type="BOOL" level="HIGH"
 label="006E0043005000410050002000500052004500530020004C0049004D0049005
4"/>

```



```

<unit class="alarm" ID="AlarmHistNcpapLow" type="BOOL" level="HIGH"
 label="004C004F00570020006E004300500041005000200050005200450053"/>
<unit class="alarm" ID="AlarmHistOcclusion" type="BOOL" level="HIGH"
 label="00430049005200430055004900540020004F00430043004C005500530049004F
 004E"/>
<unit class="alarm" ID="AlarmHistOpenSV" type="BOOL" level="HIGH"
 label="005300410046004500540059002000560041004C00560045"/>
<unit class="alarm" ID="AlarmHistPeepLow" type="BOOL" level="HIGH"
 label="004C004F005700200050004500450050"/>
<unit class="alarm" ID="AlarmHistPpeakHigh" type="BOOL" level="HIGH"
 label="004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmHistPpeakHighExt" type="BOOL" level="HIGH"
 label="0045005800540020004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmHistPpeakLow" type="BOOL" level="HIGH"
 label="004C004F005700200050007000650061006B"/>
<unit class="alarm" ID="AlarmHistPwrAcLoss" type="BOOL" level="HIGH"
 label="004C004F005300530020004F004600200041002F0043"/>
<unit class="alarm" ID="AlarmHistPwrBattLow" type="BOOL" level="HIGH"
 label="004C004F005700200042004100540054004500520059"/>
<unit class="alarm" ID="AlarmHistRateHigh" type="BOOL" level="MED"
 label="004800490047004800200052004100540045"/>
<unit class="alarm" ID="AlarmHistTest" type="BOOL" level="HIGH"
 label="0041004C00410052004D00200054004500530054"/>
<unit class="alarm" ID="AlarmHistVeHigh" type="BOOL" level="MED"
 label="0048004900470048002000560065"/>
<unit class="alarm" ID="AlarmHistVeLow" type="BOOL" level="HIGH"
 label="004C004F0057002000560065"/>
<unit class="alarm" ID="AlarmHistVteLow" type="BOOL" level="HIGH"
 label="004C004F00570020005600740065"/>
<unit class="alarm" ID="AlarmHistVtHigh" type="BOOL" level="LOW"
 label="00480049004700480020005600740065"/>
<unit class="scalar" ID="WaveAnlg0" type="WORD" scale="E+3" range="0000:2710"
 epoch="01F4" size="0032" units="" label="0041006E0061006C006F006700200030"/>
<unit class="scalar" ID="WaveAnlg1" type="WORD" scale="E+3" range="0000:2710"
 epoch="01F4" size="0032" units="" label="0041006E0061006C006F006700200031"/>
<unit class="scalar" ID="WaveFexp" type="WORD" scale="E+2" range="8AD0:7530"
 epoch="01F4" size="0032" units="004C002F006D0069006E"
 label="0046006500780070"/>
<unit class="scalar" ID="WaveFinsp" type="WORD" scale="E+2" range="8AD0:7530"
 epoch="01F4" size="0032" units="004C002F006D0069006E"
 label="00460069006E00730070"/>
<unit class="scalar" ID="WaveFlow" type="WORD" scale="E+2" range="8AD0:7530"
 epoch="01F4" size="0032" units="004C002F006D0069006E"
 label="0046006C006F0077"/>
<unit class="scalar" ID="WaveMetric" type="UWORD" range="0000:0000" epoch="01F4"
 size="0032" units="" label=""/>

```

```

<unit class="scalar" ID="WavePaw" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000610077"/>
<unit class="scalar" ID="WavePes" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000650073"/>
<unit class="scalar" ID="WavePinsp" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F"
 label="00500069006E00730070"/>
<unit class="scalar" ID="WavePtp" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000740070"/>
<unit class="scalar" ID="WavePtr" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000740072"/>
<unit class="scalar" ID="WaveVt" type="WORD" scale="E+1" range="D8F0:7530"
 epoch="01F4" size="0032" units="006D004C" label="00560074"/>
<unit class="info" ID="SysInfoConfig" type="UWORD"/>
<unit class="info" ID="SysInfoModel" type="ENUM"
 label="00420061007300650020004D006F00640065006C">
 <enum value="0000" label="0049006E00760061006C00690064"/>
 <enum value="0001" label="004100760065006100200043006F006D0070"/>
 <enum value="0002" label="00410076006500610020002B"/>
 <enum value="0003" label="0041007600650061"/>
</unit>
<unit class="info" ID="SysInfoOUI" type="TEXT"/>
<unit class="info" ID="SysInfoSerial" type="TEXT"/>
<unit class="info" ID="SysInfoSwVer" type="TEXT"/>
<unit class="info" ID="SysInfoTimeTotal" type="UINT" scale="E+2"
 range="00000000:0098967F"/>
<unit class="info" ID="PatInfoID" type="TEXT"
 label="004900440045004E00540049004600490043004100540049004F004E"/>
</profile>

```

## 12.2 Profile Message – AVEA Pediatric

```

<profile model="AVEA Comp" profileVersion="2.0" voxpVersion="3.0" textEncoding="UTF-16"
msgID="0002">
 <unit class="setting" ID="SetFiO2" type="WORD" resolution="0001" range="0015:0064"
 units="0025" label="00460069004F0032"/>
 <unit class="setting" ID="SetFlowBias" type="WORD" scale="E+1" resolution="0001"
 range="0004:0032" units="004C002F006D0069006E"
 label="004200690061007300200046006C006F0077"/>
 <unit class="setting" ID="SetFlowCycle" type="WORD" resolution="0005"
 range="0000:002D" units="0025"
 label="0046006C006F00770020004300790063006C0065"/>
 <unit class="setting" ID="SetFlowCyclePsv" type="WORD" resolution="0005"
 range="0005:002D" units="0025"
 label="0050005300560020004300790063006C0065"/>

```

```

<unit class="setting" ID="SetFlowDemand" type="WORD" resolution="0001"
 range="0000:0001" units="0030003D004F0066006600200031003D004F006E"
 label="00440065006D0061006E006400200046006C006F0077"/>
<unit class="setting" ID="SetFlowInsp" type="WORD" scale="E+1" resolution="0001"
 range="0001:004B" units="004C002F006D0069006E"
 label="005000650061006B00200046006C006F0077"/>
<unit class="setting" ID="SetPauseInsp" type="WORD" scale="E+2" resolution="0001"
 range="0000:012C" units="007300650063"
 label="0049006E00730070002000500061007500730065"/>
<unit class="setting" ID="SetPresHigh" type="WORD" resolution="0001"
 range="0000:005A" units="0063006D00480032004F"
 label="005000720065007300200048006900670068"/>
<unit class="setting" ID="SetPresInsp" type="WORD" resolution="0001"
 range="0000:005A" units="0063006D00480032004F"
 label="0049006E0073007000200050007200650073"/>
<unit class="setting" ID="SetPresLow" type="WORD" resolution="0001"
 range="0000:002D" units="0063006D00480032004F"
 label="00500072006500730020004C006F0077"/>
<unit class="setting" ID="SetPresNasalCPAP" type="WORD" resolution="0001"
 range="0002:000A" units="0063006D00480032004F"
 label="006E0043005000410050"/>
<unit class="setting" ID="SetPresPeep" type="WORD" resolution="0001"
 range="0000:0032" units="0063006D00480032004F" label="0050004500450050"/>
<unit class="setting" ID="SetPresPsv" type="WORD" resolution="0001"
 range="0000:005A" units="0063006D00480032004F" label="005000530056"/>
<unit class="setting" ID="SetRate" type="WORD" resolution="0001" range="0001:0096"
 units="00620070006D" label="0052006100740065"/>
<unit class="setting" ID="SetRiseInsp" type="WORD" resolution="0001"
 range="0001:0009" units="" label="0049006E0073007000200052006900730065"/>
<unit class="setting" ID="SetRisePsv" type="WORD" resolution="0001" range="0001:0009"
 units="" label="00500053005600200052006900730065"/>
<unit class="setting" ID="SetRiseVsync" type="WORD" resolution="0001"
 range="0001:0009" units=""
 label="005600730079006E006300200052006900730065"/>
<unit class="setting" ID="SetTimeHigh" type="WORD" scale="E+1" resolution="0001"
 range="0002:012C" units="007300650063"
 label="00540069006D006500200048006900670068"/>
<unit class="setting" ID="SetTimeHighPsv" type="WORD" resolution="0001"
 range="0000:0001" units="0030003D004F0066006600200031003D004F006E"
 label="0054002000480069006700680020005000530056"/>
<unit class="setting" ID="SetTimeHighSync" type="WORD" resolution="0001"
 range="0000:0032" units="0025"
 label="005400200048006900670068002000530079006E0063"/>
<unit class="setting" ID="SetTimeInsp" type="WORD" scale="E+2" resolution="0001"
 range="0014:01F4" units="007300650063"
 label="0049006E00730070002000540069006D0065"/>

```

```

<unit class="setting" ID="SetTimeLow" type="WORD" scale="E+1" resolution="0001"
 range="0002:012C" units="007300650063"
 label="00540069006D00650020004C006F0077"/>
<unit class="setting" ID="SetTimeLowSync" type="WORD" resolution="0001"
 range="0000:0032" units="0025"
 label="00540020004C006F0077002000530079006E0063"/>
<unit class="setting" ID="SetTmaxPsv" type="WORD" scale="E+2" resolution="0001"
 range="0014:01F4" units="007300650063"
 label="00500053005600200054006D00610078"/>
<unit class="setting" ID="SetTrigFlow" type="WORD" scale="E+1" resolution="0001"
 range="0001:00C8" units="004C002F006D0069006E"
 label="0046006C006F007700200054007200690067"/>
<unit class="setting" ID="SetTrigPres" type="WORD" scale="E+1" resolution="0001"
 range="0001:00C8" units="0063006D00480032004F"
 label="005000720065007300200054007200690067"/>
<unit class="setting" ID="SetVol" type="WORD" scale="E+1" resolution="0001"
 range="0019:01F4" units="006D004C" label="0056006F006C0075006D0065"/>
<unit class="setting" ID="SetVolAssured" type="WORD" scale="E+1" resolution="0001"
 range="0000:01F4" units="006D004C"
 label="004D00610063006800200056006F006C"/>
<unit class="setting" ID="SetVolLimit" type="WORD" scale="E+1" resolution="0001"
 range="0019:02EE" units="006D004C"
 label="0056006F006C0020004C0069006D00690074"/>
<unit class="setting" ID="SetVolSigh" type="WORD" resolution="0001" range="0000:0001"
 units="0030003D004F0066006600200031003D004F006E"
 label="0053006900670068"/>
<unit class="setting" ID="SetVolWave" type="ENUM"
 label="00570061007600650066006F0072006D">
 <enum value="0000" label="005300510055004100520045"/>
 <enum value="0001"
 label="0044004500430045004C00450052004100540049004E0047"/>
</unit>
<unit class="setting" ID="SetVsync" type="WORD" resolution="0001" range="0000:0001"
 units="0030003D004F0066006600200031003D004F006E"
 label="005600730079006E0063"/>
<unit class="setting" ID="LimitApnea" type="WORD" resolution="0001" range="0006:003C"
 units="007300650063"
 label="00410070006E0065006100200049006E00740065007200760061006C"/>
<unit class="setting" ID="LimitPeepLow" type="WORD" resolution="0001"
 range="0000:003C" units="0063006D00480032004F"
 label="004C006F007700200050004500450050"/>
<unit class="setting" ID="LimitPpeakHigh" type="WORD" resolution="0001"
 range="000A:0069" units="0063006D00480032004F"
 label="004800690067006800200050007000650061006B"/>

```

```

<unit class="setting" ID="LimitPpeakLow" type="WORD" resolution="0001"
 range="0001:0063" units="0063006D00480032004F"
 label="004C006F007700200050007000650061006B"/>
<unit class="setting" ID="LimitRateHigh" type="WORD" resolution="0001"
 range="0001:00C8" units="00620070006D"
 label="004800690067006800200052006100740065"/>
<unit class="setting" ID="LimitVeHigh" type="WORD" scale="E+2" resolution="0001"
 range="0000:012C" units="004C" label="0048006900670068002000560065"/>
<unit class="setting" ID="LimitVeLow" type="WORD" scale="E+2" resolution="0001"
 range="0000:012C" units="004C" label="004C006F0077002000560065"/>
<unit class="setting" ID="LimitVteHigh" type="WORD" scale="E+1" resolution="0001"
 range="0019:03E8" units="006D004C"
 label="00480069006700680020005600740065"/>
<unit class="setting" ID="LimitVteLow" type="WORD" scale="E+1" resolution="0001"
 range="0000:03E8" units="006D004C" label="004C006F00770020005600740065"/>
<unit class="setting" ID="SetAAC" type="BOOL" resolution="1" range="0:1" units=""
 label="004100410043000A004F006E"/>
<unit class="setting" ID="SetCircComp" type="WORD" scale="E+1" resolution="0001"
 range="0000:004B" units="006D004C002F0063006D00480032004F"
 label="004300690072006300200043006F006D0070"/>
<unit class="setting" ID="SetEttDia" type="WORD" scale="E+1" resolution="0001"
 range="0014:0064" units="006D006D"
 label="004400690061006D0065007400650072"/>
<unit class="setting" ID="SetEttLen" type="WORD" scale="E+1" resolution="0001"
 range="0014:0104" units="0063006D" label="004C0065006E006700740068"/>
<unit class="setting" ID="SetHumidifier" type="BOOL" resolution="1" range="0:1" units=""
 label="00480055004D0049004400490046004900450052"/>
<unit class="setting" ID="SetLanguage" type="ENUM"
 label="004C0061006E00670075006100670065003A">
 <enum value="0000" label="0045006E0067006C006900730068"/>
 <enum value="0001" label="004600720061006E00E7006100690073"/>
 <enum value="0002" label="0044006500750074007300630068"/>
 <enum value="0003" label="004900740061006C00690061006E006F"/>
 <enum value="0004" label="0050006F0072007400750067007500EA0073"/>
 <enum value="0005" label="004500730070006100F1006F006C"/>
 <enum value="0006" label="6C498BED"/>
 <enum value="0007" label="004E0065006400650072006C0061006E00640073"/>
</unit>
<unit class="setting" ID="SetLeakComp" type="BOOL" resolution="1" range="0:1" units=""
 label="004C00650061006B00200043006F006D0070000A004F006E"/>
<unit class="setting" ID="SetMode" type="ENUM"
 label="004D004F00440045002000530045004C004500430054">
 <enum value="0001"
 label="00410050005200560020002F00200042004900500048004100530049004
 3"/>

```

```

 <enum value="0002"
 label="00410050005200560020002F00200042004900500048004100530049004
 3"/>
 <enum value="0003" label="0050005200560043002000530049004D0056"/>
 <enum value="0004" label="005000520056004300200041002F0043"/>
 <enum value="0005" label="00430050004100500020002F0020005000530056"/>
 <enum value="0006" label="005400430050004C002000530049004D0056"/>
 <enum value="0007" label="005400430050004C00200041002F0043"/>
 <enum value="0008" label="00430050004100500020002F0020005000530056"/>
 <enum value="0009"
 label="00500052004500530053005500520045002000530049004D0056"/>
 <enum value="000A"
 label="0050005200450053005300550052004500200041002F0043"/>
 <enum value="000B" label="00430050004100500020002F0020005000530056"/>
 <enum value="000C"
 label="0056004F004C0055004D0045002000530049004D0056"/>
 <enum value="000D" label="0056004F004C0055004D004500200041002F0043"/>
 <enum value="000E" label="004E006100730061006C00200043005000410050"/>
 </unit>
<unit class="setting" ID="SetModellv" type="ENUM"
 label="0049004C00560020004D006F00640065003A">
 <enum value="0000" label="004F00660066"/>
 <enum value="0001" label="004D00610073007400650072"/>
 <enum value="0002" label="0053006C006100760065"/>
</unit>
<unit class="setting" ID="SetPatSize" type="ENUM"
 label="00500041005400490045004E0054002000530049005A0045002000530045004
 C004500430054">
 <enum value="0000" label="004E0065006F"/>
 <enum value="0001" label="005000650064"/>
 <enum value="0002" label="004100640075006C0074"/>
</unit>
<unit class="setting" ID="SetPatWt" type="WORD" scale="E+2" resolution="0001"
 range="0001:02EE" units="006B0067"
 label="005000740020005700650069006700680074"/>
<unit class="setting" ID="SetPresBaro" type="WORD" resolution="0001"
 range="0221:02F8" units="006D006D00480067"
 label="004200610072006F00200050007200650073"/>
<unit class="setting" ID="SetIncrFiO2" type="WORD" resolution="0001"
 range="0000:004F" units="0025"
 label="0049006E006300720065006100730065002000460069004F0032003A"/>
<unit class="setting" ID="SetSensitivityLowVte" type="WORD" resolution="0001"
 range="0001:0005" units=""
 label="004C006F0077002000560074006500200041006C00610072006D003A"/>
<unit class="monitor" ID="MntrAutoPEEP" type="WORD" range="0000:0032"
 units="0063006D00480032004F" label="004100750074006F0050004500450050"/>

```

```

<unit class="monitor" ID="MntrAutoPEEPdelta" type="WORD" range="0000:0032"
 units="0063006D00480032004F"
 label="0064004100750074006F0050004500450050"/>
<unit class="monitor" ID="MntrAutoPEEPesoph" type="WORD" range="0000:0032"
 units="0063006D00480032004F"
 label="004100750074006F005000450045005000650073"/>
<unit class="monitor" ID="MntrC20" type="WORD" scale="E+2" range="0000:01F4"
 units="" label="004300320030002F0043"/>
<unit class="monitor" ID="MntrCcw" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="004300630077"/>
<unit class="monitor" ID="MntrCdyn" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="004300640079006E"/>
<unit class="monitor" ID="MntrCdynNorm" type="WORD" scale="E+2" range="0000:01F4"
 units="006D004C002F0063006D00480032004F002F006B0067"
 label="004300640079006E002F006B0067"/>
<unit class="monitor" ID="MntrClung" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="0043006C0075006E0067"/>
<unit class="monitor" ID="MntrCstat" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="00430073007400610074"/>
<unit class="monitor" ID="MntrCstatNorm" type="WORD" scale="E+2" range="0000:01F4"
 units="006D004C002F0063006D00480032004F002F006B0067"
 label="00430073007400610074002F006B0067"/>
<unit class="monitor" ID="MntrFiO2" type="WORD" range="0000:0064" units="0025"
 label="00460069004F0032"/>
<unit class="monitor" ID="MntrIE" type="WORD" scale="E+1" range="FC19:03E7" units=""
 label="0049003A0045"/>
<unit class="monitor" ID="MntrLeak" type="WORD" range="0000:0064" units="0025"
 label="004C00650061006B"/>
<unit class="monitor" ID="MntrMIP" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="004D00490050"/>
<unit class="monitor" ID="MntrNcpapMeanFlow" type="WORD" scale="E+1"
 range="0000:0BB8" units="004C002F006D0069006E"
 label="004300500041005000200046006C006F0077"/>
<unit class="monitor" ID="MntrNcpapPres" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="006E0043005000410050"/>
<unit class="monitor" ID="MntrP100" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="0050003100300030"/>
<unit class="monitor" ID="MntrPair" type="WORD" range="0000:0050"
 units="0070007300690067" label="00410069007200200049006E006C00650074"/>
<unit class="monitor" ID="MntrPawDelta" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0064005000610077"/>
<unit class="monitor" ID="MntrPeep" type="WORD" range="0000:0032"
 units="0063006D00480032004F" label="0050004500450050"/>
<unit class="monitor" ID="MntrPefr" type="WORD" scale="E+1" range="0000:0BB8"
 units="004C002F006D0069006E" label="0050004500460052"/>

```

```

<unit class="monitor" ID="MntrPesDelta" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0064005000650073"/>
<unit class="monitor" ID="MntrPifr" type="WORD" scale="E+1" range="0000:0BB8"
 units="004C002F006D0069006E" label="0050004900460052"/>
<unit class="monitor" ID="MntrPmean" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0050006D00650061006E"/>
<unit class="monitor" ID="MntrPO2" type="WORD" range="0000:0050"
 units="0070007300690067" label="004F003200200049006E006C00650074"/>
<unit class="monitor" ID="MntrPpeak" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0050007000650061006B"/>
<unit class="monitor" ID="MntrPplat" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="00500070006C00610074"/>
<unit class="monitor" ID="MntrPplatPtp" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="00500074007000200050006C00610074"/>
<unit class="monitor" ID="MntrPtpPEEP" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="00500074007000200050004500450050"/>
<unit class="monitor" ID="MntrRate" type="WORD" range="0000:00C8"
 units="00620070006D" label="0052006100740065"/>
<unit class="monitor" ID="MntrRateMand" type="WORD" range="0000:00C8"
 units="00620070006D" label="004D0061006E006400200052006100740065"/>
<unit class="monitor" ID="MntrRateSpon" type="WORD" range="0000:00C8"
 units="00620070006D" label="00530070006F006E00200052006100740065"/>
<unit class="monitor" ID="MntrRimp" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="00520069006D0070"/>
<unit class="monitor" ID="MntrRlung" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="0052006C0075006E0067"/>
<unit class="monitor" ID="MntrRpeak" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="0052007000650061006B"/>
<unit class="monitor" ID="MntrRrs" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="005200720073"/>
<unit class="monitor" ID="MntrRSBIndex" type="WORD" range="0000:01F4"
 units="00620032002F006D0069006E002F004C" label="0066002F00560074"/>
<unit class="monitor" ID="MntrTe" type="WORD" scale="E+2" range="0000:270F"
 units="007300650063" label="00540065"/>
<unit class="monitor" ID="MntrTi" type="WORD" scale="E+2" range="0000:270F"
 units="007300650063" label="00540069"/>
<unit class="monitor" ID="MntrVdel" type="INT" scale="E+5" range="00000000:3B8B87C0"
 units="006D004C" label="005600640065006C"/>
<unit class="monitor" ID="MntrVeSpon" type="WORD" scale="E+2" range="0000:2706"
 units="004C" label="00530070006F006E002000560065"/>

```



```

<unit class="monitor" ID="MntrVeSponNorm" type="WORD" range="0000:03E7"
 units="006D004C002F006B0067"
 label="00530070006F006E002000560065002F006B0067"/>
<unit class="monitor" ID="MntrVeTotal" type="WORD" scale="E+2" range="0000:2706"
 units="004C" label="0054006F00740061006C002000560065"/>
<unit class="monitor" ID="MntrVeTotalNorm" type="WORD" range="0000:03E7"
 units="006D004C002F006B0067"
 label="0054006F00740061006C002000560065002F006B0067"/>
<unit class="monitor" ID="MntrVte" type="INT" scale="E+5" range="00000000:3B8B87C0"
 units="006D004C" label="005600740065"/>
<unit class="monitor" ID="MntrVteMand" type="INT" scale="E+5"
 range="00000000:3B8B87C0" units="006D004C"
 label="004D0061006E00640020005600740065"/>
<unit class="monitor" ID="MntrVteMandNorm" type="WORD" scale="E+2"
 range="0000:0BB8" units="006D004C002F006B0067"
 label="004D0061006E00640020005600740065002F006B0067"/>
<unit class="monitor" ID="MntrVteNorm" type="WORD" scale="E+2" range="0000:0BB8"
 units="006D004C002F006B0067" label="005600740065002F006B0067"/>
<unit class="monitor" ID="MntrVteSpon" type="INT" scale="E+5"
 range="00000000:3B8B87C0" units="006D004C"
 label="00530070006F006E0020005600740065"/>
<unit class="monitor" ID="MntrVteSponNorm" type="WORD" scale="E+2"
 range="0000:0BB8" units="006D004C002F006B0067"
 label="00530070006F006E0020005600740065002F006B0067"/>
<unit class="monitor" ID="MntrVti" type="INT" scale="E+5" range="00000000:3B8B87C0"
 units="006D004C" label="005600740069"/>
<unit class="monitor" ID="MntrVtiNorm" type="WORD" scale="E+2" range="0000:0BB8"
 units="006D004C002F006B0067" label="005600740069002F006B0067"/>
<unit class="monitor" ID="MntrWobImposed" type="WORD" scale="E+2"
 range="0000:07D0" units="006A006F0075006C00650073002F004C"
 label="0057004F00420069"/>
<unit class="monitor" ID="MntrWobPatient" type="WORD" scale="E+2" range="0000:07D0"
 units="006A006F0075006C00650073002F004C" label="0057004F00420070"/>
<unit class="monitor" ID="MntrWobVent" type="WORD" scale="E+2" range="0000:07D0"
 units="006A006F0075006C00650073002F004C" label="0057004F00420076"/>
<unit class="alarm" ID="AlarmActive" type="BOOL" label=""/>
<unit class="alarm" ID="AlarmActivePriority" type="ENUM" label="">
 <enum value="0001" label="0048004900470048"/>
 <enum value="0002" label="004D00450044"/>
 <enum value="0003" label="004C004F0057"/>
 <enum value="0004" label="0041004C004500520054"/>
</unit>
<unit class="alarm" ID="AlarmApnea" type="BOOL" level="HIGH"
 label="00410050004E0045004100200049004E00540045005200560041004C"/>

```

```

<unit class="alarm" ID="AlarmCircDisc" type="BOOL" level="HIGH"
 label="004300490052004300550049005400200044004900530043004F004E004E0045
 00430054"/>
<unit class="alarm" ID="AlarmFanFail" type="BOOL" level="LOW"
 label="00460041004E0020004600410049004C005500520045"/>
<unit class="alarm" ID="AlarmFiO2High" type="BOOL" level="HIGH"
 label="0048004900470048002000460069004F0032"/>
<unit class="alarm" ID="AlarmFiO2Low" type="BOOL" level="HIGH"
 label="004C004F0057002000460069004F0032"/>
<unit class="alarm" ID="AlarmIlvSlaveDisc" type="BOOL" level="HIGH"
 label="0049004C005600200044004900530043004F004E004E004500430054"/>
<unit class="alarm" ID="AlarmInop" type="BOOL" level="HIGH"
 label="00560045004E005400200049004E004F0050"/>
<unit class="alarm" ID="AlarmInvalidGasId" type="BOOL" level="MED"
 label="0049004E00560041004C004900440020004700410053002000490044"/>
<unit class="alarm" ID="AlarmLimitIE" type="BOOL" level="LOW"
 label="0049003A00450020004C0049004D00490054"/>
<unit class="alarm" ID="AlarmLimitTi" type="BOOL" level="LOW"
 label="004D0041005800200049004E00530050002000540049004D0045"/>
<unit class="alarm" ID="AlarmLimitVol" type="BOOL" level="LOW"
 label="0056004F004C0055004D00450020004C0049004D00490054"/>
<unit class="alarm" ID="AlarmLossAir" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004100490052"/>
<unit class="alarm" ID="AlarmLossGas" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004700410053"/>
<unit class="alarm" ID="AlarmLossHeliox" type="BOOL" level="HIGH"
 label="004C004F005300530020004F0046002000480045004C0049004F0058"/>
<unit class="alarm" ID="AlarmLossO2" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004F0032"/>
<unit class="alarm" ID="AlarmNcpapHigh" type="BOOL" level="HIGH"
 label="00480049004700480020006E004300500041005000200050005200450053"/>
<unit class="alarm" ID="AlarmNcpapHighPresLimit" type="BOOL" level="HIGH"
 label="006E0043005000410050002000500052004500530020004C0049004D0049005
 4"/>
<unit class="alarm" ID="AlarmNcpapLow" type="BOOL" level="HIGH"
 label="004C004F00570020006E004300500041005000200050005200450053"/>
<unit class="alarm" ID="AlarmOcclusion" type="BOOL" level="HIGH"
 label="00430049005200430055004900540020004F00430043004C005500530049004F
 004E"/>
<unit class="alarm" ID="AlarmOpenSV" type="BOOL" level="HIGH"
 label="005300410046004500540059002000560041004C00560045"/>
<unit class="alarm" ID="AlarmPeepLow" type="BOOL" level="HIGH"
 label="004C004F005700200050004500450050"/>
<unit class="alarm" ID="AlarmPpeakHigh" type="BOOL" level="HIGH"
 label="004800490047004800200050007000650061006B"/>

```

```
<unit class="alarm" ID="AlarmPpeakHighExt" type="BOOL" level="HIGH"
 label="0045005800540020004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmPpeakLow" type="BOOL" level="HIGH"
 label="004C004F005700200050007000650061006B"/>
<unit class="alarm" ID="AlarmPwrAcLoss" type="BOOL" level="HIGH"
 label="004C004F005300530020004F004600200041002F0043"/>
<unit class="alarm" ID="AlarmPwrBattLow" type="BOOL" level="HIGH"
 label="004C004F005700200042004100540054004500520059"/>
<unit class="alarm" ID="AlarmRateHigh" type="BOOL" level="MED"
 label="004800490047004800200052004100540045"/>
<unit class="alarm" ID="AlarmSilence" type="BOOL" label=""/>
<unit class="alarm" ID="AlarmTest" type="BOOL" level="HIGH"
 label="0041004C00410052004D00200054004500530054"/>
<unit class="alarm" ID="AlarmVeHigh" type="BOOL" level="MED"
 label="0048004900470048002000560065"/>
<unit class="alarm" ID="AlarmVeLow" type="BOOL" level="HIGH"
 label="004C004F0057002000560065"/>
<unit class="alarm" ID="AlarmVteLow" type="BOOL" level="HIGH"
 label="004C004F00570020005600740065"/>
<unit class="alarm" ID="AlarmVtHigh" type="BOOL" level="LOW"
 label="00480049004700480020005600740065"/>
<unit class="alarm" ID="AlarmHistApnea" type="BOOL" level="HIGH"
 label="00410050004E0045004100200049004E00540045005200560041004C"/>
<unit class="alarm" ID="AlarmHistCircDisc" type="BOOL" level="HIGH"
 label="004300490052004300550049005400200044004900530043004F004E004E0045
 00430054"/>
<unit class="alarm" ID="AlarmHistFanFail" type="BOOL" level="LOW"
 label="00460041004E0020004600410049004C005500520045"/>
<unit class="alarm" ID="AlarmHistFiO2High" type="BOOL" level="HIGH"
 label="0048004900470048002000460069004F0032"/>
<unit class="alarm" ID="AlarmHistFiO2Low" type="BOOL" level="HIGH"
 label="004C004F0057002000460069004F0032"/>
<unit class="alarm" ID="AlarmHistIlvSlaveDisc" type="BOOL" level="HIGH"
 label="0049004C005600200044004900530043004F004E004E004500430054"/>
<unit class="alarm" ID="AlarmHistInop" type="BOOL" level="HIGH"
 label="00560045004E005400200049004E004F0050"/>
<unit class="alarm" ID="AlarmHistInvalidGasId" type="BOOL" level="MED"
 label="0049004E00560041004C004900440020004700410053002000490044"/>
<unit class="alarm" ID="AlarmHistLimitIE" type="BOOL" level="LOW"
 label="0049003A00450020004C0049004D00490054"/>
<unit class="alarm" ID="AlarmHistLimitTi" type="BOOL" level="LOW"
 label="004D0041005800200049004E00530050002000540049004D0045"/>
<unit class="alarm" ID="AlarmHistLimitVol" type="BOOL" level="LOW"
 label="0056004F004C0055004D00450020004C0049004D00490054"/>
<unit class="alarm" ID="AlarmHistLossAir" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004100490052"/>
```

```

<unit class="alarm" ID="AlarmHistLossGas" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004700410053"/>
<unit class="alarm" ID="AlarmHistLossHeliox" type="BOOL" level="HIGH"
 label="004C004F005300530020004F0046002000480045004C0049004F0058"/>
<unit class="alarm" ID="AlarmHistLossO2" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004F0032"/>
<unit class="alarm" ID="AlarmHistNcpapHigh" type="BOOL" level="HIGH"
 label="00480049004700480020006E004300500041005000200050005200450053"/>
<unit class="alarm" ID="AlarmHistNcpapHighPresLimit" type="BOOL" level="HIGH"
 label="006E0043005000410050002000500052004500530020004C0049004D0049005
 4"/>
<unit class="alarm" ID="AlarmHistNcpapLow" type="BOOL" level="HIGH"
 label="004C004F00570020006E004300500041005000200050005200450053"/>
<unit class="alarm" ID="AlarmHistOcclusion" type="BOOL" level="HIGH"
 label="00430049005200430055004900540020004F00430043004C005500530049004F
 004E"/>
<unit class="alarm" ID="AlarmHistOpenSV" type="BOOL" level="HIGH"
 label="005300410046004500540059002000560041004C00560045"/>
<unit class="alarm" ID="AlarmHistPeepLow" type="BOOL" level="HIGH"
 label="004C004F005700200050004500450050"/>
<unit class="alarm" ID="AlarmHistPpeakHigh" type="BOOL" level="HIGH"
 label="004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmHistPpeakHighExt" type="BOOL" level="HIGH"
 label="0045005800540020004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmHistPpeakLow" type="BOOL" level="HIGH"
 label="004C004F005700200050007000650061006B"/>
<unit class="alarm" ID="AlarmHistPwrAcLoss" type="BOOL" level="HIGH"
 label="004C004F005300530020004F004600200041002F0043"/>
<unit class="alarm" ID="AlarmHistPwrBattLow" type="BOOL" level="HIGH"
 label="004C004F005700200042004100540054004500520059"/>
<unit class="alarm" ID="AlarmHistRateHigh" type="BOOL" level="MED"
 label="004800490047004800200052004100540045"/>
<unit class="alarm" ID="AlarmHistTest" type="BOOL" level="HIGH"
 label="0041004C00410052004D00200054004500530054"/>
<unit class="alarm" ID="AlarmHistVeHigh" type="BOOL" level="MED"
 label="0048004900470048002000560065"/>
<unit class="alarm" ID="AlarmHistVeLow" type="BOOL" level="HIGH"
 label="004C004F0057002000560065"/>
<unit class="alarm" ID="AlarmHistVteLow" type="BOOL" level="HIGH"
 label="004C004F00570020005600740065"/>
<unit class="alarm" ID="AlarmHistVtHigh" type="BOOL" level="LOW"
 label="00480049004700480020005600740065"/>
<unit class="scalar" ID="WaveAnlg0" type="WORD" scale="E+3" range="0000:2710"
 epoch="01F4" size="0032" units="" label="0041006E0061006C006F006700200030"/>
<unit class="scalar" ID="WaveAnlg1" type="WORD" scale="E+3" range="0000:2710"
 epoch="01F4" size="0032" units="" label="0041006E0061006C006F006700200031"/>

```

```

<unit class="scalar" ID="WaveFexp" type="WORD" scale="E+2" range="8AD0:7530"
 epoch="01F4" size="0032" units="004C002F006D0069006E"
 label="0046006500780070"/>
<unit class="scalar" ID="WaveFinsp" type="WORD" scale="E+2" range="8AD0:7530"
 epoch="01F4" size="0032" units="004C002F006D0069006E"
 label="00460069006E00730070"/>
<unit class="scalar" ID="WaveFlow" type="WORD" scale="E+2" range="8AD0:7530"
 epoch="01F4" size="0032" units="004C002F006D0069006E"
 label="0046006C006F0077"/>
<unit class="scalar" ID="WaveMetric" type="UWORD" range="0000:0000" epoch="01F4"
 size="0032" units="" label=""/>
<unit class="scalar" ID="WavePaw" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000610077"/>
<unit class="scalar" ID="WavePes" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000650073"/>
<unit class="scalar" ID="WavePinsp" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F"
 label="00500069006E00730070"/>
<unit class="scalar" ID="WavePtp" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000740070"/>
<unit class="scalar" ID="WavePtr" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000740072"/>
<unit class="scalar" ID="WaveVt" type="WORD" scale="E+1" range="F63C:1D4C"
 epoch="01F4" size="0032" units="006D004C" label="00560074"/>
<unit class="info" ID="SysInfoConfig" type="UWORD"/>
<unit class="info" ID="SysInfoModel" type="ENUM"
 label="00420061007300650020004D006F00640065006C">
 <enum value="0000" label="0049006E00760061006C00690064"/>
 <enum value="0001" label="004100760065006100200043006F006D0070"/>
 <enum value="0002" label="00410076006500610020002B"/>
 <enum value="0003" label="0041007600650061"/>
</unit>
<unit class="info" ID="SysInfoOUI" type="TEXT"/>
<unit class="info" ID="SysInfoSerial" type="TEXT"/>
<unit class="info" ID="SysInfoSwVer" type="TEXT"/>
<unit class="info" ID="SysInfoTimeTotal" type="UINT" scale="E+2"
 range="00000000:0098967F"/>
<unit class="info" ID="PatInfold" type="TEXT"
 label="004900440045004E00540049004600490043004100540049004F004E"/>
</profile>

```

### 12.3 Profile Message – AVEA Neonate

```

<profile model="Avea Comp" profileVersion="2.0" voxpVersion="3.0" textEncoding="UTF-16"
 msgID="0003">
 <unit class="setting" ID="SetFiO2" type="WORD" resolution="0001" range="0015:0064"
 units="0025" label="00460069004F0032"/>

```

```

<unit class="setting" ID="SetFlowBias" type="WORD" scale="E+1" resolution="0001"
 range="0004:0032" units="004C002F006D0069006E"
 label="004200690061007300200046006C006F0077"/>
<unit class="setting" ID="SetFlowCycle" type="WORD" resolution="0005"
 range="0000:002D" units="0025"
 label="0046006C006F00770020004300790063006C0065"/>
<unit class="setting" ID="SetFlowCyclePsv" type="WORD" resolution="0005"
 range="0005:002D" units="0025"
 label="0050005300560020004300790063006C0065"/>
<unit class="setting" ID="SetFlowDemand" type="WORD" resolution="0001"
 range="0000:0001" units="0030003D004F0066006600200031003D004F006E"
 label="00440065006D0061006E006400200046006C006F0077"/>
<unit class="setting" ID="SetFlowInsp" type="WORD" scale="E+1" resolution="0001"
 range="0004:012C" units="004C002F006D0069006E"
 label="005000650061006B00200046006C006F0077"/>
<unit class="setting" ID="SetPauseInsp" type="WORD" scale="E+2" resolution="0001"
 range="0000:012C" units="007300650063"
 label="0049006E00730070002000500061007500730065"/>
<unit class="setting" ID="SetPresHigh" type="WORD" resolution="0001"
 range="0000:005A" units="0063006D00480032004F"
 label="005000720065007300200048006900670068"/>
<unit class="setting" ID="SetPresInsp" type="WORD" resolution="0001"
 range="0000:0050" units="0063006D00480032004F"
 label="0049006E0073007000200050007200650073"/>
<unit class="setting" ID="SetPresLow" type="WORD" resolution="0001"
 range="0000:002D" units="0063006D00480032004F"
 label="00500072006500730020004C006F0077"/>
<unit class="setting" ID="SetPresNasalCPAP" type="WORD" resolution="0001"
 range="0002:000A" units="0063006D00480032004F"
 label="006E0043005000410050"/>
<unit class="setting" ID="SetPresPeep" type="WORD" resolution="0001"
 range="0000:0032" units="0063006D00480032004F" label="0050004500450050"/>
<unit class="setting" ID="SetPresPsv" type="WORD" resolution="0001" range="0000:0050"
 units="0063006D00480032004F" label="005000530056"/>
<unit class="setting" ID="SetRate" type="WORD" resolution="0001" range="0001:0096"
 units="00620070006D" label="0052006100740065"/>
<unit class="setting" ID="SetRiseInsp" type="WORD" resolution="0001"
 range="0001:0009" units="" label="0049006E0073007000200052006900730065"/>
<unit class="setting" ID="SetRisePsv" type="WORD" resolution="0001" range="0001:0009"
 units="" label="00500053005600200052006900730065"/>
<unit class="setting" ID="SetRiseVsync" type="WORD" resolution="0001"
 range="0001:0009" units=""
 label="005600730079006E006300200052006900730065"/>
<unit class="setting" ID="SetTimeHigh" type="WORD" scale="E+1" resolution="0001"
 range="0002:012C" units="007300650063"
 label="00540069006D006500200048006900670068"/>

```

```

<unit class="setting" ID="SetTimeHighPsv" type="WORD" resolution="0001"
 range="0000:0001" units="0030003D004F0066006600200031003D004F006E"
 label="0054002000480069006700680020005000530056"/>
<unit class="setting" ID="SetTimeHighSync" type="WORD" resolution="0001"
 range="0000:0032" units="0025"
 label="005400200048006900670068002000530079006E0063"/>
<unit class="setting" ID="SetTimeInsp" type="WORD" scale="E+2" resolution="0001"
 range="000F:012C" units="007300650063"
 label="0049006E00730070002000540069006D0065"/>
<unit class="setting" ID="SetTimeLow" type="WORD" scale="E+1" resolution="0001"
 range="0002:012C" units="007300650063"
 label="00540069006D00650020004C006F0077"/>
<unit class="setting" ID="SetTimeLowSync" type="WORD" resolution="0001"
 range="0000:0032" units="0025"
 label="00540020004C006F0077002000530079006E0063"/>
<unit class="setting" ID="SetTmaxPsv" type="WORD" scale="E+2" resolution="0001"
 range="000F:0BB8" units="007300650063"
 label="00500053005600200054006D00610078"/>
<unit class="setting" ID="SetTrigFlow" type="WORD" scale="E+1" resolution="0001"
 range="0001:00C8" units="004C002F006D0069006E"
 label="0046006C006F007700200054007200690067"/>
<unit class="setting" ID="SetTrigPres" type="WORD" scale="E+1" resolution="0001"
 range="0001:00C8" units="0063006D00480032004F"
 label="005000720065007300200054007200690067"/>
<unit class="setting" ID="SetVol" type="WORD" scale="E+1" resolution="0001"
 range="0014:0BB8" units="006D004C" label="0056006F006C0075006D0065"/>
<unit class="setting" ID="SetVolAssured" type="WORD" scale="E+1" resolution="0001"
 range="0000:0BB8" units="006D004C"
 label="004D00610063006800200056006F006C"/>
<unit class="setting" ID="SetVolLimit" type="WORD" scale="E+1" resolution="0001"
 range="0019:02EE" units="006D004C"
 label="0056006F006C0020004C0069006D00690074"/>
<unit class="setting" ID="SetVolSigh" type="WORD" resolution="0001" range="0000:0001"
 units="0030003D004F0066006600200031003D004F006E"
 label="0053006900670068"/>
<unit class="setting" ID="SetVolWave" type="ENUM"
 label="00570061007600650066006F0072006D">
 <enum value="0000" label="005300510055004100520045"/>
 <enum value="0001"
 label="0044004500430045004C00450052004100540049004E0047"/>
</unit>
<unit class="setting" ID="SetVsync" type="WORD" resolution="0001" range="0000:0001"
 units="0030003D004F0066006600200031003D004F006E"
 label="005600730079006E0063"/>

```

```

<unit class="setting" ID="LimitApnea" type="WORD" resolution="0001" range="0006:003C"
 units="007300650063"
 label="00410070006E0065006100200049006E00740065007200760061006C"/>
<unit class="setting" ID="LimitPeepLow" type="WORD" resolution="0001"
 range="0000:003C" units="0063006D00480032004F"
 label="004C006F007700200050004500450050"/>
<unit class="setting" ID="LimitPpeakHigh" type="WORD" resolution="0001"
 range="000A:0055" units="0063006D00480032004F"
 label="004800690067006800200050007000650061006B"/>
<unit class="setting" ID="LimitPpeakLow" type="WORD" resolution="0001"
 range="0001:0050" units="0063006D00480032004F"
 label="004C006F007700200050007000650061006B"/>
<unit class="setting" ID="LimitRateHigh" type="WORD" resolution="0001"
 range="0001:00C8" units="00620070006D"
 label="004800690067006800200052006100740065"/>
<unit class="setting" ID="LimitVeHigh" type="WORD" scale="E+2" resolution="0001"
 range="0000:01F4" units="004C" label="0048006900670068002000560065"/>
<unit class="setting" ID="LimitVeLow" type="WORD" scale="E+2" resolution="0001"
 range="0000:01F4" units="004C" label="004C006F0077002000560065"/>
<unit class="setting" ID="LimitVteHigh" type="WORD" scale="E+1" resolution="0001"
 range="0002:0BB8" units="006D004C"
 label="00480069006700680020005600740065"/>
<unit class="setting" ID="LimitVteLow" type="WORD" scale="E+1" resolution="0001"
 range="0000:0BB8" units="006D004C" label="004C006F00770020005600740065"/>
<unit class="setting" ID="SetAAC" type="BOOL" resolution="1" range="0:1" units=""
 label="004100410043000A004F006E"/>
<unit class="setting" ID="SetCircComp" type="WORD" scale="E+1" resolution="0001"
 range="0000:004B" units="006D004C002F0063006D00480032004F"
 label="004300690072006300200043006F006D0070"/>
<unit class="setting" ID="SetEttDia" type="WORD" scale="E+1" resolution="0001"
 range="0014:0064" units="006D006D"
 label="004400690061006D0065007400650072"/>
<unit class="setting" ID="SetEttLen" type="WORD" scale="E+1" resolution="0001"
 range="0014:0096" units="0063006D" label="004C0065006E006700740068"/>
<unit class="setting" ID="SetHumidifier" type="BOOL" resolution="1" range="0:1" units=""
 label="00480055004D0049004400490046004900450052"/>
<unit class="setting" ID="SetLanguage" type="ENUM"
 label="004C0061006E00670075006100670065003A">
 <enum value="0000" label="0045006E0067006C006900730068"/>
 <enum value="0001" label="004600720061006E00E7006100690073"/>
 <enum value="0002" label="0044006500750074007300630068"/>
 <enum value="0003" label="004900740061006C00690061006E006F"/>
 <enum value="0004" label="0050006F0072007400750067007500EA0073"/>
 <enum value="0005" label="004500730070006100F1006F006C"/>
 <enum value="0006" label="6C498BED"/>
 <enum value="0007" label="004E0065006400650072006C0061006E00640073"/>

```



```

</unit>
<unit class="setting" ID="SetLeakComp" type="BOOL" resolution="1" range="0:1" units=""
 label="004C00650061006B00200043006F006D0070000A004F006E"/>
<unit class="setting" ID="SetMode" type="ENUM"
 label="004D004F00440045002000530045004C004500430054">
 <enum value="0001"
 label="00410050005200560020002F00200042004900500048004100530049004
 3"/>
 <enum value="0002"
 label="00410050005200560020002F00200042004900500048004100530049004
 3"/>
 <enum value="0003" label="0050005200560043002000530049004D0056"/>
 <enum value="0004" label="005000520056004300200041002F0043"/>
 <enum value="0005" label="00430050004100500020002F0020005000530056"/>
 <enum value="0006" label="005400430050004C002000530049004D0056"/>
 <enum value="0007" label="005400430050004C00200041002F0043"/>
 <enum value="0008" label="00430050004100500020002F0020005000530056"/>
 <enum value="0009"
 label="00500052004500530053005500520045002000530049004D0056"/>
 <enum value="000A"
 label="0050005200450053005300550052004500200041002F0043"/>
 <enum value="000B" label="00430050004100500020002F0020005000530056"/>
 <enum value="000C"
 label="0056004F004C0055004D0045002000530049004D0056"/>
 <enum value="000D" label="0056004F004C0055004D004500200041002F0043"/>
 <enum value="000E" label="004E006100730061006C00200043005000410050"/>
</unit>
<unit class="setting" ID="SetModellv" type="ENUM"
 label="0049004C00560020004D006F00640065003A">
 <enum value="0000" label="004F00660066"/>
 <enum value="0001" label="004D00610073007400650072"/>
 <enum value="0002" label="0053006C006100760065"/>
</unit>
<unit class="setting" ID="SetPatSize" type="ENUM"
 label="00500041005400490045004E0054002000530049005A0045002000530045004
 C004500430054">
 <enum value="0000" label="004E0065006F"/>
 <enum value="0001" label="005000650064"/>
 <enum value="0002" label="004100640075006C0074"/>
</unit>
<unit class="setting" ID="SetPatWt" type="WORD" scale="E+2" resolution="0001"
 range="000A:0640" units="006B0067"
 label="005000740020005700650069006700680074"/>
<unit class="setting" ID="SetPresBaro" type="WORD" resolution="0001"
 range="0221:02F8" units="006D006D00480067"
 label="004200610072006F00200050007200650073"/>

```

```

<unit class="setting" ID="SetIncrFiO2" type="WORD" resolution="0001"
 range="0000:004F" units="0025"
 label="0049006E006300720065006100730065002000460069004F0032003A"/>
<unit class="setting" ID="SetSensitivityLowVte" type="WORD" resolution="0001"
 range="0001:0005" units=""
 label="004C006F0077002000560074006500200041006C00610072006D003A"/>
<unit class="monitor" ID="MntrAutoPEEP" type="WORD" range="0000:0032"
 units="0063006D00480032004F" label="004100750074006F0050004500450050"/>
<unit class="monitor" ID="MntrAutoPEEPdelta" type="WORD" range="0000:0032"
 units="0063006D00480032004F"
 label="0064004100750074006F0050004500450050"/>
<unit class="monitor" ID="MntrAutoPEEPesoph" type="WORD" range="0000:0032"
 units="0063006D00480032004F"
 label="004100750074006F005000450045005000650073"/>
<unit class="monitor" ID="MntrC20" type="WORD" scale="E+2" range="0000:01F4"
 units="" label="004300320030002F0043"/>
<unit class="monitor" ID="MntrCcw" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="004300630077"/>
<unit class="monitor" ID="MntrCdyn" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="004300640079006E"/>
<unit class="monitor" ID="MntrCdynNorm" type="WORD" scale="E+2" range="0000:01F4"
 units="006D004C002F0063006D00480032004F002F006B0067"
 label="004300640079006E002F006B0067"/>
<unit class="monitor" ID="MntrClung" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="0043006C0075006E0067"/>
<unit class="monitor" ID="MntrCstat" type="WORD" scale="E+2" range="0000:7530"
 units="006D004C002F0063006D00480032004F" label="00430073007400610074"/>
<unit class="monitor" ID="MntrCstatNorm" type="WORD" scale="E+2" range="0000:01F4"
 units="006D004C002F0063006D00480032004F002F006B0067"
 label="00430073007400610074002F006B0067"/>
<unit class="monitor" ID="MntrFiO2" type="WORD" range="0000:0064" units="0025"
 label="00460069004F0032"/>
<unit class="monitor" ID="MntrIE" type="WORD" scale="E+1" range="FC19:03E7" units=""
 label="0049003A0045"/>
<unit class="monitor" ID="MntrLeak" type="WORD" range="0000:0064" units="0025"
 label="004C00650061006B"/>
<unit class="monitor" ID="MntrMIP" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="004D00490050"/>
<unit class="monitor" ID="MntrNcpapMeanFlow" type="WORD" scale="E+1"
 range="0000:0BB8" units="004C002F006D0069006E"
 label="004300500041005000200046006C006F0077"/>
<unit class="monitor" ID="MntrNcpapPres" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="006E0043005000410050"/>
<unit class="monitor" ID="MntrP100" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="0050003100300030"/>

```

```

<unit class="monitor" ID="MntrPair" type="WORD" range="0000:0050"
 units="0070007300690067" label="00410069007200200049006E006C00650074"/>
<unit class="monitor" ID="MntrPawDelta" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0064005000610077"/>
<unit class="monitor" ID="MntrPeep" type="WORD" range="0000:0032"
 units="0063006D00480032004F" label="0050004500450050"/>
<unit class="monitor" ID="MntrPefr" type="WORD" scale="E+1" range="0000:0BB8"
 units="004C002F006D0069006E" label="0050004500460052"/>
<unit class="monitor" ID="MntrPesDelta" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0064005000650073"/>
<unit class="monitor" ID="MntrPifr" type="WORD" scale="E+1" range="0000:0BB8"
 units="004C002F006D0069006E" label="0050004900460052"/>
<unit class="monitor" ID="MntrPmean" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0050006D00650061006E"/>
<unit class="monitor" ID="MntrPO2" type="WORD" range="0000:0050"
 units="0070007300690067" label="004F003200200049006E006C00650074"/>
<unit class="monitor" ID="MntrPpeak" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="0050007000650061006B"/>
<unit class="monitor" ID="MntrPplat" type="WORD" range="0000:0078"
 units="0063006D00480032004F" label="00500070006C00610074"/>
<unit class="monitor" ID="MntrPplatPtp" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="00500074007000200050006C00610074"/>
<unit class="monitor" ID="MntrPtpPEEP" type="WORD" range="FFC4:0078"
 units="0063006D00480032004F" label="00500074007000200050004500450050"/>
<unit class="monitor" ID="MntrRate" type="WORD" range="0000:00C8"
 units="00620070006D" label="0052006100740065"/>
<unit class="monitor" ID="MntrRateMand" type="WORD" range="0000:00C8"
 units="00620070006D" label="004D0061006E006400200052006100740065"/>
<unit class="monitor" ID="MntrRateSpon" type="WORD" range="0000:00C8"
 units="00620070006D" label="00530070006F006E00200052006100740065"/>
<unit class="monitor" ID="MntrRimp" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="00520069006D0070"/>
<unit class="monitor" ID="MntrRlung" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="0052006C0075006E0067"/>
<unit class="monitor" ID="MntrRpeak" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="0052007000650061006B"/>
<unit class="monitor" ID="MntrRrs" type="WORD" scale="E+1" range="0000:03E8"
 units="0063006D00480032004F002F004C002F005300650063"
 label="005200720073"/>
<unit class="monitor" ID="MntrRSBIndex" type="WORD" range="0000:01F4"
 units="00620032002F006D0069006E002F004C" label="0066002F00560074"/>
<unit class="monitor" ID="MntrTe" type="WORD" scale="E+2" range="0000:270F"
 units="007300650063" label="00540065"/>

```

```

<unit class="monitor" ID="MntrTi" type="WORD" scale="E+2" range="0000:270F"
 units="007300650063" label="00540069"/>
<unit class="monitor" ID="MntrVdel" type="INT" scale="E+5" range="00000000:3B8B87C0"
 units="006D004C" label="005600640065006C"/>
<unit class="monitor" ID="MntrVeSpon" type="WORD" scale="E+2" range="0000:2706"
 units="004C" label="00530070006F006E002000560065"/>
<unit class="monitor" ID="MntrVeSponNorm" type="WORD" range="0000:03E7"
 units="006D004C002F006B0067"
 label="00530070006F006E002000560065002F006B0067"/>
<unit class="monitor" ID="MntrVeTotal" type="WORD" scale="E+2" range="0000:2706"
 units="004C" label="0054006F00740061006C002000560065"/>
<unit class="monitor" ID="MntrVeTotalNorm" type="WORD" range="0000:03E7"
 units="006D004C002F006B0067"
 label="0054006F00740061006C002000560065002F006B0067"/>
<unit class="monitor" ID="MntrVte" type="INT" scale="E+5" range="00000000:3B8B87C0"
 units="006D004C" label="005600740065"/>
<unit class="monitor" ID="MntrVteMand" type="INT" scale="E+5"
 range="00000000:3B8B87C0" units="006D004C"
 label="004D0061006E00640020005600740065"/>
<unit class="monitor" ID="MntrVteMandNorm" type="WORD" scale="E+2"
 range="0000:0BB8" units="006D004C002F006B0067"
 label="004D0061006E00640020005600740065002F006B0067"/>
<unit class="monitor" ID="MntrVteNorm" type="WORD" scale="E+2" range="0000:0BB8"
 units="006D004C002F006B0067" label="005600740065002F006B0067"/>
<unit class="monitor" ID="MntrVteSpon" type="INT" scale="E+5"
 range="00000000:3B8B87C0" units="006D004C"
 label="00530070006F006E0020005600740065"/>
<unit class="monitor" ID="MntrVteSponNorm" type="WORD" scale="E+2"
 range="0000:0BB8" units="006D004C002F006B0067"
 label="00530070006F006E0020005600740065002F006B0067"/>
<unit class="monitor" ID="MntrVti" type="INT" scale="E+5" range="00000000:3B8B87C0"
 units="006D004C" label="005600740069"/>
<unit class="monitor" ID="MntrVtiNorm" type="WORD" scale="E+2" range="0000:0BB8"
 units="006D004C002F006B0067" label="005600740069002F006B0067"/>
<unit class="monitor" ID="MntrWobImposed" type="WORD" scale="E+2"
 range="0000:07D0" units="006A006F0075006C00650073002F004C"
 label="0057004F00420069"/>
<unit class="monitor" ID="MntrWobPatient" type="WORD" scale="E+2" range="0000:07D0"
 units="006A006F0075006C00650073002F004C" label="0057004F00420070"/>
<unit class="monitor" ID="MntrWobVent" type="WORD" scale="E+2" range="0000:07D0"
 units="006A006F0075006C00650073002F004C" label="0057004F00420076"/>
<unit class="alarm" ID="AlarmActive" type="BOOL" label=""/>
<unit class="alarm" ID="AlarmActivePriority" type="ENUM" label="">
 <enum value="0001" label="0048004900470048"/>
 <enum value="0002" label="004D00450044"/>
 <enum value="0003" label="004C004F0057"/>

```

```

 <enum value="0004" label="0041004C004500520054"/>
 </unit>
 <unit class="alarm" ID="AlarmApnea" type="BOOL" level="HIGH"
 label="00410050004E0045004100200049004E00540045005200560041004C"/>
 <unit class="alarm" ID="AlarmCircDisc" type="BOOL" level="HIGH"
 label="004300490052004300550049005400200044004900530043004F004E004E0045
 00430054"/>
 <unit class="alarm" ID="AlarmFanFail" type="BOOL" level="LOW"
 label="00460041004E0020004600410049004C005500520045"/>
 <unit class="alarm" ID="AlarmFiO2High" type="BOOL" level="HIGH"
 label="0048004900470048002000460069004F0032"/>
 <unit class="alarm" ID="AlarmFiO2Low" type="BOOL" level="HIGH"
 label="004C004F0057002000460069004F0032"/>
 <unit class="alarm" ID="AlarmIlvSlaveDisc" type="BOOL" level="HIGH"
 label="0049004C005600200044004900530043004F004E004E004500430054"/>
 <unit class="alarm" ID="AlarmInop" type="BOOL" level="HIGH"
 label="00560045004E005400200049004E004F0050"/>
 <unit class="alarm" ID="AlarmInvalidGasId" type="BOOL" level="MED"
 label="0049004E00560041004C004900440020004700410053002000490044"/>
 <unit class="alarm" ID="AlarmLimitIE" type="BOOL" level="LOW"
 label="0049003A00450020004C0049004D00490054"/>
 <unit class="alarm" ID="AlarmLimitTi" type="BOOL" level="LOW"
 label="004D0041005800200049004E00530050002000540049004D0045"/>
 <unit class="alarm" ID="AlarmLimitVol" type="BOOL" level="LOW"
 label="0056004F004C0055004D00450020004C0049004D00490054"/>
 <unit class="alarm" ID="AlarmLossAir" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004100490052"/>
 <unit class="alarm" ID="AlarmLossGas" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004700410053"/>
 <unit class="alarm" ID="AlarmLossHeliox" type="BOOL" level="HIGH"
 label="004C004F005300530020004F0046002000480045004C0049004F0058"/>
 <unit class="alarm" ID="AlarmLossO2" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004F0032"/>
 <unit class="alarm" ID="AlarmNcpapHigh" type="BOOL" level="HIGH"
 label="00480049004700480020006E004300500041005000200050005200450053"/>
 <unit class="alarm" ID="AlarmNcpapHighPresLimit" type="BOOL" level="HIGH"
 label="006E0043005000410050002000500052004500530020004C0049004D0049005
 4"/>
 <unit class="alarm" ID="AlarmNcpapLow" type="BOOL" level="HIGH"
 label="004C004F00570020006E004300500041005000200050005200450053"/>
 <unit class="alarm" ID="AlarmOcclusion" type="BOOL" level="HIGH"
 label="00430049005200430055004900540020004F00430043004C005500530049004F
 004E"/>
 <unit class="alarm" ID="AlarmOpenSV" type="BOOL" level="HIGH"
 label="005300410046004500540059002000560041004C00560045"/>

```

```

<unit class="alarm" ID="AlarmPeepLow" type="BOOL" level="HIGH"
 label="004C004F005700200050004500450050"/>
<unit class="alarm" ID="AlarmPpeakHigh" type="BOOL" level="HIGH"
 label="004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmPpeakHighExt" type="BOOL" level="HIGH"
 label="0045005800540020004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmPpeakLow" type="BOOL" level="HIGH"
 label="004C004F005700200050007000650061006B"/>
<unit class="alarm" ID="AlarmPwrAcLoss" type="BOOL" level="HIGH"
 label="004C004F005300530020004F004600200041002F0043"/>
<unit class="alarm" ID="AlarmPwrBattLow" type="BOOL" level="HIGH"
 label="004C004F005700200042004100540054004500520059"/>
<unit class="alarm" ID="AlarmRateHigh" type="BOOL" level="MED"
 label="004800490047004800200052004100540045"/>
<unit class="alarm" ID="AlarmSilence" type="BOOL" label=""/>
<unit class="alarm" ID="AlarmTest" type="BOOL" level="HIGH"
 label="0041004C00410052004D00200054004500530054"/>
<unit class="alarm" ID="AlarmVeHigh" type="BOOL" level="MED"
 label="0048004900470048002000560065"/>
<unit class="alarm" ID="AlarmVeLow" type="BOOL" level="HIGH"
 label="004C004F0057002000560065"/>
<unit class="alarm" ID="AlarmVteLow" type="BOOL" level="HIGH"
 label="004C004F00570020005600740065"/>
<unit class="alarm" ID="AlarmVtHigh" type="BOOL" level="LOW"
 label="00480049004700480020005600740065"/>
<unit class="alarm" ID="AlarmHistApnea" type="BOOL" level="HIGH"
 label="00410050004E0045004100200049004E00540045005200560041004C"/>
<unit class="alarm" ID="AlarmHistCircDisc" type="BOOL" level="HIGH"
 label="004300490052004300550049005400200044004900530043004F004E004E0045
 00430054"/>
<unit class="alarm" ID="AlarmHistFanFail" type="BOOL" level="LOW"
 label="00460041004E0020004600410049004C005500520045"/>
<unit class="alarm" ID="AlarmHistFiO2High" type="BOOL" level="HIGH"
 label="0048004900470048002000460069004F0032"/>
<unit class="alarm" ID="AlarmHistFiO2Low" type="BOOL" level="HIGH"
 label="004C004F0057002000460069004F0032"/>
<unit class="alarm" ID="AlarmHistIlvSlaveDisc" type="BOOL" level="HIGH"
 label="0049004C005600200044004900530043004F004E004E004500430054"/>
<unit class="alarm" ID="AlarmHistInop" type="BOOL" level="HIGH"
 label="00560045004E005400200049004E004F0050"/>
<unit class="alarm" ID="AlarmHistInvalidGasId" type="BOOL" level="MED"
 label="0049004E00560041004C004900440020004700410053002000490044"/>
<unit class="alarm" ID="AlarmHistLimitLE" type="BOOL" level="LOW"
 label="0049003A00450020004C0049004D00490054"/>
<unit class="alarm" ID="AlarmHistLimitTi" type="BOOL" level="LOW"
 label="004D0041005800200049004E00530050002000540049004D0045"/>

```

```
<unit class="alarm" ID="AlarmHistLimitVol" type="BOOL" level="LOW"
 label="0056004F004C0055004D00450020004C0049004D00490054"/>
<unit class="alarm" ID="AlarmHistLossAir" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004100490052"/>
<unit class="alarm" ID="AlarmHistLossGas" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004700410053"/>
<unit class="alarm" ID="AlarmHistLossHeliox" type="BOOL" level="HIGH"
 label="004C004F005300530020004F0046002000480045004C0049004F0058"/>
<unit class="alarm" ID="AlarmHistLossO2" type="BOOL" level="HIGH"
 label="004C004F005300530020004F00460020004F0032"/>
<unit class="alarm" ID="AlarmHistNcpapHigh" type="BOOL" level="HIGH"
 label="00480049004700480020006E004300500041005000200050005200450053"/>
<unit class="alarm" ID="AlarmHistNcpapHighPresLimit" type="BOOL" level="HIGH"
 label="006E0043005000410050002000500052004500530020004C0049004D0049005
 4"/>
<unit class="alarm" ID="AlarmHistNcpapLow" type="BOOL" level="HIGH"
 label="004C004F00570020006E004300500041005000200050005200450053"/>
<unit class="alarm" ID="AlarmHistOcclusion" type="BOOL" level="HIGH"
 label="00430049005200430055004900540020004F00430043004C005500530049004F
 004E"/>
<unit class="alarm" ID="AlarmHistOpenSV" type="BOOL" level="HIGH"
 label="005300410046004500540059002000560041004C00560045"/>
<unit class="alarm" ID="AlarmHistPeepLow" type="BOOL" level="HIGH"
 label="004C004F005700200050004500450050"/>
<unit class="alarm" ID="AlarmHistPpeakHigh" type="BOOL" level="HIGH"
 label="004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmHistPpeakHighExt" type="BOOL" level="HIGH"
 label="0045005800540020004800490047004800200050007000650061006B"/>
<unit class="alarm" ID="AlarmHistPpeakLow" type="BOOL" level="HIGH"
 label="004C004F005700200050007000650061006B"/>
<unit class="alarm" ID="AlarmHistPwrAcLoss" type="BOOL" level="HIGH"
 label="004C004F005300530020004F004600200041002F0043"/>
<unit class="alarm" ID="AlarmHistPwrBattLow" type="BOOL" level="HIGH"
 label="004C004F005700200042004100540054004500520059"/>
<unit class="alarm" ID="AlarmHistRateHigh" type="BOOL" level="MED"
 label="004800490047004800200052004100540045"/>
<unit class="alarm" ID="AlarmHistTest" type="BOOL" level="HIGH"
 label="0041004C00410052004D00200054004500530054"/>
<unit class="alarm" ID="AlarmHistVeHigh" type="BOOL" level="MED"
 label="0048004900470048002000560065"/>
<unit class="alarm" ID="AlarmHistVeLow" type="BOOL" level="HIGH"
 label="004C004F0057002000560065"/>
<unit class="alarm" ID="AlarmHistVteLow" type="BOOL" level="HIGH"
 label="004C004F00570020005600740065"/>
<unit class="alarm" ID="AlarmHistVtHigh" type="BOOL" level="LOW"
 label="00480049004700480020005600740065"/>
```

```

<unit class="scalar" ID="WaveAnlg0" type="WORD" scale="E+3" range="0000:2710"
 epoch="01F4" size="0032" units="" label="0041006E0061006C006F006700200030"/>
<unit class="scalar" ID="WaveAnlg1" type="WORD" scale="E+3" range="0000:2710"
 epoch="01F4" size="0032" units="" label="0041006E0061006C006F006700200031"/>
<unit class="scalar" ID="WaveFexp" type="WORD" scale="E+2" range="8AD0:7530"
 epoch="01F4" size="0032" units="004C002F006D0069006E"
 label="0046006500780070"/>
<unit class="scalar" ID="WaveFinsp" type="WORD" scale="E+2" range="8AD0:7530"
 epoch="01F4" size="0032" units="004C002F006D0069006E"
 label="00460069006E00730070"/>
<unit class="scalar" ID="WaveFlow" type="WORD" scale="E+2" range="8AD0:7530"
 epoch="01F4" size="0032" units="004C002F006D0069006E"
 label="0046006C006F0077"/>
<unit class="scalar" ID="WaveMetric" type="UWORD" range="0000:0000" epoch="01F4"
 size="0032" units="" label=""/>
<unit class="scalar" ID="WavePaw" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000610077"/>
<unit class="scalar" ID="WavePes" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000650073"/>
<unit class="scalar" ID="WavePinsp" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F"
 label="00500069006E00730070"/>
<unit class="scalar" ID="WavePtp" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000740070"/>
<unit class="scalar" ID="WavePtr" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="0063006D00480032004F" label="005000740072"/>
<unit class="scalar" ID="WaveVt" type="WORD" scale="E+1" range="F63C:1D4C"
 epoch="01F4" size="0032" units="006D004C" label="00560074"/>
<unit class="info" ID="SysInfoConfig" type="UWORD"/>
<unit class="info" ID="SysInfoModel" type="ENUM"
 label="00420061007300650020004D006F00640065006C">
 <enum value="0000" label="0049006E00760061006C00690064"/>
 <enum value="0001" label="004100760065006100200043006F006D0070"/>
 <enum value="0002" label="00410076006500610020002B"/>
 <enum value="0003" label="0041007600650061"/>
</unit>
<unit class="info" ID="SysInfoOUI" type="TEXT"/>
<unit class="info" ID="SysInfoSerial" type="TEXT"/>
<unit class="info" ID="SysInfoSwVer" type="TEXT"/>
<unit class="info" ID="SysInfoTimeTotal" type="UINT" scale="E+2"
 range="00000000:0098967F"/>
<unit class="info" ID="PatInfoID" type="TEXT"
 label="004900440045004E00540049004600490043004100540049004F004E"/>
</profile>

```



## 12.4 Profile Message – VELA

```

<profile model="Vela Comprehensive" profileVersion="1.0" voxpVersion="2.0"
 textEncoding="UTF-32" msgID="0001">
 <unit class="setting" ID="SetFiO2" type="WORD" resolution="0001" range="0015:0064"
 units="00000025" label="00000046000000690000004F00000032"/>
 <unit class="setting" ID="SetFiO2IncreaseActive" type="BOOL" resolution="1" range="0:1"
 units="000000300000003D0000004F000000660000006600000020000000310000003
 D0000004F0000006E"
 label="00000031000000300000003000000025000000200000004F00000032"/>
 <unit class="setting" ID="SetFlowBias" type="WORD" scale="E+1" resolution="0001"
 range="0064:00C8" units="0000004C0000002F0000006D000000690000006E"
 label="0000004200000069000000610000007300000020000000460000006C0000006F
 00000077"/>
 <unit class="setting" ID="SetFlowCycle" type="WORD" resolution="0005"
 range="0000:001E" units="00000025"
 label="000000500000004300000020000000460000006C0000006F0000007700000020
 0000004300000079000000630000006C00000065"/>
 <unit class="setting" ID="SetFlowCyclePsv" type="WORD" resolution="0005"
 range="0005:001E" units="00000025"
 label="000000500000005300000056000000200000004300000079000000630000006C
 00000065"/>
 <unit class="setting" ID="SetFlowInsp" type="WORD" scale="E+1" resolution="0001"
 range="0064:0578" units="0000004C0000002F0000006D000000690000006E"
 label="0000005000000065000000610000006B00000020000000460000006C0000006F
 00000077"/>
 <unit class="setting" ID="SetNebulizerActive" type="BOOL" resolution="1" range="0:1"
 units="000000300000003D0000004F000000660000006600000020000000310000003
 D0000004F0000006E" label="0000004E0000006500000062"/>
 <unit class="setting" ID="SetPanelLockActive" type="BOOL" resolution="1" range="0:1"
 units="000000300000003D0000004F000000660000006600000020000000310000003
 D0000004F0000006E" label=""/>
 <unit class="setting" ID="SetPauseInsp" type="WORD" scale="E+2" resolution="0001"
 range="0000:00C8" units="000000530000006500000063"
 label="000000490000006E000000730000007000000020000000500000006100000075
 0000007300000065"/>
 <unit class="setting" ID="SetPresHigh" type="WORD" resolution="0001"
 range="0000:003C" units="000000630000006D00000048000000320000004F"
 label="0000005000000072000000650000007300000020000000480000006900000067
 00000068"/>
 <unit class="setting" ID="SetPresInsp" type="WORD" resolution="0001"
 range="0001:0064" units="000000630000006D00000048000000320000004F"
 label="000000490000006E000000730000007000000020000000500000007200000065
 00000073"/>
 <unit class="setting" ID="SetPresInspNPPV" type="WORD" resolution="0001"
 range="0001:0028" units="000000630000006D00000048000000320000004F"

```

```

label="0000004E0000005000000050000000560000002000000050000000690000006E
0000007300000070"/>
<unit class="setting" ID="SetPresLow" type="WORD" resolution="0001"
range="0000:002D" units="000000630000006D00000048000000320000004F"
label="000000500000000720000006500000073000000200000004C0000006F00000077
"/>
<unit class="setting" ID="SetPresPeep" type="WORD" resolution="0001"
range="0000:0023" units="000000630000006D00000048000000320000004F"
label="000000500000000450000004500000050"/>
<unit class="setting" ID="SetPresPsv" type="WORD" resolution="0001"
range="0000:003C" units="000000630000006D00000048000000320000004F"
label="0000005000000005300000056"/>
<unit class="setting" ID="SetPresPsvNPPV" type="WORD" resolution="0001"
range="0000:0028" units="000000630000006D00000048000000320000004F"
label="0000004E00000050000000500000005600000020000000500000005300000056
"/>
<unit class="setting" ID="SetRate" type="WORD" resolution="0001" range="0002:0050"
units="00000062000000700000006D"
label="00000052000000610000007400000065"/>
<unit class="setting" ID="SetTimeHigh" type="WORD" scale="E+1" resolution="0001"
range="0003:012C" units="000000530000006500000063"
label="00000054000000690000006D0000006500000020000000480000006900000067
00000068"/>
<unit class="setting" ID="SetTimeHighPsv" type="BOOL" resolution="1" range="0:1"
units="000000300000003D0000004F000000660000006600000020000000310000003
D0000004F0000006E"
label="0000005400000020000000480000006900000067000000680000002000000050
0000005300000056"/>
<unit class="setting" ID="SetTimeHighSync" type="WORD" resolution="0005"
range="0000:0032" units="00000025"
label="0000005400000020000000480000006900000067000000680000002000000053
000000790000006E00000063"/>
<unit class="setting" ID="SetTimeInsp" type="WORD" scale="E+2" resolution="0001"
range="0003:03E8" units="000000530000006500000063"
label="000000490000006E00000073000000700000002000000054000000690000006D
00000065"/>
<unit class="setting" ID="SetTimeLow" type="WORD" scale="E+1" resolution="0001"
range="0003:012C" units="000000530000006500000063"
label="00000054000000690000006D00000065000000200000004C0000006F0000007
7"/>
<unit class="setting" ID="SetTimeLowSync" type="WORD" resolution="0005"
range="0000:0032" units="00000025"
label="00000054000000200000004C0000006F00000077000000200000005300000079
0000006E00000063"/>
<unit class="setting" ID="SetTmaxPsv" type="WORD" scale="E+2" resolution="0001"
range="001E:012C" units="000000530000006500000063"

```

```

 label="000000500000000530000005600000020000000540000006D0000006100000078
 "/>
<unit class="setting" ID="SetTrigFlow" type="WORD" scale="E+1" resolution="0001"
 range="000A:00C8" units="0000004C0000002F0000006D000000690000006E"
 label="000000460000006C0000006F0000007700000020000000540000007200000069
 00000067"/>
<unit class="setting" ID="SetVol" type="WORD" scale="E+1" resolution="0005"
 range="0032:07D0" units="0000006D0000006C"
 label="000000560000006F0000006C000000750000006D00000065"/>
<unit class="setting" ID="SetVolAssured" type="WORD" scale="E+1" resolution="0005"
 range="01F4:4E20" units="0000006D0000006C"
 label="0000004100000073000000730000007500000072000000650000006400000020
 000000560000006F0000006C"/>
<unit class="setting" ID="SetVolLimit" type="WORD" scale="E+1" resolution="0005"
 range="01F4:61A8" units="0000006D0000006C"
 label="000000560000006F0000006C000000200000004C000000690000006D0000006
 900000074"/>
<unit class="setting" ID="SetVolSigh" type="BOOL" resolution="1" range="0:1"
 units="0000003000000003D0000004F000000660000006600000020000000310000003
 D0000004F0000006E" label="00000053000000690000006700000068"/>
<unit class="setting" ID="SetVolWave" type="ENUM"
 label="00000057000000610000007600000065000000660000006F000000720000006D
 ">
 <enum value="0000"
 label="000000530000005100000055000000410000005200000045"/>
 <enum value="0001"
 label="000000440000004500000043000000450000004C000000450000005200000
 04100000054000000490000004E00000047"/>
</unit>
<unit class="setting" ID="SetVsync" type="BOOL" resolution="1" range="0:1"
 units="0000003000000003D0000004F000000660000006600000020000000310000003
 D0000004F0000006E" label="0000005600000073000000790000006E00000063"/>
<unit class="setting" ID="LimitApnea" type="WORD" resolution="0001" range="000A:003C"
 units="000000530000006500000063"
 label="00000041000000700000006E000000650000006100000020000000490000006E
 00000074000000650000007200000076000000610000006C"/>
<unit class="setting" ID="LimitPpeakHigh" type="WORD" resolution="0001"
 range="0005:0078" units="000000630000006D00000048000000320000004F"
 label="0000004800000069000000670000006800000020000000500000007000000065
 000000610000006B"/>
<unit class="setting" ID="LimitPpeakLow" type="WORD" resolution="0001"
 range="0002:003C" units="000000630000006D00000048000000320000004F"
 label="0000004C0000006F000000770000002000000050000000700000006500000061
 0000006B"/>
<unit class="setting" ID="LimitRateHigh" type="WORD" resolution="0001"
 range="0003:0096" units="00000062000000700000006D"

```

```

 label="0000004800000069000000670000006800000020000000520000006100000074
 00000065"/>
<unit class="setting" ID="LimitVeLow" type="WORD" scale="E+2" resolution="0001"
 range="0000:270F" units="0000004C"
 label="0000004C0000006F00000077000000200000005600000065"/>
<unit class="setting" ID="SetAltitude" type="WORD" resolution="0001"
 range="FED4:0BF4"
 units="0000004D0000006500000074000000650000007200000073"
 label="000000410000006C000000740000006900000074000000750000006400000065
 "/>
<unit class="setting" ID="SetFiO2Monitoring" type="ENUM"
 label="00000046000000690000004F00000032000000200000004D0000006F0000006
 E00000069000000740000006F000000720000000A000000450000006E000000610000
 00620000006C0000006500000064">
 <enum value="0000"
 label="00000046000000690000004F00000032000000200000004D0000006F00000
 06E00000069000000740000006F000000720000000A00000044000000690000007
 300000061000000620000006C0000006500000064"/>
 <enum value="0001"
 label="00000046000000690000004F00000032000000200000004D0000006F00000
 06E00000069000000740000006F000000720000000A000000450000006E0000006
 1000000620000006C0000006500000064"/>
 </unit>
<unit class="setting" ID="SetHumidifier" type="BOOL" resolution="1" range="0:1"
 units="0000003000000003D0000004F000000660000006600000020000000310000003
 D0000004F0000006E"
 label="00000048000000750000006D0000006900000064000000690000006600000069
 0000006500000072"/>
<unit class="setting" ID="SetLanguage" type="ENUM"
 label="0000004C000000610000006E0000006700000075000000610000006700000065
 ">
 <enum value="0000"
 label="00000044000000650000007500000074000000730000006300000068"/>
 <enum value="0001"
 label="000000450000006E000000670000006C000000690000007300000068"/>
 <enum value="0002"
 label="0000004300000068000000690000006E000000650000007300000065"/>
 <enum value="0003"
 label="00000045000000730000007000000061000000F10000006F0000006C"/>
 <enum value="0004"
 label="0000004600000072000000610000006E000000E7000000610000006900000
 073"/>
 <enum value="0005"
 label="000000500000006F0000006C000000730000006B00000069"/>
 <enum value="0006" label="0000004400000075000000740000006300000068"/>

```

```
<enum value="0007"
 label="000000480000000750000006E0000006700000061000000720000006900000
 0610000006E"/>
</unit>
<unit class="setting" ID="SetLeakComp" type="BOOL" resolution="1" range="0:1"
 units="0000003000000003D0000004F000000660000006600000020000000310000003
 D0000004F0000006E"
 label="0000004C0000006B00000020000000430000006F0000006D00000070"/>
<unit class="setting" ID="SetMode" type="ENUM"
 label="0000004D0000004F00000044000000450000002000000053000000450000004
 C000000450000004300000054">
 <enum value="0001"
 label="0000004E00000020000000500000002000000050000000200000005600000
 020000000200000002000000041000000200000002F0000002000000043"/>
 <enum value="0002"
 label="0000004E00000020000000500000002000000050000000200000005600000
 0200000002000000020000000530000002000000049000000200000004D00000020
 00000056"/>
 <enum value="0003"
 label="0000004E00000020000000500000002000000050000000200000005600000
 02000000020000000200000005000000020000000530000002000000056"/>
 <enum value="0004"
 label="0000004E00000020000000500000002000000050000000200000005600000
 02000000020000000200000005000000020000000530000002000000056"/>
 <enum value="0005"
 label="0000004100000050000000520000005600000020000000420000006900000
 0500000006800000061000000730000006900000063"/>
 <enum value="0006"
 label="0000004100000050000000520000005600000020000000420000006900000
 0500000006800000061000000730000006900000063"/>
 <enum value="0007"
 label="00000050000000052000000560000004300000020000000530000004900000
 04D00000056"/>
 <enum value="0008"
 label="00000050000000052000000560000004300000020000000410000002F00000
 043"/>
 <enum value="0009"
 label="0000004300000020000000500000002000000041000000200000005000000
 02000000020000000200000005000000020000000530000002000000056"/>
 <enum value="000A"
 label="00000050000000020000000520000002000000045000000200000005300000
 0200000005300000020000000550000002000000052000000200000004500000020
 00000020000000020000000530000002000000049000000200000004D00000020000
 00056"/>
 <enum value="000B"
 label="00000050000000020000000520000002000000045000000200000005300000
```

```

0200000005300000020000000550000002000000052000000200000004500000020
0000002000000002000000041000000200000002F0000002000000043"/>
<enum value="000C"
 label="00000043000000020000000500000002000000041000000200000005000000
02000000020000000200000005000000020000000530000002000000056"/>
<enum value="000D"
 label="000000560000000200000004F000000200000004C000000200000005500000
0200000004D00000020000000450000002000000020000000200000005300000020
000000490000000200000004D0000002000000056"/>
<enum value="000E"
 label="000000560000000200000004F000000200000004C000000200000005500000
0200000004D00000020000000450000002000000020000000200000004100000020
0000002F00000002000000043"/>
</unit>
<unit class="setting" ID="SetNebulizerTime" type="WORD" resolution="0001"
 range="0001:003C"
 units="0000004D0000000690000006E00000075000000740000006500000073"
 label="0000004E00000065000000620000002000000054000000690000006D00000065
"/>
<unit class="setting" ID="SetPanelLockEnable" type="ENUM"
 label="0000004C0000006F000000630000006B000000730000000A000000450000006
E00000061000000620000006C0000006500000064">
 <enum value="0000"
 label="0000004C0000006F000000630000006B000000730000000A0000004400000
0690000007300000061000000620000006C0000006500000064"/>
 <enum value="0001"
 label="0000004C0000006F000000630000006B000000730000000A0000004500000
06E00000061000000620000006C0000006500000064"/>
</unit>
<unit class="setting" ID="SetVeLowOffEnable" type="ENUM"
 label="0000004C0000006F00000077000000200000004D000000690000006E0000002
0000000560000006F0000006C0000000A0000004F0000006600000066000000200000
00450000006E00000061000000620000006C0000006500000064">
 <enum value="0000"
 label="0000004C0000006F00000077000000200000004D000000690000006E00000
020000000560000006F0000006C0000000A0000004F00000066000000660000002
000000044000000690000007300000061000000620000006C0000006500000064"/
 >
 <enum value="0001"
 label="0000004C0000006F00000077000000200000004D000000690000006E00000
020000000560000006F0000006C0000000A0000004F00000066000000660000002
0000000450000006E00000061000000620000006C0000006500000064"/>
</unit>
<unit class="monitor" ID="MntrFiO2" type="WORD" range="0000:0064" units="00000025"
 label="000000460000000690000004F00000032"/>

```

```

<unit class="monitor" ID="MntrIE" type="WORD" scale="E+1" range="FC19:03E7" units=""
 label="000000490000003A00000045"/>
<unit class="monitor" ID="MntrPeep" type="WORD" range="0000:0063"
 units="000000630000006D00000048000000320000004F"
 label="00000050000000450000004500000050"/>
<unit class="monitor" ID="MntrPmean" type="WORD" range="0000:0063"
 units="000000630000006D00000048000000320000004F"
 label="000000500000006D00000065000000610000006E"/>
<unit class="monitor" ID="MntrPO2" type="WORD" range="0000:0050"
 units="000000700000000730000006900000067"
 label="0000004F0000003200000020000000490000006E0000006C0000006500000074
"/>
<unit class="monitor" ID="MntrPpeak" type="WORD" range="0000:008C"
 units="000000630000006D00000048000000320000004F"
 label="0000005000000007000000065000000610000006B"/>
<unit class="monitor" ID="MntrRate" type="WORD" range="0000:00FA"
 units="000000620000000700000006D"
 label="00000052000000610000007400000065"/>
<unit class="monitor" ID="MntrRateSpon" type="WORD" range="0000:00FA"
 units="000000620000000700000006D"
 label="000000530000000700000006F0000006E00000020000000520000006100000074
00000065"/>
<unit class="monitor" ID="MntrTe" type="WORD" scale="E+2" range="0000:03E7"
 units="0000005300000006500000063" label="0000005400000065"/>
<unit class="monitor" ID="MntrTi" type="WORD" scale="E+2" range="0000:03E7"
 units="0000005300000006500000063" label="0000005400000069"/>
<unit class="monitor" ID="MntrVeSpon" type="WORD" scale="E+1" range="0000:03E7"
 units="0000004C"
 label="000000530000000700000006F0000006E000000200000005600000065"/>
<unit class="monitor" ID="MntrVeTotal" type="WORD" scale="E+2" range="0000:03E7"
 units="0000004C" label="0000005600000065"/>
<unit class="monitor" ID="MntrVte" type="INT" scale="E+5" range="00000000:17D78400"
 units="0000006D0000006C" label="0000005600000007400000065"/>
<unit class="monitor" ID="MntrVteMand" type="INT" scale="E+5"
 range="00000000:17D78400" units="0000006D0000006C"
 label="0000004D000000610000006E00000064000000200000005600000074"/>
<unit class="monitor" ID="MntrVteSpon" type="INT" scale="E+5"
 range="00000000:17D78400" units="0000006D0000006C"
 label="000000530000000700000006F0000006E000000200000005600000074"/>
<unit class="monitor" ID="MntrVti" type="INT" scale="E+5" range="00000000:17D78400"
 units="0000006D0000006C" label="0000005600000007400000069"/>
<unit class="alarm" ID="AlarmSilence" type="BOOL" label=""/>
<unit class="alarm" ID="AlarmActive" type="BOOL" label=""/>
<unit class="alarm" ID="AlarmActivePriority" type="ENUM" label="">
 <enum value="0001" label="00000048000000490000004700000048"/>
 <enum value="0002" label="0000004D0000004500000044"/>

```

```

 <enum value="0003" label="0000004C0000004F00000057"/>
 <enum value="0004" label="000000410000004C000000450000005200000054"/>
</unit>
<unit class="alarm" ID="AlarmInop" type="BOOL" level="HIGH"
 label="00000056000000450000004E0000005400000020000000490000004E0000004F
 00000050"/>
<unit class="alarm" ID="AlarmMotorFault" type="BOOL" level="HIGH"
 label="0000004D0000006F000000740000006F00000072000000200000004600000061
 000000750000006C00000074"/>
<unit class="alarm" ID="AlarmApnea" type="BOOL" level="HIGH"
 label="00000041000000500000004E000000450000004100000020000000490000004E
 00000054000000450000005200000056000000410000004C"/>
<unit class="alarm" ID="AlarmHwFault" type="BOOL" level="HIGH"
 label="000000480000002F00000057000000200000004600000041000000550000004C
 00000054"/>
<unit class="alarm" ID="AlarmCircDisc" type="BOOL" level="HIGH"
 label="0000004300000049000000520000004300000055000000490000005400000020
 0000004600000041000000550000004C00000054"/>
<unit class="alarm" ID="AlarmPeepHigh" type="BOOL" level="HIGH"
 label="0000004800000049000000470000004800000020000000500000004500000045
 00000050"/>
<unit class="alarm" ID="AlarmPwrBattLow" type="BOOL" level="HIGH"
 label="0000004C0000004F000000570000002000000042000000410000005400000054
 000000450000005200000059"/>
<unit class="alarm" ID="AlarmLossO2" type="BOOL" level="HIGH"
 label="0000004F0000003200000020000000490000006E0000006C0000006500000074
 000000200000004C0000004F00000057"/>
<unit class="alarm" ID="AlarmPpeakHighExt" type="BOOL" level="HIGH"
 label="0000004800000049000000470000004800000020000000500000004900000050
 0000002C00000020000000530000005500000053000000540000002E"/>
<unit class="alarm" ID="AlarmPpeakHigh" type="BOOL" level="HIGH"
 label="0000004800000049000000470000004800000020000000500000004900000050
 "/>
<unit class="alarm" ID="AlarmPpeakLow" type="BOOL" level="HIGH"
 label="0000004C0000004F0000005700000020000000500000004900000050"/>
<unit class="alarm" ID="AlarmFiO2Range" type="BOOL" level="HIGH"
 label="000000250000004F000000320000002000000052000000410000004E00000047
 00000045000000200000004500000052000000520000004F00000052"/>
<unit class="alarm" ID="AlarmDefaults" type="BOOL" level="MED"
 label="00000044000000450000004600000041000000550000004C0000005400000053
 "/>
<unit class="alarm" ID="AlarmCheckEvents" type="BOOL" level="MED"
 label="000000430000004800000045000000430000004B000000200000004500000056
 000000450000004E0000005400000053"/>

```



```
<unit class="alarm" ID="AlarmPwrAcLoss" type="BOOL" level="MED"
 label="0000004F0000004E000000200000004200000041000000540000005400000045
 000000520000005900000020000000500000004F000000570000004500000052"/>
<unit class="alarm" ID="AlarmTransducerFault" type="BOOL" level="MED"
 label="0000005800000044000000430000005200000020000000460000004100000055
 0000004C00000054"/>
<unit class="alarm" ID="AlarmPwrBattMed" type="BOOL" level="MED"
 label="0000004D00000045000000440000002000000042000000410000005400000054
 000000450000005200000059"/>
<unit class="alarm" ID="AlarmO2PressHigh" type="BOOL" level="MED"
 label="0000004F0000003200000020000000490000004E0000004C0000004500000054
 0000002000000048000000490000004700000048"/>
<unit class="alarm" ID="AlarmFanFail" type="BOOL" level="MED"
 label="00000046000000410000004E000000200000004600000041000000490000004C
 000000550000005200000045"/>
<unit class="alarm" ID="AlarmRateHigh" type="BOOL" level="MED"
 label="0000004800000049000000470000004800000020000000520000004100000054
 00000045"/>
<unit class="alarm" ID="AlarmVeLow" type="BOOL" level="MED"
 label="0000004C0000004F00000057000000200000005600000065"/>
<unit class="alarm" ID="AlarmFiO2SensorFail" type="BOOL" level="MED"
 label="0000004F000000320000002000000053000000450000004E000000530000004F
 00000052000000200000004600000041000000490000004C000000550000005200000
 045"/>
<unit class="alarm" ID="AlarmFiO2Cal" type="BOOL" level="MED"
 label="000000430000004800000045000000430000004B000000200000004F00000032
 0000002000000043000000410000004C"/>
<unit class="alarm" ID="AlarmNoCalData" type="BOOL" level="LOW"
 label="0000004E0000004F0000002000000043000000410000004C0000002000000044
 000000410000005400000041"/>
<unit class="alarm" ID="AlarmInvalidSN" type="BOOL" level="LOW"
 label="000000490000004E00000056000000410000004C000000490000004400000020
 00000053000000450000005200000049000000410000004C000000200000004E00000
 0550000004D000000420000004500000052"/>
<unit class="alarm" ID="AlarmEEPROMFault" type="BOOL" level="ALERT"
 label="000000450000004500000050000000520000004F0000004D0000002000000046
 00000041000000550000004C00000054"/>
<unit class="alarm" ID="AlarmFlowSensorDisc" type="BOOL" level="ALERT"
 label="000000460000004C0000004F000000570000002000000053000000450000004E
 000000530000004F000000520000002000000044000000490000005300000043"/>
<unit class="alarm" ID="AlarmDirtyFilter" type="BOOL" level="ALERT"
 label="000000430000004800000045000000430000004B000000200000004600000049
 0000004C000000540000004500000052"/>
<unit class="alarm" ID="AlarmPatientDefaults" type="BOOL" level="LOW"
 label="00000050000000410000005400000049000000450000004E0000005400000020
 00000044000000450000004600000041000000550000004C0000005400000053"/>
```

```

<unit class="alarm" ID="AlarmClockBattLow" type="BOOL" level="MED"
 label="0000004C0000004F0000005700000020000000430000004C0000004F0000004
 30000004B000000200000004200000041000000540000005400000045000000520000
 0059"/>
<unit class="alarm" ID="AlarmHistInop" type="BOOL" level="HIGH"
 label="000000560000000450000004E0000005400000020000000490000004E0000004F
 00000050"/>
<unit class="alarm" ID="AlarmHistMotorFault" type="BOOL" level="HIGH"
 label="0000004D00000006F000000740000006F00000072000000200000004600000061
 0000007500000006C00000074"/>
<unit class="alarm" ID="AlarmHistApnea" type="BOOL" level="HIGH"
 label="000000410000005000000004E000000450000004100000020000000490000004E
 00000054000000450000005200000056000000410000004C"/>
<unit class="alarm" ID="AlarmHistHwFault" type="BOOL" level="HIGH"
 label="0000004800000002F00000057000000200000004600000041000000550000004C
 00000054"/>
<unit class="alarm" ID="AlarmHistCircDisc" type="BOOL" level="HIGH"
 label="0000004300000049000000520000004300000055000000490000005400000020
 0000004600000041000000550000004C00000054"/>
<unit class="alarm" ID="AlarmHistPeepHigh" type="BOOL" level="HIGH"
 label="0000004800000049000000470000004800000020000000500000004500000045
 00000050"/>
<unit class="alarm" ID="AlarmHistPwrBattLow" type="BOOL" level="HIGH"
 label="0000004C0000004F000000570000002000000042000000410000005400000054
 000000450000005200000059"/>
<unit class="alarm" ID="AlarmHistLossO2" type="BOOL" level="HIGH"
 label="0000004F0000003200000020000000490000006E0000006C0000006500000074
 000000200000004C0000004F00000057"/>
<unit class="alarm" ID="AlarmHistPpeakHighExt" type="BOOL" level="HIGH"
 label="0000004800000049000000470000004800000020000000500000004900000050
 0000002C00000020000000530000005500000053000000540000002E"/>
<unit class="alarm" ID="AlarmHistPpeakHigh" type="BOOL" level="HIGH"
 label="0000004800000049000000470000004800000020000000500000004900000050
 "/>
<unit class="alarm" ID="AlarmHistPpeakLow" type="BOOL" level="HIGH"
 label="0000004C0000004F0000005700000020000000500000004900000050"/>
<unit class="alarm" ID="AlarmHistFiO2Range" type="BOOL" level="HIGH"
 label="000000250000004F000000320000002000000052000000410000004E00000047
 00000045000000200000004500000052000000520000004F00000052"/>
<unit class="alarm" ID="AlarmHistDefaults" type="BOOL" level="MED"
 label="00000044000000450000004600000041000000550000004C0000005400000053
 "/>
<unit class="alarm" ID="AlarmHistCheckEvents" type="BOOL" level="MED"
 label="000000430000004800000045000000430000004B000000200000004500000056
 000000450000004E0000005400000053"/>

```

```
<unit class="alarm" ID="AlarmHistPwrAcLoss" type="BOOL" level="MED"
 label="0000004F0000004E000000200000004200000041000000540000005400000045
 000000520000005900000020000000500000004F000000570000004500000052"/>
<unit class="alarm" ID="AlarmHistTransducerFault" type="BOOL" level="MED"
 label="0000005800000044000000430000005200000020000000460000004100000055
 0000004C00000054"/>
<unit class="alarm" ID="AlarmHistPwrBattMed" type="BOOL" level="MED"
 label="0000004D00000045000000440000002000000042000000410000005400000054
 000000450000005200000059"/>
<unit class="alarm" ID="AlarmHistO2PressHigh" type="BOOL" level="MED"
 label="0000004F0000003200000020000000490000004E0000004C0000004500000054
 0000002000000048000000490000004700000048"/>
<unit class="alarm" ID="AlarmHistFanFail" type="BOOL" level="MED"
 label="00000046000000410000004E000000200000004600000041000000490000004C
 000000550000005200000045"/>
<unit class="alarm" ID="AlarmHistRateHigh" type="BOOL" level="MED"
 label="0000004800000049000000470000004800000020000000520000004100000054
 00000045"/>
<unit class="alarm" ID="AlarmHistVeLow" type="BOOL" level="MED"
 label="0000004C0000004F00000057000000200000005600000065"/>
<unit class="alarm" ID="AlarmHistFiO2SensorFail" type="BOOL" level="MED"
 label="0000004F000000320000002000000053000000450000004E000000530000004F
 00000052000000200000004600000041000000490000004C000000550000005200000
 045"/>
<unit class="alarm" ID="AlarmHistFiO2Cal" type="BOOL" level="MED"
 label="000000430000004800000045000000430000004B000000200000004F00000032
 0000002000000043000000410000004C"/>
<unit class="alarm" ID="AlarmHistNoCalData" type="BOOL" level="LOW"
 label="0000004E0000004F0000002000000043000000410000004C0000002000000044
 000000410000005400000041"/>
<unit class="alarm" ID="AlarmHistInvalidSN" type="BOOL" level="LOW"
 label="000000490000004E00000056000000410000004C000000490000004400000020
 00000053000000450000005200000049000000410000004C000000200000004E00000
 0550000004D000000420000004500000052"/>
<unit class="alarm" ID="AlarmHistEEPROMFault" type="BOOL" level="ALERT"
 label="000000450000004500000050000000520000004F0000004D0000002000000046
 00000041000000550000004C00000054"/>
<unit class="alarm" ID="AlarmHistFlowSensorDisc" type="BOOL" level="ALERT"
 label="000000460000004C0000004F000000570000002000000053000000450000004E
 000000530000004F000000520000002000000044000000490000005300000043"/>
<unit class="alarm" ID="AlarmHistDirtyFilter" type="BOOL" level="ALERT"
 label="000000430000004800000045000000430000004B000000200000004600000049
 0000004C000000540000004500000052"/>
<unit class="alarm" ID="AlarmHistPatientDefaults" type="BOOL" level="LOW"
 label="00000050000000410000005400000049000000450000004E0000005400000020
 00000044000000450000004600000041000000550000004C0000005400000053"/>
```

```

<unit class="alarm" ID="AlarmHistClockBattLow" type="BOOL" level="MED"
 label="0000004C0000004F0000005700000020000000430000004C0000004F0000004
 30000004B000000200000004200000041000000540000005400000045000000520000
 0059"/>
<unit class="scalar" ID="WaveFlow" type="WORD" scale="E+2" range="8AD0:7530"
 epoch="01F4" size="0032" units="0000004C00000002F0000006D000000690000006E"
 label="00000056000000020000000280000006C000000700000006D00000029"/>
<unit class="scalar" ID="WaveMetric" type="UWORD" range="0000:007F" epoch="01F4"
 size="0032" units="" label=""/>
<unit class="scalar" ID="WavePaw" type="WORD" scale="E+2" range="E890:2EE0"
 epoch="01F4" size="0032" units="000000630000006D00000048000000320000004F"
 label="000000500000006100000077"/>
<unit class="scalar" ID="WaveVt" type="WORD" scale="E+1" range="E4A8:5208"
 epoch="01F4" size="0032" units="0000006D0000006C"
 label="0000005600000074000000280000006D0000006C00000029"/>
<unit class="info" ID="SysInfoConfig" type="UWORD"/>
<unit class="info" ID="SysInfoDateTime" type="TEXT"/>
<unit class="info" ID="SysInfoModel" type="ENUM"
 label="00000056000000450000004E00000054000000200000004D0000004F00000044
 000000450000004C">
 <enum value="0000"
 label="00000056000000650000006C0000006100000020000000430000006F00000
 06D00000070000000720000006500000068000000650000006E000000730000006
 90000007600000065"/>
 <enum value="0001"
 label="00000056000000650000006C0000006100000020000000500000006C00000
 07500000073"/>
 <enum value="0002"
 label="00000056000000650000006C0000006100000020000000500000006C00000
 0750000007300000020000000490000006E0000007400000065000000720000006E
 0000006100000074000000690000006F0000006E000000610000006C"/>
 <enum value="0003"
 label="00000056000000650000006C0000006100000020000000420000006100000
 0730000006900000063"/>
</unit>
<unit class="info" ID="SysInfoOUI" type="TEXT"/>
<unit class="info" ID="SysInfoSerial" type="TEXT"/>
<unit class="info" ID="SysInfoSwVer" type="TEXT"/>
<unit class="info" ID="SysInfoSwVerBootLoader" type="TEXT"/>
<unit class="info" ID="SysInfoTimeTotal" type="UINT" scale="E+2"
 range="00000000:0098967F"
 label="00000056000000650000006E0000007400000020000000680000006F00000075
 0000007200000073"/>
<unit class="info" ID="SysInfoTurbineSerial" type="TEXT"/>

```

```

<unit class="info" ID="PatInfolD" type="TEXT"
 label="0000004900000044000000450000004E00000054000000490000004600000049
 000000430000004100000054000000490000004F0000004E"/>
</profile>

```

## 12.5 Config Message

```

<config mode="ACTIVE" msgID="0004">
 <unit class="setting" ID="SetFiO2"/>
 <unit class="setting" ID="SetFlowBias"/>
 <unit class="setting" ID="SetFlowCycle"/>
 <unit class="setting" ID="SetFlowCyclePsv"/>
 <unit class="setting" ID="SetFlowDemand"/>
 <unit class="setting" ID="SetFlowInsp"/>
 <unit class="setting" ID="SetPauseInsp"/>
 <unit class="setting" ID="SetPresHigh"/>
 <unit class="setting" ID="SetPresInsp"/>
 <unit class="setting" ID="SetPresLow"/>
 <unit class="setting" ID="SetPresNasalCPAP"/>
 <unit class="setting" ID="SetPresPeep"/>
 <unit class="setting" ID="SetPresPsv"/>
 <unit class="setting" ID="SetRate"/>
 <unit class="setting" ID="SetRiseInsp"/>
 <unit class="setting" ID="SetRisePsv"/>
 <unit class="setting" ID="SetRiseVsync"/>
 <unit class="setting" ID="SetTimeHigh"/>
 <unit class="setting" ID="SetTimeHighPsv"/>
 <unit class="setting" ID="SetTimeHighSync"/>
 <unit class="setting" ID="SetTimeInsp"/>
 <unit class="setting" ID="SetTimeLow"/>
 <unit class="setting" ID="SetTimeLowSync"/>
 <unit class="setting" ID="SetTmaxPsv"/>
 <unit class="setting" ID="SetTrigFlow"/>
 <unit class="setting" ID="SetTrigPres"/>
 <unit class="setting" ID="SetVol"/>
 <unit class="setting" ID="SetVolAssured"/>
 <unit class="setting" ID="SetVolLimit"/>
 <unit class="setting" ID="SetVolSigh"/>
 <unit class="setting" ID="SetVolWave"/>
 <unit class="setting" ID="SetVsync"/>
 <unit class="setting" ID="LimitApnea"/>
 <unit class="setting" ID="LimitPeepLow"/>
 <unit class="setting" ID="LimitPpeakHigh"/>
 <unit class="setting" ID="LimitPpeakLow"/>
 <unit class="setting" ID="LimitRateHigh"/>
 <unit class="setting" ID="LimitVeHigh"/>
 <unit class="setting" ID="LimitVeLow"/>

```

```
<unit class="setting" ID="LimitVteHigh"/>
<unit class="setting" ID="LimitVteLow"/>
<unit class="setting" ID="SetAAC"/>
<unit class="setting" ID="SetCircComp"/>
<unit class="setting" ID="SetEttDia"/>
<unit class="setting" ID="SetEttLen"/>
<unit class="setting" ID="SetHumidifier"/>
<unit class="setting" ID="SetLanguage"/>
<unit class="setting" ID="SetLeakComp"/>
<unit class="setting" ID="SetMode"/>
<unit class="setting" ID="SetModellv"/>
<unit class="setting" ID="SetPatSize"/>
<unit class="setting" ID="SetPatWt"/>
<unit class="setting" ID="SetPresBaro"/>
<unit class="setting" ID="SetIncrFiO2"/>
<unit class="setting" ID="SetSensitivityLowVte"/>
<unit class="monitor" ID="MntrAutoPEEP"/>
<unit class="monitor" ID="MntrAutoPEEPdelta"/>
<unit class="monitor" ID="MntrAutoPEEPesoph"/>
<unit class="monitor" ID="MntrC20"/>
<unit class="monitor" ID="MntrCcw"/>
<unit class="monitor" ID="MntrCdyn"/>
<unit class="monitor" ID="MntrCdynNorm"/>
<unit class="monitor" ID="MntrClung"/>
<unit class="monitor" ID="MntrCstat"/>
<unit class="monitor" ID="MntrCstatNorm"/>
<unit class="monitor" ID="MntrFiO2"/>
<unit class="monitor" ID="MntrIE"/>
<unit class="monitor" ID="MntrLeak"/>
<unit class="monitor" ID="MntrMIP"/>
<unit class="monitor" ID="MntrNcpapMeanFlow"/>
<unit class="monitor" ID="MntrNcpapPres"/>
<unit class="monitor" ID="MntrP100"/>
<unit class="monitor" ID="MntrPair"/>
<unit class="monitor" ID="MntrPawDelta"/>
<unit class="monitor" ID="MntrPeep"/>
<unit class="monitor" ID="MntrPefr"/>
<unit class="monitor" ID="MntrPesDelta"/>
<unit class="monitor" ID="MntrPifr"/>
<unit class="monitor" ID="MntrPmean"/>
<unit class="monitor" ID="MntrPO2"/>
<unit class="monitor" ID="MntrPpeak"/>
<unit class="monitor" ID="MntrPplat"/>
<unit class="monitor" ID="MntrPplatPtp"/>
<unit class="monitor" ID="MntrPtpPEEP"/>
<unit class="monitor" ID="MntrRate"/>
```

```
<unit class="monitor" ID="MntrRateMand"/>
<unit class="monitor" ID="MntrRateSpon"/>
<unit class="monitor" ID="MntrRimp"/>
<unit class="monitor" ID="MntrRlung"/>
<unit class="monitor" ID="MntrRpeak"/>
<unit class="monitor" ID="MntrRrs"/>
<unit class="monitor" ID="MntrRSBIndex"/>
<unit class="monitor" ID="MntrTe"/>
<unit class="monitor" ID="MntrTi"/>
<unit class="monitor" ID="MntrVdel"/>
<unit class="monitor" ID="MntrVeSpon"/>
<unit class="monitor" ID="MntrVeSponNorm"/>
<unit class="monitor" ID="MntrVeTotal"/>
<unit class="monitor" ID="MntrVeTotalNorm"/>
<unit class="monitor" ID="MntrVte"/>
<unit class="monitor" ID="MntrVteMand"/>
<unit class="monitor" ID="MntrVteMandNorm"/>
<unit class="monitor" ID="MntrVteNorm"/>
<unit class="monitor" ID="MntrVteSpon"/>
<unit class="monitor" ID="MntrVteSponNorm"/>
<unit class="monitor" ID="MntrVti"/>
<unit class="monitor" ID="MntrVtiNorm"/>
<unit class="monitor" ID="MntrWobImposed"/>
<unit class="monitor" ID="MntrWobPatient"/>
<unit class="monitor" ID="MntrWobVent"/>
<unit class="alarm" ID="AlarmActive"/>
<unit class="alarm" ID="AlarmActivePriority"/>
<unit class="alarm" ID="AlarmApnea"/>
<unit class="alarm" ID="AlarmCircDisc"/>
<unit class="alarm" ID="AlarmFanFail"/>
<unit class="alarm" ID="AlarmFiO2High"/>
<unit class="alarm" ID="AlarmFiO2Low"/>
<unit class="alarm" ID="AlarmIlvSlaveDisc"/>
<unit class="alarm" ID="AlarmInop"/>
<unit class="alarm" ID="AlarmInvalidGasId"/>
<unit class="alarm" ID="AlarmLimitIE"/>
<unit class="alarm" ID="AlarmLimitTi"/>
<unit class="alarm" ID="AlarmLimitVol"/>
<unit class="alarm" ID="AlarmLossAir"/>
<unit class="alarm" ID="AlarmLossGas"/>
<unit class="alarm" ID="AlarmLossHeliox"/>
<unit class="alarm" ID="AlarmLossO2"/>
<unit class="alarm" ID="AlarmNcpapHigh"/>
<unit class="alarm" ID="AlarmNcpapHighPresLimit"/>
<unit class="alarm" ID="AlarmNcpapLow"/>
<unit class="alarm" ID="AlarmOcclusion"/>
```

```
<unit class="alarm" ID="AlarmOpenSV"/>
<unit class="alarm" ID="AlarmPeepLow"/>
<unit class="alarm" ID="AlarmPpeakHigh"/>
<unit class="alarm" ID="AlarmPpeakHighExt"/>
<unit class="alarm" ID="AlarmPpeakLow"/>
<unit class="alarm" ID="AlarmPwrAcLoss"/>
<unit class="alarm" ID="AlarmPwrBattLow"/>
<unit class="alarm" ID="AlarmRateHigh"/>
<unit class="alarm" ID="AlarmSilence"/>
<unit class="alarm" ID="AlarmTest"/>
<unit class="alarm" ID="AlarmVeHigh"/>
<unit class="alarm" ID="AlarmVeLow"/>
<unit class="alarm" ID="AlarmVteLow"/>
<unit class="alarm" ID="AlarmVtHigh"/>
<unit class="alarm" ID="AlarmHistApnea"/>
<unit class="alarm" ID="AlarmHistCircDisc"/>
<unit class="alarm" ID="AlarmHistFanFail"/>
<unit class="alarm" ID="AlarmHistFiO2High"/>
<unit class="alarm" ID="AlarmHistFiO2Low"/>
<unit class="alarm" ID="AlarmHistIlvSlaveDisc"/>
<unit class="alarm" ID="AlarmHistInop"/>
<unit class="alarm" ID="AlarmHistInvalidGasId"/>
<unit class="alarm" ID="AlarmHistLimitE"/>
<unit class="alarm" ID="AlarmHistLimitTi"/>
<unit class="alarm" ID="AlarmHistLimitVol"/>
<unit class="alarm" ID="AlarmHistLossAir"/>
<unit class="alarm" ID="AlarmHistLossGas"/>
<unit class="alarm" ID="AlarmHistLossHeliox"/>
<unit class="alarm" ID="AlarmHistLossO2"/>
<unit class="alarm" ID="AlarmHistNcpapHigh"/>
<unit class="alarm" ID="AlarmHistNcpapHighPresLimit"/>
<unit class="alarm" ID="AlarmHistNcpapLow"/>
<unit class="alarm" ID="AlarmHistOcclusion"/>
<unit class="alarm" ID="AlarmHistOpenSV"/>
<unit class="alarm" ID="AlarmHistPeepLow"/>
<unit class="alarm" ID="AlarmHistPpeakHigh"/>
<unit class="alarm" ID="AlarmHistPpeakHighExt"/>
<unit class="alarm" ID="AlarmHistPpeakLow"/>
<unit class="alarm" ID="AlarmHistPwrAcLoss"/>
<unit class="alarm" ID="AlarmHistPwrBattLow"/>
<unit class="alarm" ID="AlarmHistRateHigh"/>
<unit class="alarm" ID="AlarmHistTest"/>
<unit class="alarm" ID="AlarmHistVeHigh"/>
<unit class="alarm" ID="AlarmHistVeLow"/>
<unit class="alarm" ID="AlarmHistVteLow"/>
<unit class="alarm" ID="AlarmHistVtHigh"/>
```



```
<unit class="scalar" ID="WaveAnlg0"/>
<unit class="scalar" ID="WaveAnlg1"/>
<unit class="scalar" ID="WaveFexp"/>
<unit class="scalar" ID="WaveFinsp"/>
<unit class="scalar" ID="WaveFlow"/>
<unit class="scalar" ID="WaveMetric"/>
<unit class="scalar" ID="WavePaw"/>
<unit class="scalar" ID="WavePes"/>
<unit class="scalar" ID="WavePinsp"/>
<unit class="scalar" ID="WavePtp"/>
<unit class="scalar" ID="WavePtr"/>
<unit class="scalar" ID="WaveVt"/>
<unit class="info" ID="SysInfoConfig"/>
<unit class="info" ID="SysInfoModel"/>
<unit class="info" ID="SysInfoOUI"/>
<unit class="info" ID="SysInfoSerial"/>
<unit class="info" ID="SysInfoSwVer"/>
<unit class="info" ID="SysInfoTimeTotal"/>
<unit class="info" ID="PatInfolD"/>
</config>
```

## 12.6 Data Message – AVEA

### Alarm message request from the Host to the AVEA:

&lt;link cmd="query" class="alarm" msgID="0007"/&gt;

Alarm message reply from AVEA to Host:

```
<data class="alarm" crc="059C"
 msgID="0007">00000000000000000000000000000000100000000001100000000
 010000</data>
```

Monitor message request from Host to AVEA:

&lt;link cmd="query" class="monitor" msgID="0005"/&gt;

Monitor message reply from AVEA to Host:

```
<data class="monitor" crc="3A99"
 msgID="0005">8000800080000067800005C580008000800080000015FF9D000A8
 000800080000029001E000604218000023F00080000002480008000800000080008
 00008000800000608000800002A9004502FCC25500000000013F02A43E5D02601
 70B80008000000000000000002F52A7C80008000800000098</data>
```

Scalar message requested from Host to AVEA:

&lt;link cmd="query" class="scalar" msgID="0006"/&gt;

&lt;data class="scalar" crc="B501" seq="0391"

[illegible]

### Setting message request from Host to AVEA:

&lt;link cmd="query" class="setting" msgID="0007"/&gt;

### Setting message reply from AVEA to Host:

```
<data class="setting" crc="67EA"
```

```
msgID="0007">0015001400000019000100C80000000F000F000600030006000000
0C0005000500050028000000000004B00140000004B000A001E03E8000000BB8000
```

0000000000014000300280008004B01F400320BB80000000000037009610000000  
0D000000000007402F8004F0003</data>

Info message request from Host to AVEA:

```
<link cmd="query" class="info" msgID="0004"/>
```

Info message reply from AVEA to Host:

<data class="info" crc="CA65"

```
msgID="0004">00FF00010030003000300037004100310030003000350039004200420034
00370030003000000004100420056003000310030003100370000000320000000260BF0000
</data>
```

## 12.7 Data Message - VELA

### Alarm message request from the Host to the VELA:

```
<link cmd="query" class="alarm" msgID="0003"/>
```

Alarm message reply from VELA to Host:

<data class="alarm" crc="E2E7"

[illegible]

Monitor message request from Host to VELA:

```
<link cmd="query" class="monitor" msgID="0004"/>
```

Monitor message reply from VELA to Host:

&lt;data class="monitor" crc="9B46"

```
msgID="0004">0014FFDF0000000580000016000C000001800073000001C202424
EE002424EE08000000002FAF080</data>
```

Scalar message requested from Host to VELA:

&lt;link cmd="query" class="scalar" msgID="0005"/&gt;

Scalar message reply from VELA to Host:

<data class="scalar" crc="AE76" seq="0094"

msgID="0005">0960094C0924090608E808C00898087A08480834080207DA07B2079E078A  
07620744058C0366FB46F60AF5ECF57EF3EEF394F40CF1DCF240F344F31CF2F4F394F40  
2F47AF48EF506F56AF592F5F6F63CF68CF6F0F72CF79AF7E0F844F89EF902F970F9B600  
3100310031003100310031003100310031003100310031003100310031003100310031003200320  
0320032003200320032003200320032003200320032003200320032003200320032003200320032  
00320032003200320032003200320032003200320032003207D607E907EF08020814081B08  
2D083A0840085908530866087F087F08780898089E0885081B06AA057104C904390384032  
602A9020D01D401AF017001510132012501250125011F0113010600F300F300F300DA00D4

```
00CE00C800BB00B500B500B500A80FF71046109410E1112D117711C1120A1251129712D
D1321136313A513E614261465148F148F1482142513C0135E12F112791208118A1106108C
10150F990F220EB00E410DD50D6B0D050CA10C410BE30B870B2F0AD90A860A3609E809
9E0957091308D3</data>
```

---

Setting message request from Host to VELA:

```
<link cmd="query" class="setting" msgID="0006"/>
```

Setting message reply from VELA to Host:

```
<data class="setting" crc="37C2"
 msgID="0006">00150006480010019015E008001000F000F000F000600008001800
 1000C002800000001E00140000012C00141388800161A800001000140028000300
 4B000500960001000010000E001E00000000</data>
```

---

Info message request from Host to VELA:

```
<link cmd="query" class="info" msgID="0007"/>
```

Info message reply from VELA to Host:

```
<data class="info" crc="6979"
 msgID="0007">FFC0000000320000003000000030000000340000002D000000300000003
 20000002D00000031000000300000005400000030000000340000003A000000350000003
 90000003A00000034000000300000002E0000003000000032000000340000000000000
 00030000000300000003000000037000000410000003100000030000000300000030000
 000300000004200000042000000420000003700000038000000430000000000000041000
 000450000005400000030000000320000003200000032000000300000000000000058000
 000300000002E00000030000000300000002E00000030000000300000000000000030000
 000310000002E00000030000000300000002E0000003300000030000000000000009C000
 000410000004500000030000000300000003200000034000000340000003300000000000
 00000</data>
```

## 12.8 Link Commands

Link Ping message request from the Host to AVEA:

```
<link cmd="ping" msgID="0002"/>
```

Link Ack message from AVEA to Host (i.e. for the above message – msgID="0002"):

```
<link cmd="ack" msgID="0002"/>
```

---

Link Restart message request from the Host to AVEA:

```
<link cmd="restart"/>
```

---

Link Nak message from the AVEA to Host:

```
<link cmd="nak" error="seq" msgID="0002"/>
```

---

```
<link cmd="nak" error="seq" msgID="0002"/>
```

---

Link Ping message request from the Host to VELA:

```
<link cmd="ping" msgID="0002"/>
```

Link Ack message from VELA to Host (i.e. for the above message – msgID="0002"):

```
<link cmd="ack" msgID="0002"/>
```

---

Link Restart message request from the Host to VELA:

```
<link cmd="restart"/>
```

---

Link Nak message from the VELA to Host:

```
<link cmd="nak" error="seq" msgID="0002"/>
```

```
<link cmd="nak" error="seq" msgID="0002"/>
```

## 13 Revision History

### 13.1 Rev A

This is the initial release of consolidating the VOXP Specification (91315 Rev D) with the AVEA VOXP Specialization, VELA VOXP Specialization, and the VOXP Specification Addendum (91316, 91415, and ER-1794, respectively). In addition, to Marketing's request to consolidate the various documents into a single document, additional parameters were added that will support end tidal CO<sub>2</sub>, closed loop FiO<sub>2</sub>, and SiPAP mode.

### 13.2 Rev B

10.4.2.9	Was:	"Label: "Vent hours""	
10.5.1.1	Was:	"0=Disable; 1=Enable"	
10.5.1.10	Was:	"cmH20"	
10.5.1.11	Was:	"cmH20"	
10.5.1.12	Was:	"cmH20"	
10.5.1.13	Was:	"cmH20"	
10.5.1.14	Was:	"cmH20"	
10.5.1.15	Was:	"cmH20"	
10.5.1.28	Was:	"cmH20"	
10.5.1.30	Was:	"Adult/Ped/Neo)"	
10.5.1.30	Inserted:	"Scale (Pediatric/Neo):	1"
10.5.1.30	Was:	"Resolution (Adult):	1"
10.5.1.30	Was:	"Range (Adult):	10 - 250"
10.5.1.30	Was:	"Units (Pediatric/Neo):	L"
10.5.1.31	Was:	"Scale (Adult):	2"

10.5.1.31	Was:	"Scale (Pediatric/Neo):	4"
10.5.1.31	Inserted:	"Resolution (Adult): 100"	
10.5.1.31	Was:	Range (Adult): 0 - 250	
10.5.1.31	Was:	"Units (Pediatric/Neo):	L"
10.5.1.32	Was:	"Scale (Adult): 2"	
10.5.1.32	Was:	"Scale (Pediatric/Neo):	4"
10.5.1.32	Was:	"Resolution (Adult): 1"	
10.5.1.32	Was:	Range (Adult): 10 - 250	
10.5.1.32	Was:	"Units (Pediatric/Neo):	L"
10.5.1.38	Was:	"Scale: 0; 1"	
10.5.1.38	Was:	"Range (Adult/Ped/Neo):	1-145;1-193"
10.5.1.42	Was:	"cmH20"	
10.5.1.43	Was:	"cmH20"	
10.5.1.44	Was:	"cmH20"	
10.5.1.52	Was:	"Scale (Adult): 2"	
10.5.1.52	Was:	"Scale (Pediatric/Neo):	4"
10.5.1.52	Was:	"Resolution (Adult): 1"	
10.5.1.52	Was:	"Range (Adult):	10 – 300"
10.5.1.52	Was:	"Units (Pediatric/Neo):	L"
10.5.1.53	Was:	"Scale (Adult): 2"	
10.5.1.53	Was:	"Scale (Pediatric/Neo):	4"
10.5.1.52	Was:	"Resolution (Adult): 1"	
10.5.1.53	Was:	"Range (Adult):	0 – 300"
10.5.1.53	Was:	"Units (Pediatric/Neo):	L"
10.5.1.54	Was:	"SetAAC 2.2.2.1.3"	
10.5.1.64	Was:	"Scale (Adult): 0"	
10.5.1.64	Was:	"Scale (Pediatric): 1"	
10.5.1.64	Was:	"Resolution: 1"	
10.5.1.64	Inserted:	"Resolution (Pediatric):	10"
10.5.1.64	Inserted:	"Resolution (Neonate):	1"
10.5.1.64	Was:	"Range (Adult):	1 – 300"
10.5.1.64	Was:	"Range (Pediatric):	10 – 7500"
10.5.1.72	Was:	"Enum value = label: 0=Maximum; 1=Normal; 2=APOD"	
10.5.2.6	Was:	"Range: 5 - 30"	
10.5.2.10	Was:	"Scale: 1"	
10.5.2.10	Was:	"Resolution: 1"	
10.5.2.10	Was:	"Range: 0 - 200"	
10.5.2.22	Was:	"Scale: 1"	
10.5.2.22	Was:	"Resolution: 1"	
10.5.2.22	Was:	"Range: 3 - 100"	
10.5.2.23	Was:	"Scale: 1"	
10.5.2.23	Was:	"Resolution: 1"	
10.5.2.23	Was:	"Range: 30 - 3000"	
10.5.2.25	Was:	"Range: 30 - 300"	
10.5.2.28	Was:	"Range: 50 - 2000"	
10.5.2.29	Was:	"Range: 0 - 20000"	

10.5.2.35	Was:	"6 – 150 (mmHg); 8 – 200 (kPa)"
10.5.2.40	Was:	"Scale: 1"
10.5.2.40	Was:	"Resolution: 1"
10.5.2.40	Was:	"0 - 999"
10.5.2.42	Was:	"Label: FiO2 Monitor (Enabled/Disabled)"
10.5.2.43	Was:	"Label: Humidifier (On/Off)"
10.5.2.43	Was:	"VELA GUI/Membrane: Humidifier (On/Off)"
10.5.2.45	Was:	"Label: Leak Comp (On/Off)"
10.5.2.45	Was:	"VELA GUI/Membrane: Leak Comp (On/Off)"
10.5.2.46	Was:	"1=NPPV A/C; 2=NPPV SIMV; 3=NPPVCPAP PSV-Volume; 4=NPPVCPAP PSV-Pressure; 5=APRV BiPhasic-Volume; 6=APRV BiPhasic-Pressure; 7=PRVC SIMV; 8=PRVC A/C; 9=CPAP PSV-Pressure; 10=PRESSURE SIMV; 11= PRESSURE A/C; 12=CPAP PSV-Volume; 13=VOLUME SIMV; 14=VOLUME A/C"
10.5.2.46	Was:	"Label: Mode"
10.5.2.46	Was:	"VELA GUI/Membrane: Mode"
10.5.2.49	Was:	"VELA GUI/Membrane: PANEL LOCK"
10.5.2.50	Was:	"Enum value = label: 0=Low Min Vol Off Disabled; 1=Low Min Vol Off Enabled"
10.6.1.1	Was:	"Resolution: 1"
10.6.1.2	Was:	"Resolution: 1"
10.6.1.3	Was:	"Resolution: 1"
10.6.1.4	Was:	"Resolution: 1"
10.6.1.5	Was:	"Scale: 2"
10.6.1.5	Was:	"Resolution: 100"
10.6.1.5	Was:	"Range: 30000"
10.6.1.6	Was:	"Scale: 2"
10.6.1.6	Was:	"Resolution (Adult): 100"
10.6.1.6	Was:	"Resolution (Pediatric): 100"
10.6.1.6	Was:	"Resolution (Neo): 1"
10.6.1.6	Was:	"Range: 0 - 30000"
10.6.1.6	Inserted:	"Range (Neo): 0 - 30000"
10.6.1.8	Was:	"Scale: 2"
10.6.1.8	Was:	"Resolution: 100"
10.6.1.8	Was:	"Range: 0 - 30000"
10.6.1.9	Was:	"Scale : 2"
10.6.1.9	Was:	"Resolution (Adult): 100"
10.6.1.9	Was:	"Resolution (Pediatric): 100"
10.6.1.9	Was:	"Resolution (Neo): 1"
10.6.1.9	Was:	"Range: 0 - 30000"
10.6.1.9	Inserted:	"Range (Neo): 0 - 30000"
10.6.1.11	Was:	"Scale : 1; 2"
10.6.1.11	Was:	"Resolution: 1"
10.6.1.11	Was:	"Range: 0 – 1500; 0 - 2000"
10.6.1.12	Was:	"Resolution: 1"
10.6.1.12	Was:	"Range: 0 – 100"

10.6.1.13	Was:	"Resolution:	1"
10.6.1.13	Was:	"Resolution:	1"
10.6.1.14	Was:	"Resolution:	1"
10.6.1.15	Was:	"Resolution:	1"
10.6.1.16	Was:	"Resolution:	1"
10.6.1.17	Was:	"Resolution:	1"
10.6.1.18	Was:	"Resolution:	1"
10.6.1.19	Was:	"Resolution:	1"
10.6.1.20	Was:	"Resolution:	1"
10.6.1.21	Was:	"Resolution:	1"
10.6.1.22	Was:	"Resolution:	1"
10.6.1.23	Was:	"Scale:	1"
10.6.1.23	Was:	"Resolution (Adult):	10"
10.6.1.23	Was:	"Resolution (Pediatric):	10"
10.6.1.23	Was:	"Resolution (Neo):	1"
10.6.1.23	Was:	"Range:	0 - 3000"
10.6.1.23	Inserted:	"Range (Adult/Ped):	0 - 300"
10.6.1.24	Was:	"Resolution:	1"
10.6.1.24	Was:	"Range:	should be -120 - 120"
10.6.1.25	Was:	"Scale:	1"
10.6.1.25	Was:	"Resolution (Adult):	10"
10.6.1.25	Was:	"Resolution (Pediatric):	10"
10.6.1.25	Was:	"Resolution (Neo):	1"
10.6.1.25	Was:	"Range:	0 - 3000"
10.6.1.25	Inserted:	"Range (Adult/Ped):	0 - 300"
10.6.1.26	Was:	"Resolution:	1"
10.6.1.27	Was:	"Resolution:	1"
10.6.1.28	Was:	"Resolution:	1"
10.6.1.29	Was:	"Resolution:	1"
10.6.1.30	Was:	"Resolution:	1"
10.6.1.31	Was:	"Resolution:	1"
10.6.1.31	Was:	"Resolution:	1"
10.6.1.31	Was:	"Scale:	0; 1"
10.6.1.31	Was:	"Range:	545 – 760; 727 – 1013"
10.6.1.32	Was:	"Resolution:	1"
10.6.1.33	Was:	"Resolution:	1"
10.6.1.34	Was:	"Resolution:	1"
10.6.1.35	Was:	"Resolution:	1"
10.6.1.36	Was:	"Resolution:	1"
10.6.1.37	Was:	"Resolution:	1"
10.6.1.38	Was:	"ID:	MntrRpeak)"
10.6.1.38	Was:	"Resolution:	1"
10.6.1.39	Was:	"Resolution:	1"
10.6.1.40	Was:	"Resolution:	1"
10.6.1.41	Was:	"Resolution:	1"
10.6.1.42	Was:	"Resolution:	1"



10.6.1.42	Was:	"Resolution:	1"
10.6.1.43	Was:	"Resolution:	1"
10.6.1.44	Was:	"Resolution:	1"
10.6.1.45	Was:	"Resolution:	1"
10.6.1.46	Was:	"Scale (Adult):	8"
10.6.1.46	Was:	"Scale (Pediatric/Neo):	8"
10.6.1.46	Inserted:	"Scale (Neo):	1"
10.6.1.46	Was:	"Resolution (Adult):	1000000"
10.6.1.46	Was:	"Resolution (Pediatric):	100000"
10.6.1.46	Was:	"Resolution (Neo):	10000"
10.6.1.46	Was:	"Range (Adult):	0 - 4000000000"
10.6.1.46	Was:	"Range (Pediatric):	0 - 1999000000"
10.6.1.46	Was:	"Range (Neo):	0 - 999900000"
10.6.1.46	Was:	"Units (Pediatric/Neo):	L"
10.6.1.47	Was:	"Resolution:	1"
10.6.1.48	Was:	"Resolution:	1"
10.6.1.49	Was:	"Resolution:	1"
10.6.1.50	Was:	"Scale:	2"
10.6.1.50	Was:	"Resolution (Adult/Ped):	10"
10.6.1.50	Was:	"Resolution (Neo):	1"
10.6.1.50	Inserted:	"Range (Adult/Ped):	0 - 999"
10.6.1.50	Was:	"Range:	0 - 9990"
10.6.1.51	Was:	"Resolution:	1"
10.6.1.52	Inserted:	"Scale (Adult/Ped):	2"
10.6.1.52	Was:	"Scale:	2"
10.6.1.52	Was:	"Resolution (Adult):	10"
10.6.1.52	Was:	"Resolution (Pediatric):	10"
10.6.1.52	Was:	"Resolution (Neo):	1"
10.6.1.52	Inserted:	"Range (Adult/Ped):	0 - 999"
10.6.1.52	Was:	"Range:	0 - 9990"
10.6.1.53	Was:	"Resolution:	1"
10.6.1.54	Was:	"Resolution:	1"
10.6.1.55	Was:	"Scale (Adult):	8"
10.6.1.55	Was:	"Scale (Pediatric/Neo):	8"
10.6.1.55	Inserted:	"Scale (Neo):	1"
10.6.1.55	Was:	"Resolution (Adult):	1000000"
10.6.1.55	Was:	"Resolution (Pediatric):	100000"
10.6.1.55	Was:	"Resolution (Neo):	10000"
10.6.1.55	Was:	"Range (Adult):	0 - 4000000000"
10.6.1.55	Was:	"Range (Pediatric):	0 - 1999000000"
10.6.1.55	Was:	"Range (Neo):	0 - 999000000"
10.6.1.55	Was:	"Units (Pediatric/Neo):	L"
10.6.1.56	Was:	"Scale (Adult):	8"
10.6.1.56	Was:	"Scale (Pediatric/Neo):	8"
10.6.1.56	Inserted:	"Scale (Neo):	1"
10.6.1.56	Was:	"Resolution (Adult):	1000000"

10.6.1.56	Was:	"Resolution (Pediatric):	100000"
10.6.1.56	Was:	"Resolution (Neo):	10000"
10.6.1.56	Was:	"Range (Adult):	0 - 4000000000"
10.6.1.56	Was:	"Range (Pediatric):	0 - 1999000000"
10.6.1.56	Was:	"Range (Neo):	0 - 999000000"
10.6.1.56	Was:	"Units (Pediatric/Neo):	L"
10.6.1.57	Was:	"Scale:	2"
10.6.1.57	Was:	"Resolution:	1"
10.6.1.58	Was:	"Scale:	2"
10.6.1.58	Was:	"Resolution:	1"
10.6.1.59	Was:	"Scale (Adult):	8"
10.6.1.59	Was:	"Scale (Pediatric/Neo):	8"
10.6.1.59	Inserted:	"Scale (Neo):	1"
10.6.1.59	Was:	"Resolution (Adult):	1000000"
10.6.1.59	Was:	"Resolution (Pediatric):	100000"
10.6.1.59	Was:	"Resolution (Neo):	10000"
10.6.1.59	Was:	"Range (Adult):	0 - 4000000000"
10.6.1.59	Was:	"Range (Pediatric):	0 - 1999000000"
10.6.1.59	Was:	"Range (Neo):	0 - 999000000"
10.6.1.59	Was:	"Units (Pediatric/Neo):	L"
10.6.1.60	Was:	"Scale:	2"
10.6.1.60	Was:	"Resolution:	1"
10.6.1.61	Was:	"Scale (Adult):	8"
10.6.1.61	Was:	"Scale (Pediatric/Neo):	8"
10.6.1.61	Inserted:	"Scale (Neo):	1"
10.6.1.61	Was:	"Resolution (Adult):	1000000"
10.6.1.61	Was:	"Resolution (Pediatric):	100000"
10.6.1.61	Was:	"Resolution (Neo):	10000"
10.6.1.61	Was:	"Range (Adult):	0 - 4000000000"
10.6.1.61	Was:	"Range (Pediatric):	0 - 1999000000"
10.6.1.61	Was:	"Range (Neo):	0 - 999000000"
10.6.1.61	Was:	"Units (Pediatric/Neo):	L"
10.6.1.62	Was:	"Scale:	2"
10.6.1.62	Was:	"Resolution:	1"
10.6.1.63	Was:	"Resolution:	1"
10.6.1.64	Was:	"Resolution:	1"
10.6.1.65	Was:	"Resolution:	1"
10.6.2.1	Was:	"Resolution:	1"
10.6.2.2	Was:	"Resolution:	1"
10.6.2.3	Was:	"Resolution:	1"
10.6.2.4	Was:	"Resolution:	1"
10.6.2.5	Was:	"Resolution:	1"
10.6.2.6	Was:	"Resolution:	1"
10.6.2.6	Was:	"Range:	40 - 85"
10.6.2.7	Was:	"Resolution:	1"
10.6.2.8	Was:	"Resolution:	1"

10.6.2.9	Was:	"Resolution:	1"
10.6.2.10	Was:	"Resolution:	1"
10.6.2.11	Was:	"Resolution:	1"
10.6.2.12	Was:	"Resolution:	1"
10.6.2.13	Was:	"Scale:	2"
10.6.2.13	Was:	"Range:	0 - 9990"
10.6.2.14	Was:	"Scale:	2"
10.6.2.14	Was:	"Resolution:	1"
10.6.2.14	Was:	"Range:	0 - 9990"
10.6.2.15	Was:	"Scale:	5"
10.6.2.15	Was:	"Resolution:	100000"
10.6.2.15	Was:	"Range:	0 - 4000000000"
10.6.2.16	Was:	"Scale:	5"
10.6.2.16	Was:	"Resolution:	100000"
10.6.2.16	Was:	"Range:	0 - 4000000000"
10.6.2.17	Was:	"Scale:	5"
10.6.2.17	Was:	"Resolution:	100000"
10.6.2.17	Was:	"Range:	0 - 4000000000"
10.6.2.18	Was:	"Resolution:	100000"
10.6.2.18	Was:	"Scale:	5"
10.6.2.18	Was:	"Range:	0 - 4000000000"
10.7.1.2	Was:	"Enum value = label: 1=HIGH; 2=MED; 3=LOW;"	
10.7.1.25	Was:	"LOW;"	
10.7.1.25	Was:	"VOL LIMIT"	
10.7.1.31	Was:	"Label:	nCPAP PRESSURE LIMIT"
10.7.1.31	Was:	"AVEA GUI/Membrane:	nCPAP PRESSURE LIMIT"
10.7.1.91	Was:	"LOW"	
10.7.1.134	Was:	"AlarmHistPpeakHighExt	
Description:		Indication that alarm has asserted in the past, is no longer active and has not been reset. See AlarmPpeakHighExt description above.	
Type:		BOOL	
Level:		High	
Label:		EXT HIGH Ppeak	
AVEA GUI/Membrane:		EXT HIGH Ppeak"	
10.7.2.3	Was:	"Label:	APNEA"
10.7.2.3	Was:	"VELA GUI/Membrane:	APNEA"
10.7.2.6	Was:	"Label:	CIRC FAULT"
10.7.2.6	Was:	"VELA GUI/Membrane:	CIRC FAULT"
10.7.2.7	Was:	"Level:	HIGH"
10.7.2.7	Was:	"Label:	LOW BATTERY"
10.7.2.7	Was:	"VELA GUI/Membrane:	LOW BATTERY"
10.7.2.19	Was:	"Label:	FAN FAULT"
10.7.2.19	Was:	"VELA GUI/Membrane:	FAN FAULT"
10.7.2.20	Was:	"Label:	O2 RANGE ERROR"
10.7.2.20	Was:	"VELA GUI/Membrane:	O2 RANGE ERROR"
10.7.2.21	Was:	"Label:	FLOW SENSOR DISCONNECT"

10.7.2.21	Was:	"VELA GUI/Membrane:	FLOW SENSOR DISCONNECT"
10.7.2.22	Was:	"Label:	HW FAULT"
10.7.2.22	Was:	"VELA GUI/Membrane:	HW FAULT"
10.7.2.24	Was:	"Label:	INVALID CONFIGURATION"
10.7.2.24	Was:	"VELA GUI/Membrane:	INVALID CONFIGURATION"
10.7.2.25	Was:	"Label:	LOW O2"
10.7.2.25	Was:	"VELA GUI/Membrane:	LOW O2"
10.7.2.28	Was:	"Label:	CHK O2 CAL"
10.7.2.28	Was:	"VELA GUI/Membrane:	CHK O2 CAL"
10.7.2.29	Was:	"Label:	HIGH O2"
10.7.2.29	Was:	"VELA GUI/Membrane:	HIGH O2"
10.7.2.30	Was:	"Label:	O2 SEN FAIL"
10.7.2.30	Was:	"VELA GUI/Membrane:	O2 SEN FAIL"
10.7.2.31	Was:	"Level:	MED"
10.7.2.34	Was:	"Label:	SUSTAINED HIGH PIP"
10.7.2.34	Was:	"VELA GUI/Membrane:	SUSTAINED HIGH PIP"
10.7.2.35	Was:	"Label:	LOW PRESS"
10.7.2.35	Was:	"VELA GUI/Membrane:	LOW PRESS"
10.7.2.36	Was:	"Label:	BATTERY ON"
10.7.2.36	Was:	"VELA GUI/Membrane:	BATTERY ON"
10.7.2.39	Was:	"Label:	HIGH BREATH"
10.7.2.39	Was:	"VELA GUI/Membrane:	HIGH BREATH"
10.7.2.41	Was:	"Level:	MED"
10.7.2.41	Was:	"Label:	LOW Ve"
10.7.2.41	Was:	"VELA GUI/Membrane:	LOW Ve"
10.7.2.42	Was:	"Label:	APNEA"
10.7.2.42	Was:	"VELA GUI/Membrane:	APNEA"
10.7.2.44	Was:	"Label:	CIRC FAULT"
10.7.2.44	Was:	"VELA GUI/Membrane:	CIRC FAULT"
10.7.2.45	Was:	"Label:	LOW BATTERY"
10.7.2.45	Was:	"VELA GUI/Membrane:	LOW BATTERY"
10.7.2.58	Was:	"Label:	FAN FAULT"
10.7.2.58	Was:	"VELA GUI/Membrane:	FAN FAULT"
10.7.2.59	Was:	"Label:	O2 RANGE ERROR"
10.7.2.59	Was:	"VELA GUI/Membrane:	O2 RANGE ERROR"
10.7.2.60	Was:	"Label:	FLOW SENSOR DISCONNECT"
10.7.2.60	Was:	"VELA GUI/Membrane:	FLOW SENSOR DISCONNECT"
10.7.2.61	Was:	"Label:	HW FAULT"
10.7.2.61	Was:	"VELA GUI/Membrane:	HW FAULT"
10.7.2.63	Was:	"Label:	INVALID CONFIGURATION"
10.7.2.63	Was:	"VELA GUI/Membrane:	INVALID CONFIGURATION"
10.7.2.64	Was:	"Label:	LOW O2"
10.7.2.64	Was:	"VELA GUI/Membrane:	LOW O2"
10.7.2.67	Was:	"Label:	CHK O2 CAL"
10.7.2.67	Was:	"VELA GUI/Membrane:	CHK O2 CAL"
10.7.2.68	Was:	"Label:	HIGH O2"

10.7.2.68	Was:	"VELA GUI/Membrane:	HIGH O2"
10.7.2.69	Was:	"Label:	O2 SEN FAIL"
10.7.2.69	Was:	"VELA GUI/Membrane:	O2 SEN FAIL"
10.7.2.70	Was:	"Level:	MED"
10.7.2.70	Was:	"Label:	PATIENT DEFAULTS"
10.7.2.70	Was:	"VELA GUI/Membrane:	PATIENT DEFAULTS"
10.7.2.73	Was:	"Label:	SUSTAINED HIGH PIP"
10.7.2.73	Was:	"VELA GUI/Membrane:	SUSTAINED HIGH PIP"
10.7.2.74	Was:	"Label:	LOW PRESS"
10.7.2.74	Was:	"VELA GUI/Membrane:	LOW PRESS"
10.7.2.75	Was:	"Label:	BATTERY ON"
10.7.2.75	Was:	"VELA GUI/Membrane:	BATTERY ON"
10.7.2.78	Was:	"Label:	HIGH BREATH"
10.7.2.78	Was:	"VELA GUI/Membrane:	HIGH BREATH"
10.7.2.80	Was:	"Level:	MED"
10.7.2.80	Was:	"Label:	LOW Ve"
10.7.2.80	Was:	"VELA GUI/Membrane:	LOW Ve"
Added:	10.7.2.81	ID:	AlarmSilence
		Description	Active/Inactive state of the capability to locally silence the audible ventilator alarms.
		Type:	BOOL
		Label:	ALARM SILENCE
		VELA GUI/Membrane:	ALARM SILENCE"
10.8.2.4	Was:	"Label:	Vt(ml)"
10.8.2.4	Was:	"VELA GUI/Membrane:	Vt(ml)"

### 13.3 Rev C

1.4	Removed the following references: VIASYS Open XML Protocol Specification Addendum – ER1794 Rev. A Specification, VOXP Avea Specialization – P/N 91316 Rev. C Specification, VOXP Vela Specialization – P/N 91415 Rev. A
5.5	Clarified the description for scalar example.
7.6	Added "TrendNcpapMeanFlow" and "Nasal Continuous Positive Airway Pressure Mean Flow."
7.6	was "TrendnCPAP"
10.1.1	Removed references to "Avea".
10.1.2	Added a limit to the number of waveforms and baud rate.
10.3	Removed references to class tables. The tables have been replaced with the text descriptions found in sections 10.4 – 10.8 prior to Rev A.
10.3	Added class priority scheme.
10.5.1.54	was "Label: ACC (On/Off)"
10.5.1.58	was "Label: Active (On/Off)"
10.5.1.60	was "Label: Leak Comp (On/Off)"
10.8.1.13	Added the parameter WaveVt.

Throughout the document, minor formatting and grammatical corrections per Marketing's request. For example, "Avea" was globally changed to "AVEA"